томатоло и и ортодонти

Технические характеристики

По вопросам продаж и поддержки обращайтесь:

Алматы (7273)495-231 Ангарск (3955)60-70-56 Архангельск (8182)63-90-72 Астрахань (8512)99-46-04 Барнаул (3852)73-04-60 Белгород (4722)40-23-64 Благовещенск (4162)22-76-07 Брянск (4832)59-03-52 Владивосток (423)249-28-31 Владикавказ (8672)28-90-48 Владимир (4922)49-43-18 Волоград (844)278-03-48 Волоград (8172)26-41-59 Воронеж (473)204-51-73 Екатеринбург (343)384-55-89

Иваново (4932)77-34-06 Ижевск (3412)26-03-58 Иркутск (395)279-98-46 Казань (843)206-01-48 Калининград (4012)72-03-81 Калуга (4842)92-23-67 Кемерово (3842)65-04-62 Киров (8332)68-02-04 Коломна (4966)23-41-49 Кострома (4942)77-07-48 Краснодар (861)203-40-90 Красноярск (391)204-63-61 Курск (4712)77-13-04 Курган (3522)50-90-47 Липецк (4742)52-20-81 Магнитогорск (3519)55-03-13 Москва (495)268-04-70 Мурманск (8152)59-64-93 Набережные Челны (8552)20-53-41 Нижний Новгород (831)429-08-12 Новокузнецк (3843)20-46-81 Ноябрьск (3496)41-32-12 Новосибирск (383)227-86-73 Омск (3812)21-46-40 Орел (4862)44-53-42 Оренбург (3532)37-68-04 Пенза (8412)22-31-16 Петрозаводск (8142)55-98-37 Псков (8112)59-10-37 Пермь (342)205-81-47 Ростов-на-Дону (863) 308-18-15 Рязань (4912) 46-61-64 Самара (846) 206-03-16 Санкт-Петербург (812) 309-46-40 Саратов (845) 249-38-78 Севастополь (8692) 22-31-93 Саранск (8342) 22-96-24 Симферополь (3652) 67-13-56 Смоленск (4812) 29-41-54 Сочи (862) 225-72-31 Ставрополь (8652) 20-65-13 Сургут (3462) 77-98-35 Сыктывкар (8212) 25-95-17 Тамбов (4752) 50-40-97 Тверь (4822) 63-31-35 Тольятти (8482)63-91-07 Томск (3822)98-41-53 Тула (4872)33-79-87 Тюмень (3452)66-21-18 Ульяновск (8422)24-23-59 Улан-Удь (3012)59-97-51 Уфа (347)229-48-12 Хабаровск (4212)92-98-04 Чебоксары (8352)28-53-07 Челябинск (351)202-03-61 Череповец (8202)49-02-64 Чита (3022)38-34-83 Якутск (4112)23-90-97 Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +7(7172)727-132

Киргизия +996(312)96-26-47

эл.почта: tam@nt-rt.ru || сайт: https://3m.nt-rt.ru/







3M[™] Polyether Adhesive Refill, 30600

3M[™] Scotchbond[™] Dispensing Tips, 41274

3M™ ESPE™ Express™ XT Putty Soft VPS Impression Material Refill, 36972







3M[™] Express[™] VPS Bite Registration Material Putty, 6160J

3M™ Vanish™ XT Extended Contact Varnish, 12148, 10 g Clicker™ Dispenser

3M™ Elipar™ DeepCure-S Li Ion Battery, 76985







3M™ Penta™ Mixing Tips Refill, Red (50), 77949

3M™ Filtek™ Z250 XT Nano Hybrid Universal Restorative Introductory Kit, 1470IK, Syringe

3M™ Unitek™ Stainless Steel Permanent Bicuspid Crowns, 902600, Kit







3M™ Pentamix™ 3 Wall Mount Kit, 77601

3M™ Imprint™ 4 VPS Bite Registration Material, 71529

3M™ Ketac™ Conditioner Refill, 37470







3M™ ESPE™ Stainless Steel Primary Molar Crowns

3M™ RelyX™ Unicem Self-Adhesive Universal Resin Cement Maxicap™ Capsule... Refill

3M™ Vitremer™ Core Buildup/Restorative, 3303







3M™ Polycarbonate Crowns

3M™ Filtek™ One Bulk Fill Restorative

3M™ RelyX™ Universal Resin Cement Refills







3M™ Filtek™ Easy Match Universal Restorative

3M™ Unitek ™ Anterior Band Regular Temper Narrow Upper Left/Right Central/Lateral

3M™ Imprint™ 4 Preliminary Penta™ Cartridge for Pentamix™ 3, 71520 Series








3M™ RelyX™ Universal Resin Cement Kits

3M[™] Unitek[™] Steiner Rotating Spring Wedge

3M™ Unitek™ Clear Plastic Button, 224-800, 10/Pk







3M™ Unitek™ Lingual Mathieu Pliers



3M™ Unitek™ Molar Band Seater, 811-001, 1/Ea








3M™ Unitek™ Single-Ended Ligature Director, 803-003, 1/Ea

3M™ AlastiK™ Color ID Dots, 406-829, 1/Pk

3M[™] Typodont with Stand, 3M[™] Clarity[™] Advanced Ceramic Brackets, 600-236,... 1/Ea







3M™ Unitek™ Coaxial Wire

3M[™] Concise[™] Orthodontic Chemical Cure Adhesive Enamel Bond Resin A, 196-00... 6 ml, 1/Btl

3M™ Unitek™ Preformed Lingual Arch, Horizontal







3M™ Unitek™ Molded Retainer Cases

3M™ Unitek™ Bracket Positioning Gauge

3M™ Clarity™ Retainers







3M™ Clinpro™ Clear Fluoride Treatment

3M™ Clinpro™ Sealant Refill, 12600 Series

3M[™] Unitek[™] Adjustable Sliding Yoke, 310-200, 10/Pk






3M[™] Unitek[™] Prestige Lingual Lock-Forming Pliers, 900-706, 1/Ea

3M™ AlastiK™ QuiK-StiK™ Separators, 406-410 Series

704-045 Unitek™ Disposable Spatula 100/Pk







3M™ Unitek™ T Spring

3M™ Unitek™ Eyelet, Large, 471-003, 1/Ea

3M™ Unitek™ Lingual Light Wire Flush Cut Distal End Cutter, 803-550, 1/Ea







3M™ Deluxe AlastiK™ Force Modules Needle Holder, 803-535, 1/Ea

3M™ Unitek™ Prestige Posterior Band Removing Pliers

3M™ Unitek™ Single-Ended Ligature Director, 803-003, 1/Ea







3M[™] AlastiK[™] Color Palette, 406-825, 1/Ea

3M[™] Victory Series[™] 2nd Bicuspid Tubes

3M[™] Adper[™] Scotchbond[™] Multi-Purpose Primer








3M™ Ketac™ Cem Maxicap™ Permanent Glass Ionomer Luting Cement Intro Kit, 56015

3M™ Unitek™ Bicuspid Bands Kit

3M™ RelyX™ Universal Micro Mixing Tips







3M™ Digital Bonding Tray

3M™ Unitek™ Braided Archwire

3M™ Concise™ Orthodontic Bonding Paste B, 18 g, 196-003, 1/Ea







3M™ Series 5 Facebows

3M™ Unitek™ Multi-Cure Glass Ionomer Band Cement Powder, 712-051, 35 g, 1/Btl

3M™ Unitek™ Multi-Cure Glass Ionomer Orthodontic Band Cement Dispensing Spoon,... 712-053, 1/Ea







3M™ Forsus™ Fatigue Resistant Device, Split Crimps, 885-203, 50/Pk

3M™ Forsus™ Fatigue Resistant Device, 885-100, L-Pin Spring Modules, 5/Pk

3M™ Unitek™ General Purpose Molar Bands Kit








3M™ Unitek™ Lingual Anterior Debonding Pliers, 803-532, 1/Ea


3M[™] Unitek[™] Lingual Retainer Kit, 709-060, 1 Kit

3M™ Unitek™ Posterior Band Removing Pliers, 801-713, 1/Ea





ск 000000000000



3M[™] Unitek[™] Auxiliary Ligature Tie Hook 390-145, 0.014, 100/Pk

3M™ AlastiK™ Continuous C Module Chain Reel

3M™ Unitek™ Bracket Positioning Gauge Five, 900-840, 6.0-6.5 mm, 1/Ea







956-208 CNTRD STD ME U7

3M[™] Sof-Lex[™] Diamond Polishing Spiral, 5091

3M™ Unitek™ Permachrome Resilient Archwire







Lava™ Frame Shading Liquid

3M™ Paradigm™ MZ100 Blocks for CEREC®

3M™ Polycarbonate Crowns, C-180, Anterior, Kit







3M™ Penta™ Cartridge for Pentamix™ Lite, 77944

3M™ Scotchbond™ Etchant, 7423

3M[™] Adper[™] Scotchbond[™] 1 XT Adhesive








3M™ ESPE™ Ketac™ Molar Easymix Glass Ionomer Filling Material Refill, 56625

3M™ Protemp™ Crown Temporization Material Kit

3M™ Clinpro™ Tooth Créme 0.21% Sodium Fluoride Anti-Cavity Toothpaste







3M™ Unitek™ Twist Wire

3M™ Concise™ Orthodontic Bonding Paste A, 18 g, 196-002, 1/Ea

3M™ Victory Series™ First Molar Bands Kit







3M[™] Penta[™] Universal Cartridge for Pentamix[™] 3, 71545

3M™ Filtek™ Restoratives Dispenser, 5850RD, 1 Dispenser

3M™ Forsus™ Fatigue Resistant Device, 885-200, Reactivation Spacer, 10/Pk







3M[™] Victory Series[™] V-Slot Brackets Kit

3M™ Garant™ Intraoral Tips, Intraoral Dispensing

3M[™] Transbond[™] LR Light Cure Adhesive







3M™ Impregum™ Super Quick Polyether Impression Material

3M™ RelyX™ Ultimate Clicker™ Adhesive Resin Cement Kit

3M™ Ortholux™ Luminous Curing Light Guide, 704-456, 1/Ea



3M™ Astringent Retraction Paste



3M™ CoJet™ Sand, 68411, 3 Bottles



3M[™] Ortholux[™] Luminous Curing Light Glasses, S0360X, 1/Ea



3M™ Forsus™ Fatigue Resistant Device, 885-110, L-Pins Spring Module, 10/Pk




3M™ Unitek™ General Purpose Molar Bands



3M™ Unitek™ Lingual Distal End Cutter, 803-552, 1/Ea







3M™ Unitek™ How Utility Offset Pliers, 801-701, 1/Ea

3M™ Unitek™ Prestige Series Omega Loop Forming Pliers Replacement Tip, 900-708, 1/Ea

3M[™] Unitek[™] Ligature Wire, 390-160, Preformed Coated, Natural Colored, .012 in., 100/Pk







3M)™ AlastiK™ Plastic Box With Dividers, 406-890, 1/Ea

3M™ Filtek™ Flowable Dispensing Tips

3M™ Photac™ Fil Quick Aplicap™ Light-Cured Glass Ionomer Restorative, 61000... Series







3M™ Victory Series™ Second Molar Bands

3M™ Ketac™ Cem Glass Ionomer Luting Cement Refill, 56000 Series

3M™ Ketac™ Nano Light-Curing Glass Ionomer Quick Mix Restorative Trial Kit, 3305TK







3M[™] RelyX[™] Fiber Post 3D Glass Fiber Post

3M™ Unitek™ Beta III Titanium Archwire

453-012 Unitek™ Wire Dispenser Ea







3M™ SmartClip™ SL3 Self-Ligating Brackets Kit

3M™ AlastiK™ B Modules

3M™ Unitek™ Sickle-Style Cement Remover and Bracket Positioner, 900-822, 1/Ea







3M™ Imprint™ 4 Penta™ Heavy Cartridge for Pentamix™ 3, 77784

3M™ Unitek™ Prestige End Cutter, 900-750, 1/Ea

3M™ Penta™ Elastomer Syringe Knurled Hub Cap, 77720







3M™ ESPE™ Intraoral Syringe Value Pack

3M[™] Applicator Handles and Disposable Applicator Brush Tips, 1919, (3 Handles/120 Tips)

3M™ Crown Scissor







3M™ Lava™ Ultimate CAD/CAM Restorative for CEREC®



3M[™] Sof-Lex[™] Pre-Polishing Spiral, 5090, 15ct



3M[™] Unitek[™] Nitinol Heat-Activated Archwire



3M™ Stainless Steel Permanent Molar Crowns



3M™ Imprint™ II Monophase VPS Impression Material







3M™ Imprint™ Bite Registration Material Refill, 36850

3M[™] Adper[™] Single Bond Plus Adhesive Intro Kit, 51101

3M™ RelyX™ Veneer Cement Syringe Refill







3M[™] Durelon[™] Carboxylate Luting Cement Pack

3M™ Clinpro™ Sealant Introductory Kit, 12600 Series

3M™ Unitek™ Miniature Twin Metal Brackets Kit







3M™ Clinpro™ Sealant & 3M™ Adper™ Prompt™ L-Pop™ Self-Etch Adhesive Trial Kit,... 12626TP, 1.2 mL/Syringe, 1 Syringe/Kit

3M™ Paradigm™ Temporization Material Refill

3M™ Lava™ Ultimate CAD/CAM Restorative Intro Kit for CEREC®, 2914SK







3M™ Strip Crown Form - Pedo

3M™ Unitek™ Stainless Steel Primary Molar Crowns

3M™ Adper™ Single Bond 2 Intro Kit, 51201







3M™ RelyX™ Try-In Paste Syringe Refill

3M[™] Durelon[™] Carboxylate Luting Cement Maxicap Refill, 56310

3M™ RelyX™ Ceramic Primer







3M™ Victory Series™ Superior Fit Buccal Tubes

3M[™] Victory Series[™] Active Self-Ligating Brackets

3M[™] Victory Series[™] V-Slot Brackets








3M™ Transbond™ MIP Moisture Insensitive Primer Kit, 712-021, 1/Kit

3M™ Transbond™ Plus Color Change Adhesive

3M™ Mixing Sticks (50), 1994S







3M™ Disposable Applicator Brush Tips, 1919B, 60ct

3M™ RelyX™ Fiber Post 3D Intro Kit, 56957

3M™ Unitek™ Ligature Wire







3M™ Imprint™ 4 Penta™ Cartridge for Pentamix™ 3, 77780 Series

3M™ Unitek™ Stainless Steel Permanent Molar Crowns

3M™ Impregum™ Penta™ Cartridge for Pentamix™ 3







3M)™ AlastiK™ Plastic Box With Dividers, 406-890, 1/Ea

3M™ Filtek™ Flowable Dispensing Tips

3M™ Photac™ Fil Quick Aplicap™ Light-Cured Glass Ionomer Restorative, 61000... Series







3M™ Victory Series™ Second Molar Bands

3M™ Ketac™ Cem Glass Ionomer Luting Cement Refill, 56000 Series

3M™ Ketac™ Nano Light-Curing Glass Ionomer Quick Mix Restorative Trial Kit, 3305TK







3M™ RelyX™ Fiber Post 3D Glass Fiber Post

3M™ Unitek™ Beta III Titanium Archwire

453-012 Unitek™ Wire Dispenser Ea



3M™ SmartClip™ SL3 Self-Ligating Brackets Kit



3M™ AlastiK™ B Modules



3M[™] Unitek[™] Sickle-Style Cement Remover and Bracket Positioner, 900-822, 1/Ea



3M™ Imprint™ 4 Penta™ Heavy Cartridge for Pentamix™ 3, 77784



3M[™] Unitek[™] Prestige End Cutter, 900-750, 1/Ea



3M™ Penta™ Elastomer Syringe Knurled Hub Cap, 77720



3M™ ESPE™ Intraoral Syringe Value Pack




3M[™] Applicator Handles and Disposable Applicator Brush Tips, 1919, (3 Handles/120 Tips)



3M™ Crown Scissor







3M™ Lava™ Ultimate CAD/CAM Restorative for CEREC®



3M[™] Sof-Lex[™] Pre-Polishing Spiral, 5090, 15ct



3M[™] Unitek[™] Nitinol Heat-Activated Archwire



3M™ Stainless Steel Permanent Molar Crowns



3M™ Imprint™ II Monophase VPS Impression Material







3M™ Imprint™ Bite Registration Material Refill, 36850

3M[™] Adper[™] Single Bond Plus Adhesive Intro Kit, 51101

3M™ RelyX™ Veneer Cement Syringe Refill







3M™ Durelon™ Carboxylate Luting Cement Pack

3M™ Clinpro™ Sealant Introductory Kit, 12600 Series

3M™ Unitek™ Miniature Twin Metal Brackets Kit







3M™ Clinpro™ Sealant & 3M™ Adper™ Prompt™ L-Pop™ Self-Etch Adhesive Trial Kit,... 12626TP, 1.2 mL/Syringe, 1 Syringe/Kit

3M™ Paradigm™ Temporization Material Refill

3M[™] Lava[™] Ultimate CAD/CAM Restorative Intro Kit for CEREC®, 2914SK







3M™ Strip Crown Form - Pedo

3M™ Unitek™ Stainless Steel Primary Molar Crowns

3M[™] Adper[™] Single Bond 2 Intro Kit, 51201







3M™ RelyX™ Try-In Paste Syringe Refill

3M™ Durelon™ Carboxylate Luting Cement Maxicap Refill, 56310

3M™ RelyX™ Ceramic Primer







3M[™] Victory Series[™] Superior Fit Buccal Tubes

3M™ Victory Series™ Active Self-Ligating Brackets

3M™ Victory Series™ V-Slot Brackets








3M™ Transbond™ MIP Moisture Insensitive Primer Kit, 712-021, 1/Kit

3M[™] Transbond[™] Plus Color Change Adhesive

 $3M^{TM}$ Mixing Sticks (50), 1994S







3M[™] Disposable Applicator Brush Tips, 1919B, 60ct

3M™ RelyX™ Fiber Post 3D Intro Kit, 56957

3M™ Unitek™ Ligature Wire








3M™ Imprint™ 4 Penta™ Cartridge for Pentamix™ 3, 77780 Series

3M™ Unitek™ Stainless Steel Permanent Molar Crowns

3M™ Impregum™ Penta™ Cartridge for Pentamix™ 3







3M™ Unitek™ Gurin Intraoral Hook

3M™ Unitek™ Lingual Heavy Wire Flush Cut Distal End Cutter, 803-551, 1/Ea

3M™ Transbond™ Plus Self Etching Primer







3M™ Transbond™ Plus Disposable Tips, 712-082, 25/Pack

3M™ Unitek™ Lingual Sheaths

3M[™] Unitek[™] Prestige Hard Wire Cutter, 900-751, 1/Ea







3M[™] Unitek[™] Prestige Wire Cutter, 900-753, 1/Ea

3M™ Chairside Zirconia for CEREC®

3M[™] Dual Arch Impression Tray







3M™ Ortholux™ Luminous Curing Light Eye Shield, 704-457, 1/Ea

3M™ Rocatec™ Junior Bonding System Refill

3M™ Victory Series™ Brackets







3M[™] Unitek[™] Multi-Cure Glass Ionomer Band Cement, 712-050, 1/Kit

3M[™] Victory Series[™] Brackets Kit

3M[™] Clarity[™] Ultra Self-Ligating Brackets Kit







3M™ Impregum™ Penta™ Super Quick Medium Body/Light Body Polyether Impression... Material Intro Kit

3M™ Fast Release 5% Sodium Fluoride Varnish, 12150 Series

3M[™] Silane Coupling Agent, 68300, 1 per Pack







3M™ Elastomer Syringe Washer Refill, 71050



3M™ Mixing Pad Refill







3M[™] Fiber Tip Applicators, 1919F, 100ct

3M[™] Vitremer[™] Core Buildup/Restorative, 3303

3M[™] Unitek[™] Prestige Hollow Chop Pliers, 900-716, 1/Ea







3M™ ESPE™ Ketac™ Cem Radiopaque Glass Ionomer Luting Cement Intro Kit, 37201

3M™ Transbond™ Plus Light Cure Band Adhesive Kit

3M™ Unitek™ Permachrome Posted Archwire







3M[™] Penta[™] Elastomer Nozzle, 71225

3M™ RotoMix™ Capsule Mixing Unit, 76310, 120 V/60 Hz

3M[™] Unitek[™] MIA Mobile Intraoral Arch, 080-810, U/L, Curved, 1/Ea







3M™ Lava™ Zirconia Blocks Refill, 68583, 20mm

3M™ Unitek™ Stainless Steel Primary Anterior Crowns

3M™ RelyX™ Unicem Self-Adhesive Universal Resin Cement Aplicap™ Capsule Kit







3M™ RelyX™ Temp Zinc Oxide Non-Eugenol Temporary Cement Refill

3M™ TheraMints™ 100% Xylitol Sweetened Mints, 12119M, Mint Flavor, 520 Pieces Per Bottle

3M™ Victory Series™ First Molar Bands







3M™ Elastomer Standard Spatula Refill, 71233

3M™ Elastomer Syringe Screw Refill, 71076

3M™ Brush Holders







3M™ Elastomer Syringe and Filling Device WE/EE, 71209

3M™ Disposable Mini-Sponge Applicators

3M[™] Vitremer[™] Core Buildup/Restorative Primer, 3303P, 6.5 mL







3M™ Transbond™ Supreme LV Adhesive

3M™ Filtek™ Supreme Ultra Universal Restorative Single Shade Intro Kit

3M™ Sof-Lex™ Diamond Polishing System, 5092










3M™ Sof-Lex™ Square Eyelet Discs, 1958

3M™ Unitek™ Gurin Lock

3M[™] Unitek[™] Action Line Tubing








3M™ MBT™ Appliance System Bracket Positioning Guide, 701-725, Plastic, 1/Ea

3M™ Unitek™ Combination Lingual Hook/Seating Lugs

3M™ Unitek™ Prestige Jaraback Pliers, 900-704, 1/Ea







3M™ Unitek™ Gurin Intraoral Hook

3M™ Unitek™ Lingual Heavy Wire Flush Cut Distal End Cutter, 803-551, 1/Ea

3M™ Transbond™ Plus Self Etching Primer







3M™ Transbond™ Plus Disposable Tips, 712-082, 25/Pack

3M™ Unitek™ Lingual Sheaths

3M™ Unitek™ Prestige Hard Wire Cutter, 900-751, 1/Ea







3M™ Unitek™ Prestige Wire Cutter, 900-753, 1/Ea

3M™ Chairside Zirconia for CEREC®

3M™ Dual Arch Impression Tray







3M™ Ortholux™ Luminous Curing Light Eye Shield, 704-457, 1/Ea

3M™ Rocatec™ Junior Bonding System Refill

3M[™] Victory Series [™] Brackets







3M[™] Unitek[™] Multi-Cure Glass Ionomer Band Cement, 712-050, 1/Kit

3M™ Victory Series™ Brackets Kit

3M™ Clarity™ Ultra Self-Ligating Brackets Kit







3M™ Impregum™ Penta™ Super Quick Medium Body/Light Body Polyether Impression... Material Intro Kit

3M™ Fast Release 5% Sodium Fluoride Varnish, 12150 Series

3M[™] Silane Coupling Agent, 68300, 1 per Pack



3M™ Victory Series™ Active Self-Ligating Brackets Kit



3M™ Elastomer Syringe Washer Refill, 71050



3M™ Mixing Pad Refill







3M[™] Fiber Tip Applicators, 1919F, 100ct

3M[™] Vitremer[™] Core Buildup/Restorative, 3303

3M™ Unitek™ Prestige Hollow Chop Pliers, 900-716, 1/Ea








3M™ ESPE™ Ketac™ Cem Radiopaque Glass Ionomer Luting Cement Intro Kit, 37201

3M[™] Transbond[™] Plus Light Cure Band Adhesive Kit

3M[™] Unitek[™] Permachrome Posted Archwire







3M[™] Penta[™] Elastomer Nozzle, 71225

3M™ RotoMix™ Capsule Mixing Unit, 76310, 120 V/60 Hz

3M™ Unitek™ MIA Mobile Intraoral Arch, 080-810, U/L, Curved, 1/Ea







3M[™] Lava[™] Zirconia Blocks Refill, 68583, 20mm

3M™ Unitek™ Stainless Steel Primary Anterior Crowns

3M™ RelyX™ Unicem Self-Adhesive Universal Resin Cement Aplicap™ Capsule Kit







3M™ RelyX™ Temp Zinc Oxide Non-Eugenol Temporary Cement Refill

3M[™] TheraMints[™] 100% Xylitol Sweetened Mints, 12119M, Mint Flavor, 520 Pieces Per Bottle

3M™ Victory Series™ First Molar Bands







3M™ Elastomer Standard Spatula Refill, 71233

3M™ Elastomer Syringe Screw Refill, 71076

3M™ Brush Holders







3M™ Elastomer Syringe and Filling Device WE/EE, 71209

3M™ Disposable Mini-Sponge Applicators

3M[™] Vitremer[™] Core Buildup/Restorative Primer, 3303P, 6.5 mL







3M™ Transbond™ Supreme LV Adhesive

3M™ Filtek™ Supreme Ultra Universal Restorative Single Shade Intro Kit

3M[™] Sof-Lex[™] Diamond Polishing System, 5092










3M™ Sof-Lex™ Square Eyelet Discs, 1958

3M™ Unitek™ Gurin Lock

3M™ Unitek™ Action Line Tubing







3M™ MBT™ Appliance System Bracket Positioning Guide, 701-725, Plastic, 1/Ea

3M™ Unitek™ Combination Lingual Hook/Seating Lugs

3M™ Unitek™ Prestige Jaraback Pliers, 900-704, 1/Ea







3M[™] Iso-Form[™] Temporary Crowns

3M™ Unitek™ Prestige Universal Cutter

3M™ Impregum™ Soft Polyether Impression Material







3M[™] Scotchbond[™] Universal Adhesive Unit Dose Intro Pack

3M™ RelyX™ Unicem 2 Automix Self-Adhesive Resin Cement Value Pack

3M™ PerioMed™ 0.63% Stannous Fluoride Oral Rinse







3M[™] Filtek[™] Universal Restorative, 6500 Series

3M™ Unitek™ 804-155 Open-Close Instrument

3M™ Clarity™ Aligners







3M™ Unitek™ Miniature Twin Metal Brackets

3M[™] Concise[™] Orthodontic Chemical Cure Adhesive Bonding System, 196-001, 1/Kit

3M™ Unitek™ Miniature Twin V-Slot Hook Pins, 461-140, 100/Pk







3M™ Sof-Lex™ XT Pop-on Extra Thin Disc Refills - Orange Series, 30 x 12.7mm, Superfin... 4931SF

3M™ Sof-Lex™ XT Pop-on Extra Thin Disc Refills - Orange Series, 30 x 12.7mm, Coars... 4931C

3M™ Just For Kids™ 0.4% Stannous Fluoride Brush-On Gel Refill, 12107







3M™ Unitek™ Multi-Cure Glass Ionomer Band Cement Liquid, 712-052, 25 g, 1/Btl

3M™ Scotchbond™ Universal Plus Adhesive Unit Doses L-Pop

3M™ SmartClip™ Appliance Wire Disengagement Hand Instrument, 804-162, 1/Ea







3M[™] Unitek[™] Hi-T II Stainless Steel Coil Retraction Spring

3M™ Unitek™ Level Anchorage Palatal Bar

3M™ Unitek™ Low Profile Lingual Cleats







3M[™] Unitek[™] Molar Band Seater, 811-001, 1/Ea

3M™ Unitek™ Lingual Cinch Back Pliers, 803-531, 1/Ea

3M™ APC™ Reusable 7x7 Set-Up Tray, 2600-005, 1/Ea











3M™ AlastiK™ Loose Modules

3M™ AlastiK™ Rotation Wedge

3M™ Transbond™ Plus Self Etching Primer, 712-093, Easy Roller, 1/Ea







3M™ AlastiK™ Mini-StiK Patient Sticks

3M™ Unitek™ Stainless Steel Primary Anterior Crowns, 908100, Kit

3M™ SmartClip™ Appliance Wire Disengagement Hand Instrument Slide Reorder,... 804-170, 5/Pk







3M™ Unitek™ Etching Gel System (Blue Tint) Kit

3M™ Strip Crown Form - Pedo, 915100, Kit

3M™ Elastomer Syringe Knurled Hub Cap Refill, 71043







3M[™] Adper[™] Scotchbond[™] Multi-Purpose 3 Well Dispenser, 7544

3M™ Delivery Tips With Pistons, 3303DT, 50ct

3M[™] XS Applicator Brushes, 41904, 50 ct







3M™ Sof-Lex™ Mandrel Refill

3M[™] Scotchbond[™] Universal Plus Adhesive Kits

3M™ Unitek™ Retraction Spring







3M[™] Wax-Pak[™] Sticks, 603-063, 100/Pk

3M™ Unitek™ Bondable Lingual Cleat, 709-050, 10/Pk

3M™ Unitek™ Double-Ended Ligature Director, 803-530, 1/Ea










3M™ Unitek™ Lingual Weingart Utility Pliers, 803-530, 60 Degree Bend, 1/Ea

3M™ Unitek™ Bracket Placement Instrument

3M™ Unitek™ Nitinol Open Coil Spring







3M™ AlastiK™ Separator Modules, 406-080 Series

3M™ AlastiK™ C Modules

3M™ Filtek™ Supreme Flowable Restorative, Syringe







3M™ Unitek™ Preformed Vertical Lingual Archwire with Loop

3M[™] Adper[™] Single Bond Plus Adhesive Refill, 51102

3M[™] Concise[™] Orthodontic Chemical Cure Adhesive Enamel Bond Resin B, 196-01... 6 ml, 1/Btl







3M™ APC™ Inventory Dispensing System, Back-Up Storage Units, 7 x 7, 2600-03... 1/Ea

3M™ AlastiK™ QuiK-StiK™ Ligatures

3M™ Victory Series™ Mini Tubes







3M[™] Single Bond Universal Adhesive Vial Refill

3M[™] Chairside Zirconia for CEREC® Kit, 69956

3M[™] Polyether Contact Tray Adhesive Refill, 69408







3M™ Filtek™ Bulk Fill Flowable Restorative, Capsule

3M™ Ortholux™ Luminous Curing Light Rechargeable Battery, 704-455, 1/Ea

3M™ RelyX™ Ultimate Adhesive Resin Cement Refill







3M™ Unitek™ Prestige Arch Bending Pliers, 900-721, 1/Ea

3M™ Victory Series™ Low Profile Brackets Kit

3M™ Unitek™ Nitinol Classic Archwire








3M[™] Unitek[™] Clarity[™] Ultra SL Open-Close Instrument, 804-155, 1/Ea

3M™ Unitek™ Gurin Wrench

3M[™] Unitek[™] Gurin Loop Spring







3M™ OrthoForm™ Operatory Arch Form Templates, White, 701-724, Ea

3M™ Unitek™ Raised Arm Lingual Cleats

3M[™] Unitek[™] Prestige How Utility Pliers, Os Tip, 900-702, 1/Ea







3M™ Unitek™ Prestige How Utility Pliers

3M™ APC™ Inventory Dispensing System, 2600-022, Drawer Dust Cover Set, 1/Ea

3M™ Unitek™ Nitinol Super-Elastic Archwire







3M™ AlastiK™ Lingual Ligatures, 406-993, Navy Blue, 12/Stick, 12/Pk

3M™ AlastiK™ Bobbin Chain Dispenser, 406-927, 1/Ea

3M™ Unitek™ Ball Clasp







3M™ Transbond™ MIP Moisture Insensitive Primer, 712-025, 6 ml, 1/Btl

3M™ Unitek™ Nitinol Closed Coil Spring

3M[™] Clarity[™] Advanced Ceramic Brackets Kit



3M™ Unitek™ Stainless Steel Permanent Molar Crowns, 902350, Kit



3M™ Typodonts



3M™ Mini Uni-Twin™ Brackets







3M™ Filtek™ Universal Restorative Essentials Kit

3M[™] RelyX[™] Unicem Elongation Tips for 3M[™] Aplicap[™] Capsules, 77550

3M™ Penta™ Plunger Discs, 77780







3M™ Vitrebond™ Plus Light Cure Glass Ionomer Liner/Base Refill, 7570 Series

3M[™] Elipar[™] S10 Eye Shield, 76984, 1 - Eye Shield

3M[™] Elipar[™] DeepCure-S Light Guide, 76981, 1 - 10mm Guide







3M[™] Giz Accessory Pack Maxicap, 37180

3M™ Ketac™ Cem EasyMix Triple Kit EE, 56906

3M™ Ketac™ Conditioner Refill, 37470







3M™ Ketac™ Nano Glass Ionomer Primer, 3304P, 6.5ml

3M™ Paradigm™ DeepCure Light Guide, 76982, 1 - 10mm Light Guide

3M™ ESPE™ Stainless Steel Permanent Molar Crowns, PO-96, Kit







3M™ Unitek™ Miniature Twin V-Slot Metal Brackets

3M™ Filtek™ Supreme Ultra Universal Restorative Master's Kit, 6029M, Capsule

3M™ Sof-Lex™ Finishing and Polishing Discs, Kit, 1980








3M[™] Impregum[™] Polyether Impression Materials

3M[™] Scotchbond[™] Universal DCA Dual Cure Activator - Vial, 41264, 5ml

3M™ RelyX™ U200 Self-Adhesive Resin Cement Refill, 56878, A2







3M™ RelyX™ Luting Plus Cement

3M[™] Garant[™] Dispenser 4:1/10:1 Plunger Refill, 77591

3M[™] Bendable Applicator Brushes, 3402, 120ct







3M™ Giz Accesory Pack Aplicap, 37160

3M™ Vitrebond™ Light Cure Glass Ionomer Liner/Base Refill, 7510 Series

3M™ Paradigm™ DeepCure LED Curing Light, 76974, 120 V NA







3M™ Paradigm™ 10 mm Light Guide, 76968

3M™ Disposable Wells

3M™ Lava™ Milling Bur







3M™ Filtek™ One Bulk Fill Restorative

3M™ Polycarbonate Crowns

3M™ Pentamix™ Mounting Kit Twist Grip, 950916







3M[™] Penta[™] Universal Cartridge for Pentamix[™] 2, 76994

3M™ Express™ VPS Impression Materials

3M™ Ketac™ Cem Plus Luting Cement Refills, 56930







3M™ Durelon™ Carboxylate Luting Cement Intro Kit, 38019

3M™ Protemp™ Plus Temporization Material, 46962, Shade A2, Introductory Kit

3M™ Vanish™ 5% Sodium Fluoride White Varnish, 12100 Series







3M[™] Clinpro[™] 5000 1.1% Sodium Fluoride Anti-Cavity Toothpaste, 12000 Series

712-022 Unitek™ Light Sensitive Wells 4/Pk

3M™ Incognito™ Appliance System



3M™ Victory Series™ Buccal Tubes



3M™ Lava™ Esthetic Fluorescent Full-Contour Zirconia Disc



3M™ Gold Anodized Temporary Crowns, 940000 Series







3M™ Pentamix™ 3 Automatic Mixing Unit, 77876

3M™ Express™ STD VPS Impression Material Putty Refill

3M™ Clinpro™ XT Varnish Durable Fluoride-releasing Coating







3M™ OMNI Gel™ 0.4% Stannous Fluoride Brush-On Gel Refill, 12106

3M™ Lava™ Plus High Translucency Zirconia Disc

3M[™] Filtek[™] P60 Posterior Restorative







3M™ Paradigm™ Nano Hybrid Universal Restorative, 1480 Series

3M™ Sof-Lex™ Extra-Thin Contouring and Polishing Discs Kit, 2380

3M™ Pentamix™ Lite Automatic Mixing Unit - 120V/60Hz NA, 77907







3M™ Imprint™ 4 Preliminary VPS Impression Material

3M™ Adper™ Prompt™ L-Pop™ Self-Etch Adhesive Refill, 41925, 40ct

3M™ RelyX™ Luting Plus Cement Kit







3M™ Protemp™ Plus Temporization Material Refill

3M™ Peridex™ Chlorhexidine Gluconate 0.12% Oral Rinse

3M™ Transbond™ IDB Pre-Mix Chemical Cure Adhesive Kit, 712-120, 1/Kit







3M™ Sondhi™ Rapid-Set Indirect Bonding Adhesive, 712-070, 1/Kit

3M™ Victory Series™ Low Profile Bracket System

3M™ SmartClip™ SL3 Self-Ligating Brackets







3M™ Imprint™ 3 VPS Impression Material

3M™ Crown Pliers

3M™ Ketac™ Silver Aplicap™ Silver Glass Ionomer Restorative, 37000 Series







3M™ Unitek™ Silver Solder

3M[™] Unitek[™] Pliers/Cutter Combo, 800-705, 1/Ea

3M™ Unitek™ MIA Mobile Intraoral Arch Angle, 340-018, 5/Pk







3M™ Arch-Making Pencil Kit, 603-026, 1/Kit

3M[™] Unitek[™] Twin Torque[™] Brackets

3M™ Ketac™ Nano Light-Curing Glass Ionomer Quick Mix Capsule Refill, 3305







3M™ Band Removing Pliers, Replacement Pads

3M™ Typodont Hinge Replacements, 600-209, 1/Ea

3M™ SmartClip™ Appliance Detail Model (15X), 600-226, 1 Ea







3M[™] Unitek[™] Prestige Bond Removing Pliers, Replacement Pads, 6 mm., 900-719, 4/Pk

3M™ AlastiK™ Separator Radiopaque, 406-080 Series

3M™ Cavit™ Temporary Filling Material







3M™ Pentamix™ Lite Starter Kit
— Paradigm™ Penta™

3M™ Victory Series™ Low Profile Brackets

3M[™] Unitek[™] Prestige Omega Loop Forming Pliers, 900-707, 1/Ea







3M™ Unitek™ Gemini Clear Ceramic Brackets

3M™ Paradigm™ VPS Impression Material

3M™ Ketac™ Fil Plus Aplicap™ Glass Ionomer Restorative Refill, 55000 Series







3M[™] Dispensing Wells, 712-074, 2-cavity, 4/Pk

3M™ Unitek™ Archwire Hooks and Stops Wire Split Crimpable

3M™ Pentamix™ Lite Starter Kit
— Position™ Penta™







3M™ Pentamix™ Lite Starter Kit
— Imprint™ 4 Penta™

3M™ RelyX™ Ultimate Adhesive Resin Cement Mixing Tips

3M[™] Ketac[™] Universal Aplicap[™] Glass lonomer Restorative, 61000 Series







3M™ Visio™ Beta Vario Replacement Light, 77170

3M™ Pentamix™ 3 Precision Starter Kit — Imprint™ 4 Penta™

3M[™] Gold Anodized Crowns Kit, 942001 Series







3M™ Protemp™ Crown Temporization Material

3M™ Dimension™ VPS Impression Materials

3M™ Grommets for Eye Shield







Ketac™ Molar Quick Aplicap™ Glass Ionomer Restorative

3M™ Pentamix™ 3 Starter Kit — Impregum™ Polyether

3M™ Sinfony™ Indirect Lab Composite Master System, 49962







3M™ Lava™ Sintering Beads, 68594, 100g

3M™ ESPE™ Sof-Lex™ Finishing Strips Refill, Fine/Superfine, 1956

3M™ Unitek™ Prestige Bond Removing Pliers, Replacement Blade, 900-718, 1/Ea







3M[™] Unitek[™] Lingual Retainer Wire, 260-032, 7 in., .032, Bondable, 25/Pk

3M[™] Forsus[™] EZ2 Spring Modules

3M™ Unitek™ Lateral Development Archwire Arch Evaluator, 701-726, 5/Pk







3M[™] Unitek[™] Archwire Hooks and Stops Wire Crimpable

3M™ Forsus™ Push Rod Measurement Guide, 807-014, 5/Pk

3M™ Light-Tight Box, 712-075, 1/Ea







3M™ Unitek™ Prestige Bond Removing Pliers, 900-717, 1/Ea

3M™ Indirect Bonding Technique DVD

3M™ Visio™ Vacuum Pump Oil Refill, 77150







3M™ Unitek™ Gemini Brackets Kit

3M™ Lava™ Plus High Translucency Zirconia Effect Shade

3M[™] Adper[™] Scotchbond[™] Multi-Purpose Adhesive System Intro Kit, 7540S







3M™ Ketac™ Molar Aplicap™ Glass Ionomer Restorative Refill, 56400 Series

3M™ Vitremer™ Core Buildup/Restorative Kit, 3303

3M™ Elipar™ I Ellipoid Lamp, 77160







3M™ Sinfony™ Material Shade Dispenser Refills

3M™ Unitek™ Convertible Cap Removing Pliers Tip 805-036 5/Pack

3M™ Permadyne™ 2:1 Syringe Impression Material Refill Pack







3M™ Transbond™ XT Light Cure Adhesive

3M™ RelyX™ Unicem Self-Adhesive Universal Resin Cement Maxicap™ Capsule... Intro Kit

3M™ RelyX™ Luting Cement Intro Kit, 3505







3M[™] OrthoForm[™] Diagnosis Arch Form Templates, Clear, 701-723, 1/Ea

3M™ ESPE™ Stainless Steel Primary Molar Crowns, ND-96, Kit

3M™ Clarity™ Ultra Self-Ligating Brackets







3M[™] Unitek[™] Prestige Tweed Loop Forming Pliers, 900-720, 1/Ea

3M™ Pentamix™ Lite Starter Kit — Imprint™ 4 Penta™ Putty

3M[™] Sof-Lex[™] Contouring and Polishing Discs, 1982







3M™ Garant™ Dispenser

3M™ Elipar™ S10 Light Guide

3M™ Unitek™ Modular Band Kit Organizer







3M™ Unitek™ Elastics

3M™ Filtek™ Z250 Universal Restorative

3M[™] Sof-Lex[™] Contouring and Polishing Discs, 1981







3M™ Adper™ Prompt™ L-Pop™ Self-Etch Adhesive

3M™ Unitek™ Bonding Mixing Pad, 704-030, 5/Pk

3M™ Pentamix™ Lite Starter Kit — Impregum™ Penta™







3M™ Protemp™ II Temporization Material Refill

3M™ Elipar™ Eye Shield, 76971

Forsus™ Class II Correction System







3M™ AlastiK™ Bobbin Chain Dispensing System Replacement Blade, 406-92... 2/Pk

3M™ Light Wire Arch Computer, 701-720, 1/Ea

3M™ Ketac™ Fil Plus Glass Ionomer Filling Material, 55200 Series







3M™ RelyX™ Fiber Post Drill

3M™ Pentamix™ 3 Precision Starter Kit — Impregum™ Penta™

3M™ Unitek™ Archwire Hooks and Stops Wire, Weldable







3M™ Unitek™ Stainless Steel Straight Length Archwire

3M™ Unitek™ Gemini Clear Ceramic Brackets Kit

3M™ Unitek™ Gemini Brackets







3M™ Bonding Prep Tray w/Mixing Wells (Tacky Surface)

3M™ Iso-Form™ Temporary Crown Kit

3M™ Garant™ Intraoral Tips -Yellow (50), 71462







3M[™] Penta[™] Elastomer Syringe, 71210

3M™ RelyX™ Unicem 2 Automix Mix Tips Refill

3M™ Ketac™ Molar Easymix Glass Ionomer Filling Material, 56600 Series







3M™ Ketac™ Bond Aplicap™ Glass Ionomer Base Material Refill, 37350

3M™ Ketac™ Nano Light-Curing Glass Ionomer Restorative Kit, 3304

3M™ Pentamix™ 3 Starter Kit — Imprint™ 4 Preliminary Penta™







3M™ Ketac™ Nano Light-Curing Glass Ionomer Restorative, 3304

3M™ Pentamix™ 3 Precision Starter Kit — Imprint™ 4 Penta™ Putty

3M[™] Sof-Lex[™] Extra-Thin Contouring and Polishing Discs, 2381







3M™ RelyX™ Fiber Post Refill

3M™ Scotchbond™ Universal Etchant - Etching Gel, 41263, 2-3 ml

3M™ Pentamix™ Lite Starter Kit — Imprint™ 4 Preliminary Penta™







3M™ Sof-Lex™ Extra-Thin Contouring and Polishing Discs Refill, 2382

3M™ Ramitec™ Polyether Bite Registration Material Refill, 33000 Series

3M™ RelyX™ Ultimate Adhesive Resin Cement Kit







3M[™] Clarity[™] Advanced Ceramic Brackets

3M™ Stretch Block, Iso-Form™ Temporary Crown

3M™ Etching Liquid

3M™ Etching Liquid

Details

Acid etching of enamel and dentin Clear

Typical properties

Brands: 3M™

Delivery System: Bottle

Factory ISO Certification: No ISO Certification information available

Package Contents: Refill includes 9ml bottle eClass 14 Classification Group: 34110204

3M™ Stretch Block, Iso-Form™ Temporary Crown

Details

Assists in re-establishing the cervical contour of a crown if distorted by scissors while trimming

Bicuspid and molar selections for use with posterior crowns

Typical properties

Brands: Iso-Form™

Catalog Number: SBB45,SMB67

Factory ISO Certification: No ISO Certification information available

3M™ Clarity™ Advanced Ceramic Brackets

Details

Small, low profile bracket enhances patient comfort without compromising strength, while minimizing occlusal interference

Predictable and consistent debonding with 3M's proprietary stress concentrator on the base

Available precoated with APC™ Flash-Free Adhesive, for a fast and efficient bonding experience

Translucent ceramic bracket blends with tooth color, while resisting staining and discoloration during treatment

Smooth bracket surface and rounded corners minimize soft tissue contact and abrasion

Ample under tie-wing space, allowing for flexible treatment and ligation options Water-soluble horizontal and vertical reference markers assist with efficient bracket placement

Color-coded ID marks assist with accurate bracket identification

Patient can "blend in or be bold" by choosing from more than twenty clear or colorful 3M™ AlastiK™ Easy-to-Tie Ligatures

Extensive practice marketing support portfolio to support patient recruitment

Typical properties

Adhesive Precoat: APC™ Flash-Free, APC™ II, APC™ PLUS, Non-coated

Angulation: +11°, +2°, +3°, +4°, +5°, +7°, +8°, +9°, 0°

Arch: Lower, Upper Bracket Style: Twin Brands: Clarity™

Factory ISO Certification: No ISO Certification information available

Hook Included: No, Yes Hook Type: Distal

In/Out (Imperial): 0.02 in, 0.021 in, 0.026 in, 0.029 in, 0.03 in, 0.034 in, 0.035 in,

0.045 in

In/Out (Metric): 0.51 mm, 0.53 mm, 0.66 mm, 0.74 mm, 0.76 mm, 0.86 mm, 0.87

mm, 0.89 mm, 1.14 mm

Material: Ceramic

Prescription: High Torque Side: Left, Left/Right, Right

Slot Size (Imperial): 0.018 in, 0.022 in Slot Size (Metric): 0.457 mm, 0.559 mm

Tooth: Anterior, Bicuspid, Central Incisor, Cuspid, First Bicuspid, Lateral Incisor,

Second Bicuspid

Torque: +10°, +12°, +14°, +17°, +22°, +7°, +8°, -11°, -12°, -17°, -1°, -22°, -6°, -7°, 0°

Units per Pack: 5

3M™ RelyX™ Ultimate Adhesive Resin Cement Kit

Details

High bond strength

Easy-to-use, two-component system

Versatile for total-etch, selective etch or self-etch techniques

Excellent color stability with natural fluorescence

Offers wear resistance

If you like 3M™ RelyX™ Ultimate Cement, you'll love 3M™ RelyX™ Universal

Resin Cement

Typical properties

Brands: RelyX™

Catalog Number: 56892,56894

Cure Type: Dual Cure

Delivery System: Automix Syringe

Factory ISO Certification: No ISO Certification information available

Material: Adhesive Resin

Package Contents: 1 – 5ml Vial 3M[™] Scotchbond[™] Universal Adhesive; 1 – 3ml Etchant Syringe; 25 Dispensing Tips for Etchant Syringe; 50 Microbrushes; 1 – 8.5g Syringe cement Translucent; 10 Mixing Tips (regular); 1 Bag of 10 each: Mixing Tips (wide), Intra-oral Tips; Technique Guide; Instructions for Use,1 – 5ml Vial 3M[™] Scotchbond[™] Universal Adhesive; 1 – 3ml Etchant Syringe; 25 Dispensing Tips for Etchant Syringe; 50 Microbrushes; 1 – 8.5g Syringe Cement Translucent; 10 Mixing Tips (regular); 1 Bag of 5 each: Mixing Tips (wide), Intra-

oral Tips; Technique Guide; Instructions for Use

Shade: A1, TR

3M™ Ramitec™ Polyether Bite Registration Material Refill, 33000 Series

3M ID B00007635

Details

Hydrophilic polyether characteristics provide highly precise, detailed bite registrations

Direct application onto the occlusal for simple handling Thixotropic properties offer occlusal application without slumping Provides high final hardness combined with easy trimming

Typical properties

Brands: Ramitec™

Delivery System: Hand Mix

Impression Type: Bite Registration

Material: Polyether Set Type: Regular Set Setting Time: 5 Minute Work Time: 2 min

3M[™] Sof-Lex[™] Extra-Thin Contouring and Polishing Discs Refill, 2382

Details

Easy-to-use system designed to create high-luster restorations that offer lifelike esthetics

Discs are color coded to allow for easy identification of grit sequence Pop-on disc technology designed to save you time

Extra-thin (XT) discs allow for easy access to interproximal surfaces while offering a slightly stiffer disc

Typical properties

Brands: Sof-Lex™

Factory ISO Certification: No ISO Certification information available

Grade: Coarse, Fine, Medium, Super Fine

Model: XT

Package Contents: Refill includes 85 3M[™] Sof-Lex[™] Extra-Thin Contouring and Polishing Discs, 1/2" Diameter, Superfine grit,Refill includes 85 3M[™] Sof-Lex[™] Extra-Thin Contouring and Polishing Discs, 1/2" Diameter, Medium grit,Refill includes 85 3M[™] Sof-Lex[™] Extra-Thin Contouring and Polishing Discs, 1/2" Diameter, Fine grit,Refill includes 85 Sof-Lex[™] Extra-Thin Contouring and

Polishing Discs, 1/2" Diameter, Coarse grit

Product Color: Dark Orange, Light Orange, Medium Orange, Yellow

Product Type: Extra-thin Contouring & Polishing Disc

eClass 14 Classification Group: 34111407

Dimensions and Classifications

Overall Diameter (Imperial): 0.5 in Overall Diameter (Metric): 12.7 mm

3M™ Pentamix™ Lite Starter Kit — Imprint™ 4 Preliminary Penta™

3M ID B10397638

Typical properties

Brands: Pentamix™

3M™ Scotchbond™ Universal Etchant - Etching Gel, 41263, 2-3 ml

3M ID 7100048585 UPC 04035077010030

Details

Low-viscosity, gel-like consistency Gel has a distinct blue color to assist in complete rinsing away of the etchant Flexible dispensing tip can be individually adjusted before application Improved, easy-to-place etchant

Typical properties

Brands: Scotchbond™ Delivery System: Syringe

Factory ISO Certification: No ISO Certification information available

Package Contents: Refill includes: 2 - 3ml etchant syringes and 50 dispensing

tips

3M™ RelyX™ Fiber Post Refill

Details

Elasticity of RelyX Fiber Post is similar to human dentin for a lower risk of root fracture due to the "wedge effect"

Unique resin composition and equally dispersed parallel fibers lead to superior mechanical strength

Microporous surface allows for additional bond strength through microretentive effects

Excellent light translucency for efficient light cure deep into the root canal Improved radiopacity is among the highest of available glass fiber posts RelyX™ Unicem/Unicem 2 cements make cementing a post fast, easy and reliable

Typical properties

Apical Width (Imperial): 0.024 in, 0.028 in, 0.031 in, 0.035 in Apical Width (Metric): 0.6 mm, 0.7 mm, 0.8 mm, 0.9 mm

Brands: RelyX™

Catalog Number: 56863,56872,56861,56862

Coronal Width (Imperial): 0.043 in, 0.051 in, 0.063 in, 0.075 in Coronal Width (Metric): 1.1 mm, 1.3 mm, 1.6 mm, 1.9 mm

Factory ISO Certification: No ISO Certification information available

Material: Fiber Post

Package Contents: 10 – 3M[™] RelyX[™] Fiber Posts Size 2 (Red),10 – 3M[™] RelyX[™] Fiber Posts Size 0 (White),10 – 3M[™] RelyX[™] Fiber Posts Size 3

(Blue),10 – 3M[™] RelyX[™] Fiber Posts Size 1 (Yellow)

Product Color: Blue, Red, White, Yellow

Product Type: Fiber Post

Size: 0, 1, 2, 3

3M™ Sof-Lex™ Extra-Thin Contouring and Polishing Discs, 2381

Typical properties

Brands: Sof-Lex™

Factory ISO Certification: No ISO Certification information available

Grade: Coarse, Fine, Medium, Super Fine

Model: XT

Package Contents: Refill includes 85 3M™ Sof-Lex™ Extra-Thin Contouring and Polishing Discs, 3/8" Diameter, Fine grit, Refill includes 85 3M™ Sof-Lex™ Extra-Thin Contouring and Polishing Discs, 3/8" Diameter, Superfine grit, Refill includes 85 3M[™] Sof-Lex[™] Extra-Thin Contouring and Polishing Discs, 3/8" Diameter, Coarse grit, Refill includes 85 3M™ Sof-Lex™ Extra-Thin Contouring and

Polishing Discs, 3/8" Diameter, Medium grit

Product Color: Dark Orange, Light Orange, Medium Orange, Yellow

Product Type: Extra-thin Contouring & Polishing Disc

eClass 14 Classification Group: 34111407

Dimensions and Classifications

Overall Diameter (Imperial): 0.375 in Overall Diameter (Metric): 9.5 mm

3M™ Pentamix™ 3 Precision Starter Kit — Imprint™ 4 Penta™ Putty

3M ID B10397637

Typical properties

Brands: Pentamix™

3M™ Ketac™ Nano Light-Curing Glass Ionomer Restorative, 3304

Details

The first paste/paste, resin-modified glass ionomer (RMGI) restorative based on bonded nanofiller technology

Esthetic glass ionomer restorative

High fluoride release

Improved wear resistance

Light cure on demand

Radiopaque

Automatic mixing for fast dispensing from the unit-dose Quick Mix Capsule

Also available in our 3M™ Clicker™ Dispenser

Typical properties

Brands: Ketac™

Catalog Number: 3304A3.5,3304B2,3304A1,3304A2,3304A3

Cure Type: Light Cure

Factory ISO Certification: No ISO Certification information available

Package Contents: Refill includes 1 - 12g 3M[™] Clicker[™] Dispenser in shade B2,Refill includes 1 - 12g 3M[™] Clicker[™] Dispenser in shade A1,Refill includes 1

- 12g 3M™ Clicker™ Dispenser in shade A35,Refill includes 1 - 12g 3M™

Clicker™ Dispenser in shade A3,Refill includes 1 - 12g 3M™ Clicker™ Dispenser

in shade A2

Shade: A1, A2, A3, A3.5, B2 Total Units per Case: 10 per Case

3M™ Pentamix™ 3 Starter Kit — Imprint™ 4 Preliminary Penta™

3M ID B10397639

Typical properties

Brands: Pentamix™

3M™ Ketac™ Nano Light-Curing Glass Ionomer Restorative Kit, 3304

Details

Quick Mix Capsule for intraoral placement
The first paste/paste, resin-modified glass ionomer (RMGI) restorative based on
bonded nanofiller technology
Esthetic glass ionomer restorative
High fluoride release
Improved wear resistance
Light cure on demand
Radiopaque

Typical properties

Brands: Ketac™

Catalog Number: 3304I,3304TK

Cure Type: Light Cure

Factory ISO Certification: No ISO Certification information available

Package Contents: Trial kit includes 1 - 12g 3M™ Clicker™ Dispenser A2;

Primer (65ml); 10 Delivery Tips with pistons; Technique Guide, Intro kit includes 1

- 12g 3M[™] Clicker[™] Dispenser A1; 1 - 12g 3M[™] Clicker[™] Dispenser A3; Primer (65ml); 50 Delivery Tips with pistons; Mixing Pad; 50 Fiber Tips; 1

Applicator Handle; 48 Disposable Mixing Wells; Shade Guide; Technique Guide

Shade: A2, Multiple Shades

3M™ Ketac™ Bond Aplicap™ Glass Ionomer Base Material Refill, 37350

3M ID B00034543

Details

Adheres molecularly to enamel and dentine for marginal integrity
Offers high, long-term fluoride release
High radiopacity
Aplicap™ Capsule provides consistent quality and simplifies practical handling

Typical properties

Brands: Ketac™

Catalog Number: 37350 Delivery System: Capsule

3M™ Ketac™ Molar Easymix Glass Ionomer Filling Material, 56600 Series

3M ID B00034593

Details

High compressive strength

Minimal abrasion

Excellent handling characteristics: plastic contouring and condensing properties similar to amalgam

Aplicap capsule provides consistent quality, simplified handling and a nozzle for direct application of the material

Available in four shades

Typical properties

Brands: Ketac™

3M™ RelyX™ Unicem 2 Automix Mix Tips Refill

Details

For easy, automatic mixing of RelyX Unicem 2 Automix cement

Typical properties

Brands: RelyX™

Catalog Number: 56851, Refill - 56919, Refill - 56918

Factory ISO Certification: No ISO Certification information available

Package Contents: 15 Intra-oral Tips; 15 Mixing Tips (wide),30 Mixing Tips

(regular),15 Endo Tips; 15 Mixing Tips (wide)

Product Type: Tip

3M™ Penta™ Elastomer Syringe, 71210

3M ID 7100215236

Details

Allows for easy, intraoral delivery of 3M[™] Penta[™] impression materials Accessory for use with 3M[™] Pentamix[™] mixing uints Easy to use

Typical properties

Brands: Penta™

Catalog Number: 71210

Factory ISO Certification: No ISO Certification information available

Package Contents: 1 - Elastomer Syringe; 4 - Syringe Barrels; 5 - Intraoral Tips

Product Type: Syringe Set

3M™ Garant™ Intraoral Tips - Yellow (50), 71462

3M ID 7100215419

Details

For intraoral dispensing of dental materials using the cartridge dispenser 1:1/2:1

Typical properties

Brands: Garant™

Catalog Number: 71462

Factory ISO Certification: No ISO Certification information available

Package Contents: 50 - Intraoral Tips, yellow

Product Type: Tip

3M™ Iso-Form™ Temporary Crown Kit

Typical properties

Arch: Upper/Lower Brands: Iso-Form™ Catalog Number: MC64

Factory ISO Certification: No ISO Certification information available

Side: Left/Right Tooth: Molar

3M[™] Bonding Prep Tray w/Mixing Wells (Tacky Surface)

Details

7x7 upper and lower

Typical properties

Brands: 3M™

Factory ISO Certification: No ISO Certification information available

3M™ Unitek™ Gemini Brackets

Details

Cuspid and bicuspid brackets are available with patient-friendly integral hooks Color ID marking provides easy identification

Compatible in-out relationships with our 3M™ Victory Series™ Buccal Tubes means this system works for you

Typical properties

Adhesive Precoat: APC™ II, Non-coated

Angulation: +10°, +2°, +3°, +4°, +5°, +6°, +8°, 0°

Arch: Lower, Upper Bracket Style: Twin Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Hook Included: Yes Hook Type: Distal

Material: Stainless Steel

Package Contents: 10 Bkts/Pk,5 Bkts/Pk

Prescription: MBT™, Roth* Side: Left, Left/Right, Right

Slot Size (Imperial): 0.018 in, 0.022 in Slot Size (Metric): 0.457 mm, 0.559 mm

Tooth: Anterior, Bicuspid, Central Incisor, Cuspid, First Bicuspid, Lateral Incisor,

Second Bicuspid

Torque: +10°, +14°, +17°, +7°, -11°, -12°, -17°, -1°, -22°, -6°, -7°, 0°

Units per Pack: 5, 10

3M™ Unitek™ Gemini Clear Ceramic Brackets Kit

Typical properties

Adhesive Precoat: Non-coated

Arch: Upper/Lower Bracket Style: Twin Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Hook Included: Yes Hook Type: Distal Material: Ceramic

Package Contents: 20 Bkts/Kit,12 Bkts/Kit

Prescription: MBT™

Slot Size (Imperial): 0.022 in Slot Size (Metric): 0.559 mm

Units per Pack: 12, 20

3M™ Unitek™ Stainless Steel Straight Length Archwire

Details

Round Wire - 14 in. (35.6 cm) straight lengths

Rectangular and square wires have a slightly radiused corners facilitate insertion into bracket slot with no appreciable loss of torque

Square wire for control with optimum sliding mechanics

Tighter diagonal tolerance on .016 x .016 to avoid rolling in slot

Typical properties

Arch Shape: Standard

Archwire Size (Imperial): 0.01 in, 0.012 in, 0.014 in, 0.015 in, 0.016 in, 0.016 x 0.016 in, 0.016 x 0.022 in, 0.017 x 0.022 in, 0.017 x 0.025 in, 0.018 in, 0.018 in, 0.018 x 0.018 in, 0.018 x 0.018 in, 0

Archwire Size (Metric): 0.254 mm, 0.305 mm, 0.356 mm, 0.381 mm, 0.406 mm, 0.406 x 0.406 mm, 0.406 x 0.559 mm, 0.432 x 0.559 mm, 0.432 x 0.635 mm, 0.457 mm, 0.457 x 0.457 mm, 0.457 x 0.635 mm, 0.457 x 0.635 mm, 0.508 mm, 0.508 x 0.635 mm, 0.533 x 0.635 mm, 0.533 x 0.635 mm, 0.559 x 0.711 mm, 0.61 mm, 0.635 mm, 0.666 mm, 0.711 mm, 0.762 mm, 0.813 mm, 0.914 mm, 1.016 mm, 1.143 mm, 1.295 mm, 1.397 mm

Brands: Unitek™

Cross Section Shape: Rectangular, Round, Square

Factory ISO Certification: No ISO Certification information available

Material: Stainless Steel

Units per Inner Pack: 10 per Pack, 15 per Pack, 20 per Pack, 25 per Pack, 30 per

Pack, 35 per Pack, 40 per Pack

3M™ Unitek™ Archwire Hooks and Stops Wire, Weldable

Details

Medium hard temper
One or two welds to permanently anchor
.006 in. (.15 mm) wall
Choice of two lengths: Round 2 mm or 1 mm; Rectangular 2 mm

Typical properties

Brands: Unitek™

Compatible Archwire Size (Imperial): 0.018 - 0.028 in, 0.02 in, 0.022 in, 0.03 in,

0.04 in, 0.045 in

Factory ISO Certification: No ISO Certification information available

Shape: Rectangular, Round

3M™ Pentamix™ 3 Precision Starter Kit — Impregum™ Penta™

3M ID B10397636

Typical properties

Brands: Pentamix™

3M™ RelyX™ Fiber Post Drill

Details

Available in four sizes plus universal and are color-coded to ensure an accurate match and prevent errors

Unique resin composition and equally dispersed parallel fibers lead to superior mechanical strength

Microporous surface allows for additional bond strength through microretentive effects

Excellent light translucency for efficient light cure deep into the root canal Improved radiopacity is among the highest of available glass fiber posts RelyX™ Unicem system makes cementing a post fast, easy and reliable

Typical properties

Brands: RelyX™

Catalog Number: 56864,56865,56866,56873,56867

Color Code: Blue, Red, Yellow

Factory ISO Certification: No ISO Certification information available

Package Contents: 1 Drill – Reamer/Universal,1 Drill – Size 2 (1.6mm dia.) Red,1 Drill – Size 3 (1.9mm dia.) Blue,1 Drill – Size 0 (1.1mm dia.) White,1 Drill – Size 1

(1.3mm dia.) Yellow

Product Type: Fiber Post Drill

3M™ Ketac™ Fil Plus Glass Ionomer Filling Material, 55200 Series

3M ID B00034561

Details

High sustained fluoride release Radiopaque Self-cure One-step placement

Typical properties

Brands: Ketac™

Delivery System: Hand Mix

3M™ Light Wire Arch Computer, 701-720, 1/Ea

3M ID 70202088848

Details

Used for fast determination of proper archwire size

Typical properties

Brands: 3M™

Factory ISO Certification: No ISO Certification information available

Units per Inner Pack: 1 per Pack

3M[™] AlastiK[™] Bobbin Chain Dispensing System Replacement Blade, 406-928, 2/Pk

3M ID 7100287023

Details

Not made with natural rubber latex Contains two replacement blades

Typical properties

Brands: AlastiK™

Factory ISO Certification: No ISO Certification information available

Units per Inner Pack: 2 per Pack

Forsus™ Class II Correction System

3M ID B5005051009

Details

A simple one-step installation into the headgear tube
A telescoping coaxial spring designed to allow for enhanced patient jaw opening
A comfortable stainless steel mechanism that does not bow out into the cheek
Low spring stresses to give the Forsus spring a longer, more fatigue-resistant life
Open spring design to facilitate good oral hygiene
One appointment installation, chairside, with no lab work required
Independent left and right components; install only one side if desired

Typical properties

Brands: Forsus™

3M™ Elipar™ Eye Shield, 76971

Part Number 76971 3M Product Number 76971 3M ID 7100003375

Details

For use with the 3M™ Elipar™ LED Curing Light

Typical properties

Brands: Elipar™

Catalog Number: 76971

Factory ISO Certification: No ISO Certification information available

Product Type: Eye Shield

3M™ Protemp™ II Temporization Material Refill

3M ID B00006517

Details

Easy-to-extrude delivery system
High biocompatibility
Excellent esthetics due to great color stability, high polishability and the choice of three intensive shades: A1 extra light, A3 light and B3 yellow. Shades may be mixed to allow better shade matching

Typical properties

Brands: Protemp™

3M™ Pentamix™ Lite Starter Kit — Impregum™ Penta™

3M ID B00007779

Typical properties

Brands: Pentamix™

3M™ Unitek™ Bonding Mixing Pad, 704-030, 5/Pk

3M ID 7100154482

Details

Liquid resistant surface 3" X 2 1/2" 5/pack

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

eClass 14 Classification Group: 34110110

Dimensions and Classifications

Overall Length (Imperial): 3 in Overall Length (Metric): 76.2 mm Overall Width (Imperial): 2.5 in Overall Width (Metric): 63.5 mm

3M™ Adper™ Prompt™ L-Pop™ Self-Etch Adhesive

Details

Bonds to dentin, cut and uncut enamel without a separate etch step

Clinically proven

Very low incidence of post-operative sensitivity

Desensitizes hypersensitive root surfaces

Simplifies sealant technique

Cures with any light source - halogen, laser, plasma or LED

Ensures hygienic delivery of the adhesive material

Extra small applicator works well with very small preperations and in hard-to-

reach areas of the mouth

Typical properties

Adhesive Type: Self-etch

Brands: Adper™

Catalog Number: 41928,41926

Cure Type: Light Cure Delivery System: L-Pop™

Factory ISO Certification: No ISO Certification information available Package Contents: 500 L-Pop™ applicators,100 L-Pop™ applicators

Product Type: Dental Adhesive

Units per Inner Pack: 100 per Pack, 500 per Pack

3M™ Sof-Lex™ Contouring and Polishing Discs, 1981

Typical properties

Brands: Sof-Lex™

Factory ISO Certification: No ISO Certification information available

Grade: Coarse, Fine, Medium, Super Fine

Model: Standard

Package Contents: Refill includes 85 3M[™] Sof-Lex[™] Contouring and Polishing Discs, 3/8" Diameter, Superfine grit,Refill includes 85 3M[™] Sof-Lex[™] Contouring and Polishing Discs, 3/8" Diameter, Medium grit,Refill includes 85 3M[™] Sof-Lex[™] Contouring and Polishing Discs, 3/8" Diameter, Coarse grit,Refill includes 85 3M[™] Sof-Lex[™] Contouring and Polishing Discs, 3/8" Diameter, Fine grit

Product Color: Black, Dark Blue, Light Blue, Medium Blue

Product Type: Contouring & Polishing Disc eClass 14 Classification Group: 34111407

Dimensions and Classifications

Overall Diameter (Imperial): 0.375 in Overall Diameter (Metric): 9.5 mm

3M™ Filtek™ Z250 Universal Restorative

Details

Versatile use for both anterior and posterior restorations

Reduced shrinkage for less post-operative sensitivity, low wear and high fracture toughness

Good handling with resistance to slumping

Polishes easily for long-lasting shine

Available in 12 shades for natural esthetics

Complete your premium restoration with the gold standard of adhesion, 3M[™] Scotchbond[™] Universal Plus Adhesive

Typical properties

Brands: Filtek™ Catalog Number:

1370A3.5,6020B1,6020D3,6021B2,1370A3,6020B3,6021B1,6021D3,1370A2,60 20B2,6021B3,1370A4,1370A1,1370UD,370A3.5,6020I,6021A3.5,6020UD,6020C 2,6021A1,6020A2,6021C2,6021UD,6020A1,6021A3,6020A4,6021A2,6020A3,60 21A4,6020A3.5,1370B2,1370D3,1370B1,6021I

Composite Type: Universal Cure Type: Light Cure

Delivery System: Capsule, Syringe Factory ISO Certification: ISO 13485 Inner Packs per Case: 10 Packs per Case Material Structure: Microhybrid Composite

Net Weight (Metric): 0.2 g, 4 g

Opacity: Body, Dentin

Package Contents: Refill includes 20 - 02g capsules in shade Universal Dentin,Refill includes 1 - 4g syringe in shade A4,Refill includes 1 - 4g syringe in shade Universal Dentin,Refill includes 1 - 4g syringe in shade A35,Refill includes 20 - 02g capsules in shade A4,Refill includes 20 - 02g capsules in shade B2,Refill includes 20 - 02g capsules in shade A1,Refill includes 20 - 02g capsules in shade C2,Refill includes 20 - 02g capsules in shade B1,Refill includes 20 - 02g capsules in shade D3,Refill includes 20 - 02g capsules in shade A3,Refill includes 20 - 02g capsules in shade A2,Refill includes 20 - 02g capsules in shade Incisal,Refill includes 20 - 02g capsules in shade B3,Refill includes 1 - 4g syringe in shade B3,Refill includes 1 - 4g syringe in shade A3,Refill includes 1 - 4g syringe in shade C2,Refill includes 1 - 4g syringe in shade B1,Refill includes 1 - 4g syringe in shade C2,Refill includes 1 - 4g syringe in shade Incisal,Refill includes 1 - 4g syringe in shade D3,Refill includes 1 - 4g syringe in shade D3,Refill includes 1 - 4g syringe in shade A1

Product Form: Refill

Shade: A1, A2, A3, A3.5, A4, B1, B2, B3, C2, D3, I, UD

Total Units per Case: 10 per Case

Units per Inner Pack: 1 per Pack, 20 per Pack

Units per Pack: 1, 20

3M™ Unitek™ Elastics

Details

Made with natural rubber latex Coated with cornstarch Precision-cut, surgical-grade latex Consistent range and force

No color additives

Time and cost effective – no need to count, separate, and repackage

Easy, positive identification

Size and force rating are easily identified

Cross indexed by color and symbol

Easy to dispense - prepackaged in quantities of 100/pk

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Force (Imperial): 3.5 ozf,2.0 ozf,4.0 ozf,6.0 ozf

Force (Metric): 99.223 gf,56.699 gf,170.097 gf,113.398 gf

Force Rating: Heavy, Light, Medium, Medium-light

Inner Diameter (Imperial): 0.125 in, 0.187 in, 0.188 in, 0.25 in, 0.312 in, 0.313 in,

0.375 in, 0.625 in, 0.75 in

Inner Diameter (Metric): 3.175 mm, 3.2 mm, 4.6 mm, 4.762 mm, 6.35 mm, 6.4

mm, 7.9 mm, 7.937 mm, 9.525 mm, 15.875 mm, 19.05 mm

Packaging Color: Agua, Blue, Green, Orange, Purple, Red, Yellow

3M™ Unitek™ Modular Band Kit Organizer

Details

Modular

Durable plastic containers enclose sliding kit drawers with clear plastic hinged covers to protect your inventory

Color coded labels identify band types

Typical properties

Arch: Lower, Upper, Upper/Lower

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Product Type: Band Side: Left, Right

Tooth: Bicuspid, First Molar, First/Second Molar, Molar, Pedodontic Molar,

Second Molar

3M™ Elipar™ S10 Light Guide

Details

Includes one 3mm light guide

Typical properties

Brands: Elipar™

Catalog Number: 76958

Factory ISO Certification: No ISO Certification information available

Package Contents: Includes one 3mm light guide, Includes one 10mm light

guide, Includes one 8mm light guide

Product Type: Light Guide

3M™ Garant™ Dispenser

Details

Consistent material characteristics with void-free mixing

Typical properties

Brands: Garant™

Catalog Number: 77580,77581

Factory ISO Certification: No ISO Certification information available

Package Contents: 1 - Cartridge Dispenser 4:1/10:1,1 - Cartridge Dispenser

1:1/2:1 Dispenser w/plunger Product Type: Dispenser

3M™ Sof-Lex™ Contouring and Polishing Discs, 1982

Details

Easy-to-use system allows you to create high-luster restorations that offer lifelike esthetics

Discs are color coded to allow for easy identification of grit sequence Pop-on disc technology makes changing discs easier to save you time

Typical properties

Brands: Sof-Lex™

Factory ISO Certification: No ISO Certification information available

Grade: Coarse, Fine, Medium, Super Fine

Model: Standard

Package Contents: Refill includes 85 3M[™] Sof-Lex[™] Contouring and Polishing Discs, 1/2" Diameter, Coarse grit,Refill includes 85 3M[™] Sof-Lex[™] Contouring and Polishing Discs, 1/2" Diameter, Medium grit,Refill includes 85 3M[™] Sof-Lex[™] Contouring and Polishing Discs, 1/2" Diameter, Fine grit,Refill includes 85 3M[™] Sof-Lex[™] Contouring and Polishing Discs, 1/2" Diameter, Superfine grit

Product Color: Black, Dark Blue, Light Blue, Medium Blue

Product Type: Contouring & Polishing Disc eClass 14 Classification Group: 34111407

Dimensions and Classifications

Overall Diameter (Imperial): 0.5 in Overall Diameter (Metric): 12.7 mm

3M™ Pentamix™ Lite Starter Kit — Imprint™ 4 Penta™ Putty

3M ID B10397653

Typical properties

Brands: Pentamix™

3M™ Unitek™ Prestige Tweed Loop Forming Pliers, 900-720, 1/Ea

3M ID 7000004393

Details

Efficient handling, bending, and torquing of square and rectangular wires 3M™ Unitek™ Prestige Series

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

3M™ Clarity™ Ultra Self-Ligating Brackets

Details

Translucent and stain-resistant, Clarity Ultra brackets are an excellent fullyesthetic choice for you and your patients

Innovative, smooth ceramic bracket surfaces and edges are designed for patient comfort so your patients can live confidently

Clarity™ Ultra Self-Ligating Brackets are passive, but can be activated with the proprietary oval under tie-wing ligation design

A wide ceramic door assists rotational control with a broad mesial-distal span Reliable door mechanism designed to resist unplanned opening and closing Unique stress concentrator for predictable squeeze debonding

Clarity Ultra brackets do not require elastomeric ligatures which may stain during treatment

Color-coded with identification dots, vertical and horizontal reference marks for easy bracket selection, and positioning

Compatible with many 3M appliances to give orthodontists custom control and achieve results based on treatment preferences and the unique needs of the patient

Typical properties

Adhesive Precoat: APC™ Flash-Free, Non-coated

Angulation: +2°, +3°, +4°, +8°, 0°

Arch: Lower, Upper Bracket Style: Twin Brands: Clarity™

Factory ISO Certification: No ISO Certification information available

Hook Included: No, Yes Hook Type: Distal

In/Out (Imperial): 0.02 in, 0.021 in, 0.026 in, 0.029 in, 0.03 in, 0.034 in, 0.035 in,

0.045 in

In/Out (Metric): 0.51 mm, 0.53 mm, 0.66 mm, 0.74 mm, 0.76 mm, 0.86 mm, 0.89

mm, 1.14 mm Material: Ceramic

Side: Left, Left/Right, Right Slot Size (Imperial): 0.022 in Slot Size (Metric): 0.559 mm

Tooth: Anterior, Bicuspid, Central Incisor, Cuspid, First Bicuspid, Lateral Incisor,

Second Bicuspid

Torque: +10°, +17°, -12°, -17°, -6°, -7°, 0°

Units per Pack: 5

3M™ ESPE™ Stainless Steel Primary Molar Crowns, ND-96, Kit

3M ID B5005089040 UPC 30605861011236

Typical properties

Arch: Upper/Lower Brands: 3M™

Catalog Number: ND-96

Dentition: Primary

Factory ISO Certification: No ISO Certification information available

Number of Crowns: 96 Patient: Pediatric Size: Assorted

Tooth: Assorted Molars, First Molar

3M[™] OrthoForm[™] Diagnosis Arch Form Templates, Clear, 701-723, 1/Ea

3M ID 7000025648

Details

Designed to aid in arch form symmetry determination

Typical properties

Brands: OrthoForm™

Factory ISO Certification: No ISO Certification information available

Units per Inner Pack: 1 per Pack

3M™ RelyX™ Luting Cement Intro Kit, 3505

3M ID 7000128788 UPC 30605861012554

Details

Easy clean-up
Sustained fluoride release
Virtually no post-operative sensitivity
If you like 3M™ RelyX™ Luting Cement, you'll love 3M™ RelyX™ Universal
Resin Cement

Typical properties

Brands: RelyX™

Catalog Number: 3505 Cure Type: Self-cure Delivery System: Hand Mix

Factory ISO Certification: No ISO Certification information available

Inner Packs per Case: 10 Packs per Case Material: Resin-modified Glass Ionomer

Package Contents: 1 – 16g Powder, Spoon; 1 – 9ml Liquid; Mixing Pads (1961P)

Shade: White Opaque

Total Units per Case: 10 per Case Units per Inner Pack: 1 per Pack

3M™ RelyX™ Unicem Self-Adhesive Universal Resin Cement Maxicap™ Capsule Intro Kit

3M ID B00007441

Details

Saves time by eliminating etching, priming and bonding steps
Virtually no post-operative sensitivity
Less discoloration
Strong and moisture-tolerant
If you like 3M™ RelyX™ Unicem Self-Adhesive Resin Cement, you'll love 3M™
RelyX™ Universal Resin Cement

Typical properties

Brands: RelyX™

Catalog Number: 56830 Cure Type: Dual Cure

Delivery System: Maxicap™ Capsule

3M™ Transbond™ XT Light Cure Adhesive

Details

Viscosity designed to prevent adhesive run-on and bracket skating, saving you money while reducing adhesive waste

Versatile — offers efficient bonding of ceramic and metal brackets

Excellent adhesive handling properties — easy control of bracket positioning and easy excess clean-up prior to curing

Immediate bond strength allows for archwire tie-in right after light cure Available in syringes and capsules

Typical properties

Adhesive Type: Direct Brands: Transbond™

Catalog Number: 712-031,712-030,712-035,712-036

Cure Type: Light Cure

Delivery System: Capsule, Syringe

Factory ISO Certification: No ISO Certification information available Package Contents: 4 - 4g syringes,2 - 4g syringes,25 - 0.2g capsules Product Type: Light Cure Band Adhesive Kit, Orthodontic Adhesive

Units per Inner Pack: 1 /KT, 4 per Pack, 25 per Pack

3M™ Permadyne™ 2:1 Syringe Impression Material Refill Pack

Details

Hand-mix material for one-step heavy/light impression techniques

Polyether properties

Hydrophilicity for moisture control and accuracy

Thixotropic behavior and good flow properties throughout the entire working time

Typical properties

Brands: Permadyne™

Catalog Number: 30312,30103,30751,30413,30711

Delivery System: Garant™ Dispenser, Garant™ Hand Dispenser, Hand Mix,

Pentamix™ Mixing Unit

Factory ISO Certification: No ISO Certification information available

Impression Type: Final Material: Polyether

Package Contents: 1 - base paste, 300ml; 1 - catalyst 60ml,2 - base paste, 120ml each; 2 - catalyst, 15ml each; 1 - mixing pad,4 - 50ml cartridges; 20 - 3M™ Garant™ Mixing Tips - White,2 - base paste, 120 ml each; 2 - tubes catalyst - 15ml each; 1 - mixing pad,2 - base paste, 300ml each; 2 - catalyst,

60ml each

Set Time: 5.5 min, 6 Minutes

Set Type: Regular Set

Setting Time: 5.5 Minute, 6 Minute

Viscosity: Heavy, Light

Work Time: 2 min, 2 min 30 sec, 3 min 30 sec eClass 14 Classification Group: 34110100

3M™ Unitek™ Convertible Cap Removing Pliers Tip 805-036 5/Pack

3M ID 7000004388

Details

Replacement tips are five per pack

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

3M™ Sinfony™ Material Shade Dispenser Refills

3M ID B00007998

Details

Excellent esthetics/translucency Long-term proven wear resistance Ease of finishing and polishing Plaque and stain resistance

Proven in tests to be as strong or stronger than other popular brands and to possess inherent advantages in flexibility, impact strength and color stability Handles like porcelain

Indirect restorations created from Sinfony indirect lab composite offer natural vitality, amber opalescent effect and fluorescence

Sinfony Magic Material is available in 11 shades: I1 Ocean Blue, I2 Atlantis, I3 Chestnut, I4 Havanna, I5 Orange, I6 Khaki, I7 Vanilla, I8 Honey Yellow, I9 Gingiva, I1 Red, I11 Terracotta.itality, amber opalescent effect and fluorescence

Typical properties

Brands: Sinfony™

3M™ Elipar™ I Ellipoid Lamp, 77160

3M ID 7000054660

Details

Allows for direct and functional layering of light-curedcomposite material Built-in proximity sensor automatically controls the light source once the object is positioned at the correct distance from the lamp

Typical properties

Brands: Elipar™

Factory ISO Certification: No ISO Certification information available

Product Type: Replacement Light

3M™ Vitremer™ Core Buildup/Restorative Kit, 3303

Details

Versatile system can be light-cured or self-cured to meet your application needs Unique dark-cure chemistry allows for bulk placement which eliminates the need for layering and saves time

Displays a true glass ionomer cure for long-term fluoride release and adhesion to tooth structure which provides better patient protection and increased doctor confidence

As a restorative material, it provides good handling properties and extended working time in the mouth for easier use

Typical properties

Brands: Vitremer™

Catalog Number: 3303I,3303MP-A3,3303MPP,3303MPB

Cure Type: Light Cure Delivery System: Hand Mix

Factory ISO Certification: No ISO Certification information available

Package Contents: Trial kit includes 5g bottle of powder in Pedo shade, 25ml liquid, 2ml primer, 2ml finishing gloss, and accessories, Trial kit includes 5g bottle of powder in shade A3, 25ml liquid, 2ml primer, 2ml finishing gloss, delivery tips with pistons (10 each), spoon and mixing pad, Trial kit includes 5g bottle of powder in Blue shade, 25ml liquid, 2ml primer, delivery tips with pistons (10 each), spoon, and mixing pad, Intro kit includes 5g powder each of A3, A4, C2, and C4 shades, 2 - 8ml liquid, 65ml primer, 65ml finishing gloss, shade guide and accessories

Shade: A3, Blue, Multiple Shades, PEDO

Total Units per Case: 10 per Case

3M™ Ketac™ Molar Aplicap™ Glass Ionomer Restorative Refill, 56400 Series

Details

Excellent mechanical properties for posterior applications
Bulk fill convenience
High fluoride release
Contouring and condensing properties similar to amalgam
3M™ Aplicap™ capsule provides consistent quality, simplified handling and a nozzle for direct application of the material

Typical properties

Brands: Ketac™

Catalog Number: 56430,56420,56410 Delivery System: Aplicap™ Capsule

Factory ISO Certification: No ISO Certification information available

Package Contents: Refill includes 50 capsules in shade A1,Refill includes 50 capsules in shade A3,Refill includes 50 capsules in assorted shades (10 each of

A1, A4, DYO; 20 of A3)

Shade: A1, A3, Multiple Shades

3M[™] Adper[™] Scotchbond[™] Multi-Purpose Adhesive System Intro Kit, 7540S

3M ID B00007555 UPC 30605861066908

Details

Improve patient outcomes with a total-etch system that offers excellent bond strength to a variety of dental surfaces

Creates a complete, versatile dental bonding system when used with 3M[™] Adper[™] Scotchbond[™] Multi-Purpose Primer

Obtain strong and reliable bonds both to wet and dry surfaces, even when complete isolation from moisture is difficult

Reduces microleakage that can contribute to post-operative sensitivity, helping to ensure a positive patient experience

Typical properties

Adhesive Type: Total Etch

Brands: Adper™

Catalog Number: 7540S Cure Type: Light Cure Delivery System: Vial

Factory ISO Certification: No ISO Certification information available

Product Type: Adhesive System Total Units per Case: 10 per Case

3M™ Lava™ Plus High Translucency Zirconia Effect Shade

Details

Available in 8 shades Chrome and cobalt-free

Typical properties

Brands: Lava™

CAD/CAM System: Lava

Factory ISO Certification: No ISO Certification information available

Material: Zirconia

Shade: Brown, Fluorescence, Gray, Orange, Pink, Purple, White, Yellow eClass 14 Classification Group: 22320201

3M™ Unitek™ Gemini Brackets Kit

Details

Cuspid and bicuspid brackets are available with patient-friendly integral hooks Color ID marking provides easy identification

Compatible in-out relationships with our 3M[™] Victory Series[™] Buccal Tubes means this system works for you

Typical properties

Adhesive Precoat: Non-coated

Arch: Upper/Lower Bracket Style: Twin Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Hook Included: Yes Hook Type: Distal Material: Stainless Steel

Package Contents: 20 Bkts/Kit Prescription: MBT™, Roth*

Slot Size (Imperial): 0.018 in, 0.022 in Slot Size (Metric): 0.457 mm, 0.559 mm

Torque: -11°, -6° Units per Pack: 20

3M™ Visio™ Vacuum Pump Oil Refill, 77150

3M ID 7100005840 UPC 04035077000123

Details

Ideal for polymerization of veneers made from 3M[™] Sinfony[™] Indirect Lab Composite to guarantee veneer strength and color stability

Mirrored polymerization chamber provides complete illumination for deep and thorough polymerization

Vacuum technique removes oxygen-inhibited layers for improved surface hardness

Includes three additional programs for versatile use, e.g., for curing MultiTray Light-Cured Maxillary and Mandibular Baseplates

Typical properties

Brands: Visio™

Factory ISO Certification: No ISO Certification information available

Product Type: Oil Refill

3M™ Indirect Bonding Technique DVD

3M ID B5005126014

Details

DVD

Typical properties

Brands: 3M™

Catalog Number: 012-147

Units per Inner Pack: 1 per Pack eClass 14 Classification Group: 34110205

3M™ Unitek™ Prestige Bond Removing Pliers, 900-717, 1/Ea

3M ID 7000025672

Details

Replacement blade 3M™ Unitek™ Prestige Series

Typical properties

Brands: 3M™

Factory ISO Certification: No ISO Certification information available

3M™ Light-Tight Box, 712-075, 1/Ea

3M ID 70202097799

Details

Holds 2 arches with light sensitive adhesive Plastic box with latch Black

Typical properties

Brands: 3M™

Catalog Number: 712-075

Factory ISO Certification: No ISO Certification information available

Units per Inner Pack: 1 per Pack

3M™ Forsus™ Push Rod Measurement Guide, 807-014, 5/Pk

3M ID 7000026037

Details

For 5-patient and 20-patient kits

Typical properties

Brands: Forsus™

Factory ISO Certification: No ISO Certification information available

Units per Inner Pack: 5 per Pack

3M™ Unitek™ Archwire Hooks and Stops Wire Crimpable

Details

Light loads 2 mm length Wall: round .015 in. (.38 mm) Ideal with Nitinol wire

Typical properties

Brands: Unitek™

Compatible Archwire Size (Imperial): 0.016 - 0.026 in, 0.018 - 0.025 in, 0.022 -

0.028 in, 0.022 in, 0.025 in, 0.036 in

Factory ISO Certification: No ISO Certification information available

Shape: Rectangular, Round

3M™ Unitek™ Lateral Development Archwire Arch Evaluator, 701-726, 5/Pk

3M ID 7100008869

Details

Helps estimate the necessary lateral development
Allows you to identify upper cuspid "Sweet Spot" position on wires and confirm
proper upper cuspid bracket positioning
Helps evaluate progression of lateral development
Use on diagnostic model as an early guide for archwire sizing
Can be used with wax bite to assess lateral development
Confirms arch form size prior to insertion
Tabletop guide for restoring shape following adjustment bends
Reference for customizing extra-large or extra-small arch forms in extreme cases

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Units per Inner Pack: 5 per Pack

3M™ Forsus™ EZ2 Spring Modules

Typical properties

Brands: Forsus™

Factory ISO Certification: No ISO Certification information available

Module Style: EZ2

Product Type: Fatigue Resistant Device

Side: Left, Right

Units per Inner Pack: 5 per Pack eClass 14 Classification Group: 34110690

3M[™] Unitek[™] Lingual Retainer Wire, 260-032, 7 in., .032, Bondable, 25/Pk

Part Number 260-032 3M ID 7000003978

Typical properties

Archwire Size (Imperial): 0.032 in Archwire Size (Metric): 0.813 mm

Brands: Unitek[™]

Factory ISO Certification: No ISO Certification information available

Units per Inner Pack: 25 per Pack

3M[™] Unitek[™] Prestige Bond Removing Pliers, Replacement Blade, 900-718, 1/Ea

3M ID 7000025673

Details

Replacement pads 3M™ Unitek™ Prestige Series

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

3M™ ESPE™ Sof-Lex™ Finishing Strips Refill, Fine/Superfine, 1956

3M ID 7100066455 UPC 30605861059054

Details

Designed with ease for use and to smooth proximal surfaces
Strips are gapped at their centers for easy interproximal insertion
Precision-coated with an aluminum-oxide grit to produce a smooth restoration
surface that helps prevent interproximal plaque buildup
Unique narrow strip is available for easier use in short interproximal areas

Typical properties

Brands: Sof-Lex™

Catalog Number: 1954,1956

Factory ISO Certification: No ISO Certification information available

Grade: Coarse/Medium, Fine/Super Fine

Package Contents: Refill includes 120 fine/superfine strips (7" x 5/32")

eClass 14 Classification Group: 34110000

Dimensions and Classifications

Overall Length (Imperial): 7 in Overall Length (Metric): 177.8 mm Overall Width (Imperial): 0.156 in Overall Width (Metric): 3.962 mm

3M™ Lava™ Sintering Beads, 68594, 100g

Part Number 68594 3M Product Number 68594 3M ID 7000145345

Details

Includes 3M™ Lava™ Sintering Beads

Typical properties

Brands: Lava™

CAD/CAM System: Lava Catalog Number: 68594

Factory ISO Certification: No ISO Certification information available

Product Type: Sintering Bead

3M™ Sinfony™ Indirect Lab Composite Master System, 49962

3M ID 70201112896

Details

Excellent esthetics/translucency Long-term proven wear resistance Ease of finishing and polishing Plaque and stain resistance

Proven in tests to be as strong or stronger than other popular brands and to possess inherent advantages in flexibility, impact strength and color stability Handles like porcelain

Indirect restorations created from Sinfony indirect lab composite offer natural vitality, amber opalescent effect and fluorescence

Typical properties

Factory ISO Certification: No ISO Certification information available

3M™ Pentamix™ 3 Starter Kit — Impregum™ Polyether

3M ID B10397648

Typical properties

Brands: Pentamix™

Ketac™ Molar Quick Aplicap™ Glass Ionomer Restorative

Details

Fast mixing and setting time

Excellent mechanical properties for posterior applications

Bulk fill convenience

High fluoride release

Available in five shades

3M[™] Aplicap[™] capsule provides consistent quality, simplified handling and a nozzle for direct application of the material

Aplicap capsule provides consistent quality, simplified handling and a nozzle for direct application of the material

Excellent handling characteristics: plastic contouring and condensing properties similar to amalgam

Minimal abrasion

High compressive strength

Typical properties

Brands: Ketac™

Catalog Number: 56369,56364,56366,56363

Delivery System: Aplicap™ Capsule

Factory ISO Certification: No ISO Certification information available

Package Contents: Refill includes 50 capsules in assorted shades, Refill includes 50 capsules in shade A1, Refill includes 50 capsules in shade A2, Refill includes

50 capsules in shade A3

Shade: A1, A2, A3, Multiple Shades

3M™ Grommets for Eye Shield

3M ID 7000057034

Details

Replacement part for eye shield. 10 pieces

Typical properties

Brands: 3M™

Factory ISO Certification: No ISO Certification information available

Product Type: Eye Shield

3M™ Dimension™ VPS Impression Materials

Details

Accurate impression results
Void-free mix through the 3M™ Pentamix™ Automatic Mixing Unit
Good detail reproduction
Hydrophilic properties

Typical properties

Brands: Dimension™

Catalog Number: 36913,36715,36954,36515

Delivery System: Garant™ Dispenser, Pentamix™ Mixing Unit Factory ISO Certification: No ISO Certification information available

Impression Type: Final

Material: VPS

Package Contents: 2 - base paste, 300 ml each; 2 - catalyst, 60 ml each, 4 - 50ml

cartridges; 20 3M™ Garant™ Mixing Tips - Yellow

Product Color: Green, Light Green, Pink Set Time: 4 min, 5.5 min, 6 Minutes Set Type: Quick Set, Regular Set

Setting Time: 4 Minute, 5.5 Minute, 6 Minute

Viscosity: Heavy, Light

Work Time: 1 Minute 30 Seconds, 2 min eClass 14 Classification Group: 34110100

3M™ Protemp™ Crown Temporization Material

Details

Nine preformed sizes with malleable walls for a custom fit Natural-looking esthetics Low intra-oral heat during fitting process

Typical properties

Arch: Lower, Upper, Upper/Lower

Brands: Protemp™

Catalog Number: 50617,50616,50618,50611,50610,50613,50612,50615,50614

Dentition: Permanent

Factory ISO Certification: No ISO Certification information available

Material: Composite Number of Crowns: 5 Pre-formed: Yes Side: Left/Right

Tooth: Bicuspid, Cuspid, Molar

Units per Case: 5.0

3M™ Gold Anodized Crowns Kit, 942001 Series

Details

Made from medium-hard aluminum for good durability Molar crown selection to cover adult posterior applications

Typical properties

Arch: Upper/Lower Brands: 3M™

Catalog Number: 942301,942501

Dentition: Permanent

Factory ISO Certification: No ISO Certification information available

Material: Aluminum

Number of Crowns: 72, 80

Side: Left/Right

Tooth: Assorted Bicuspids, Assorted Molars eClass 14 Classification Group: 34111100

3M™ Pentamix™ 3 Precision Starter Kit — Imprint™ 4 Penta™

3M ID B10397646

Typical properties

Brands: Pentamix™

3M™ Visio™ Beta Vario Replacement Light, 77170

Part Number 77170 3M Product Number 77170 3M ID 7000054716

Details

Ideal for polymerization of veneers made from 3M[™] Sinfony[™] Indirect Lab Composite to guarantee veneer strength and color stability

Mirrored polymerization chamber provides complete illumination for deep and thorough polymerization

Vacuum technique removes oxygen-inhibited layers for improved surface hardness

Includes three additional programs for versatile use, e.g., for curing MultiTray Light-Cured Maxillary and Mandibular Baseplates

Typical properties

Brands: Visio™

Catalog Number: 77170

Factory ISO Certification: No ISO Certification information available

Product Type: Light Unit

3M[™] Ketac[™] Universal Aplicap[™] Glass Ionomer Restorative, 61000 Series

Details

Enables fast, one-step, stress-bearing restorative solutions
Low stickiness allows an easy placement
Extended indications for use compared to other glass ionomer materials
No need for conditioner or coating
Bulk fill—no need for layering
Continuously releases fluoride for 24 months
Self-adhesive
Self-cure

Typical properties

Brands: Ketac™

Catalog Number: 61089,61088,61098,61087,61097,61086,61096,61085,61095

Cure Type: Self-cure

Delivery System: Aplicap™ Capsule, Capsule

Factory ISO Certification: No ISO Certification information available

Package Contents: Refill includes 300 capsules in A3 shade, Refill includes 50 capsules in assorted shades, Refill includes 50 capsules in A2 shade, Refill includes 50 capsules in A3.5 chade Refill includes 50 capsules in A3.5 chade Refill includes 50 capsules in A3 shade Refill includes 50 capsules in A3.5

shade, Refill includes 50 capsules in A3 shade, Refill includes 50 capsules in A1

shade,Refill includes 50 capsules in A4 shade

Shade: A1, A2, A3, A3.5, A4, Multiple Shades, White

3M™ RelyX™ Ultimate Adhesive Resin Cement Mixing Tips

Details

For easy, automatic mixing of RelyX Ultimate Adhesive Resin cement

Easy to use

Two component system - one bottle, one syringe

Total-etch, selective-etch or self-etch

Dual-cure system

High wear resistance and tooth-like fluorescence

Color stability

Covers entire spectrum of indirect indications

Typical properties

Brands: RelyX™

Catalog Number: 56885,56886,56883

Factory ISO Certification: No ISO Certification information available

Package Contents: 1 Bag - 10 ea Wide & Endo Mix Tips,1 Bag - 10 ea Wide &

Intra Oral Mix Tips, 30- Regular Tips

Product Type: Tip

3M™ Pentamix™ Lite Starter Kit — Imprint™ 4 Penta™

3M ID B10397641

Typical properties

Brands: Pentamix™

3M™ Pentamix™ Lite Starter Kit — Position™ Penta™

3M ID B10397692

Typical properties

Brands: Pentamix™

3M™ Unitek™ Archwire Hooks and Stops Wire Split Crimpable

Details

Light loads 2 mm length Wall: round .015 in. (.38 mm) Ideal with Nitinol wire

Typical properties

Brands: Unitek™

Compatible Archwire Size (Imperial): 0.016 - 0.018 in, 0.018 - 0.025 in, 0.022 -

0.028 in, 0.022 in, 0.03 in, 0.036 in, 0.04 in

Factory ISO Certification: No ISO Certification information available

Shape: Rectangular, Round

3M™ Dispensing Wells, 712-074, 2-cavity, 4/Pk

3M ID 7100240361

Details

2-cavity dispensing wells 4/pack

Typical properties

Brands: 3M™

Factory ISO Certification: No ISO Certification information available

3M™ Ketac™ Fil Plus Aplicap™ Glass Ionomer Restorative Refill, 55000 Series

Details

Provides molecular bond to enamel and dentin to guarantee a tight marginal seal

Typical properties

Brands: Ketac™

Catalog Number: 55050,55060,55010,55020,55030,55040

Cure Type: Self-cure

Delivery System: Aplicap™ Capsule

Factory ISO Certification: No ISO Certification information available

Package Contents: Refill includes 50 capsules in assorted shades: A1, A2, A3, A35 and A4,Refill includes 50 capsules in shade A1,Refill includes 50 capsules in shade A2,Refill includes 50 capsules in shade A35,Refill includes 50 capsules in

shade A3,Refill includes 50 capsules in shade A4 Shade: A1, A2, A3, A3.5, A4, Multiple Shades eClass 14 Classification Group: 34110000

3M™ Paradigm™ VPS Impression Material

Details

Easy extrusion and easy to use Hydrophilic for void-free impressions High tear strength and recovery from deformation Great flow, yet does not drip or slump

Typical properties

Brands: Paradigm™ Catalog Number:

36716,36717,5311,5312,5511,5313,5313J,5315,5513,5513J,5515 Delivery System: Garant™ Dispenser, Pentamix™ Mixing Unit Factory ISO Certification: No ISO Certification information available

Impression Type: Final

Inner Packs per Case: 10 Packs per Case

Material: VPS

Package Contents: 2 - 50ml cartridges, 2 - base pastes, 300ml each; 2 - catalysts,

60ml each,30 - 50ml cartridges

Product Color: Blue, Green, Light Blue, Lilac, Mauve, Orange, Purple, Yellow

Set Time: 2.5 min, 4 min

Set Type: Fast Set, Regular Set Setting Time: 2.5 Minute, 4 Minute

Viscosity: Heavy, Light, Medium, Regular

Work Time: 1 Minute 15 Seconds, 1 min, 2 min, 40 sec

3M™ Unitek™ Gemini Clear Ceramic Brackets

Details

True-Twin Design – improved rotational control
Ample Under Tie-Wing Area – easy ligation
Torque-in-Base – level slot line-up
Base Flange – easier bracket placement
Micro-crystalline Base – reliable bond strength
3M™ APC™ II Adhesive – fewer bonding steps saves time

Typical properties

Adhesive Precoat: Non-coated Angulation: +2°, +3°, +4°, +8°, 0°

Arch: Lower, Upper Bracket Style: Twin Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Hook Included: Yes Hook Type: Distal Material: Ceramic

Package Contents: 10 Bkts/Pk,1/Ea

Prescription: MBT™

Side: Left, Left/Right, Right Slot Size (Imperial): 0.022 in Slot Size (Metric): 0.559 mm

Tooth: Anterior, Bicuspid, Central Incisor, Cuspid, First Bicuspid, Lateral Incisor,

Second Bicuspid

Torque: +10°, +17°, -12°, -17°, -6°, -7°, 0°

Units per Pack: 10

3M™ Unitek™ Prestige Omega Loop Forming Pliers, 900-707, 1/Ea

3M ID 7000025662

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

3M™ Victory Series™ Low Profile Brackets

Details

MBT™ System for treatment control and consistent patient outcomes

Uncompromised mesial-distal dimension

Micro-etch mesh base surface with three point contact design and anatomical contouring

Vertical scribe line for accurate bracket placement

Uncompromising under tie-wing area for ease of ligation

Smooth radius surface

Color-coded quadrant ID system for bracket identification

Typical properties

Adhesive Precoat: APC™ Flash-Free, APC™ II, APC™ PLUS, Non-coated

Angulation: +10°, +2°, +3°, +4°, +5°, +6°, +8°, 0°

Arch: Lower, Upper Base Type: Weldable Bracket Style: Twin Brands: Victory Series™

Factory ISO Certification: No ISO Certification information available

Hook Included: Yes Hook Type: Distal Material: Stainless Steel

Offset Base: Yes

Prescription: MBT™, Roth* Side: Left, Left/Right, Right

Slot Size (Imperial): 0.018 in, 0.022 in Slot Size (Metric): 0.457 mm, 0.559 mm

Tooth: Anterior, Bicuspid, Central Incisor, Cuspid, First Bicuspid, Lateral Incisor,

Second Bicuspid

Torque: +10°, +12°, +17°, +22°, +7°, -10°, -12°, -17°, -1°, -22°, -3°, -6°, -7°, 0°

Units per Pack: 5

3M™ Pentamix™ Lite Starter Kit — Paradigm™ Penta™

3M ID B00007786

Typical properties

Brands: Pentamix™

3M™ Cavit™ Temporary Filling Material

Details

Proven filling material for temporary fillings
Simple to apply with filling instrument
Slight expansion of the filling material ensures a well-sealed margin

Typical properties

Brands: Cavit™

Cure Type: Self-cure

Delivery System: Jar, Tube

Factory ISO Certification: No ISO Certification information available

Package Contents: Refill includes 10, 7g tubes, Refill includes 1, 28g jar, Refill

includes 4, 7g tubes Shade: Gray, Pink, White

3M™ AlastiK™ Separator Radiopaque, 406-080 Series

Details

Not made with natural rubber latex Available in Dusky Blue Loose modules, 100/pack or 1,000/pack

Typical properties

Brands: AlastiK™

Factory ISO Certification: No ISO Certification information available

Product Color: Dusky Blue

eClass 14 Classification Group: 34110603

Dimensions and Classifications

Overall Diameter (Imperial): 0.19 in Overall Diameter (Metric): 0.19 mm

3M™ Unitek™ Prestige Bond Removing Pliers, Replacement Pads, 6 mm., 900-719, 4/Pk

Part Number 900-720 3M ID 7000025674

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

3M™ SmartClip™ Appliance Detail Model (15X), 600-226, 1 Ea

3M ID 7000024229

Typical properties

Bracket Style: Self-ligating Brands: SmartClip™

Factory ISO Certification: No ISO Certification information available

3M™ Typodont Hinge Replacements, 600-209, 1/Ea

3M ID 7100154547

Details

1 set

Typical properties

Brands: 3M™

Factory ISO Certification: No ISO Certification information available

3M™ Band Removing Pliers, Replacement Pads

Details

Replacement pads 3M™ Unitek™ Prestige Series

Typical properties

Brands: 3M™

Factory ISO Certification: No ISO Certification information available

Product Type: Band

3M™ Ketac™ Nano Light-Curing Glass Ionomer Quick Mix Capsule Refill, 3305

Details

The first paste/paste, resin-modified glass ionomer (RMGI) restorative based on bonded nanofiller technology
Esthetic glass ionomer restorative
High fluoride release
Improved wear resistance
Light cure on demand
Radiopaque

Typical properties

Brands: Ketac™

Catalog Number: 3305A4,3305A3.5,3305A1,3305B2,3305A2,3305A3

Cure Type: Light Cure Delivery System: Capsule

Factory ISO Certification: No ISO Certification information available

Package Contents: Refill includes 20 - A3 Quick Mix Capsules,Refill includes 20 - A2 Quick Mix Capsules,Refill includes 20 - A4 Quick Mix Capsules,Refill includes 20 - A35 Quick Mix Capsules,Refill includes 20 - B2 Quick Mix Capsules,Refill

includes 20 - A1 Quick Mix Capsules Shade: A1, A2, A3, A3.5, A4, B2 Total Units per Case: 10 per Case

3M™ Unitek™ Twin Torque™ Brackets

Typical properties

Angulation: 0° Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Torque: -11°, -17°, -22°, -7° eClass 14 Classification Group: 34110601

3M™ Arch-Making Pencil Kit, 603-026, 1/Kit

3M ID 7000025589

Details

1 kit

Typical properties

Brands: 3M™

Factory ISO Certification: No ISO Certification information available

3M™ Unitek™ MIA Mobile Intraoral Arch Angle, 340-018, 5/Pk

Part Number 340-018 3M ID 7100268466

Details

Is removable and easily activated

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Units per Inner Pack: 5 per Pack

3M™ Unitek™ Pliers/Cutter Combo, 800-705, 1/Ea

Part Number 800-734 3M Product Number 800-734 3M ID 7000025651

Details

Lightweight stainless steel pliers provide easy access to difficult-to-reach areas Features springback action

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

3M™ Unitek™ Silver Solder

Details

Dead soft Contains no toxic cadmium 2-DWT package

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

3M™ Ketac™ Silver Aplicap™ Silver Glass Ionomer Restorative, **37000 Series**

3M ID B5005237012

Details

High sustained fluoride release Radiopaque Self-cure One-step placement

Typical properties

Brands: Ketac™

Delivery System: Aplicap™ Capsule eClass 14 Classification Group: 34110000

3M™ Crown Pliers

Details

Ideal for enhancing crown form contours to improve interproximal contacts and gingival margins for stainless steel and temporary crowns

Smaller size to meet your procedural needs

Typical properties

Brands: 3M™

Catalog Number: 800112,800421,800417

Factory ISO Certification: No ISO Certification information available

Pliers Type: Contouring, Crimping

3M™ Imprint™ 3 VPS Impression Material

Details

Hydrophilicity and flow behavior allow for void-free impressions
Elastic recovery helps prevent distortion for a more accurate impression
Resists tearing upon removal thanks to material toughness
Recovery from deformation helps to reduce retakes
Available in a full range of viscosities, deliveries and setting times
Control moisture at every step and gently retract gingival tissue for more accurate impressions with 3M™ Astringent Retraction Paste

Typical properties

Brands: Imprint™ Catalog Number:

10771J,10771PP,10971J,10771PH,10778,10978,10779,10979,10781,10971PH,

10771,10971,36955,36956

Delivery System: Garant™ Dispenser, Pentamix™ Mixing Unit Factory ISO Certification: No ISO Certification information available

Impression Type: Final

Material: VPS

Package Contents: 4 - 50ml cartridges; 5 - 3M[™] Intra-Oral Syringes - Green; 5 - 3M[™] Garant[™] Mixing Tips - Yellow,4 - 50ml cartridges; 5 - 3M[™] Intra-Oral Syringes - Green; 5 - 3M[™] Garant[™] Mixing Tips - Yellow,4 - 50ml cartridges; 5 - 3M[™] Intra-Oral Syringes -Green; 5 - 3M[™] Garant[™] Mixing Tips - Yellow,4 - 50ml cartridges; 10 - 3M[™] Garant[™] Mixing Tips - Green,4 - 50ml cartridges; 20 - 3M[™] Garant[™] Mixing Tips - Green,4 - 50ml cartridges; 10 - 3M[™] Garant[™] Mixing Tips - Yellow; 5 - 3M[™] Intra-Oral Syringes - Green,30 - 50ml cartridges,2 - base paste, 300ml each; 2 - catalyst, 60ml each

Physical Strength: Higher, Medium

Product Color: Blue Green, Caramel, Green, Light Blue, Maroon, Pink, Royal

Blue, Teal, Yellow

Set Time: 2.5 min, 3 Minutes, 3.5 min, 4 min

Set Type: Quick Set, Regular Set

Setting Time: 2.5 Minute, 3 Minute, 3.5 Minute, 4 Minute

Total Units per Case: 10 per Case

Viscosity: Heavy, Light, Medium, Regular, Ultra-regular

Work Time: 1 Minute 15 Seconds, 1 Minute 30 Seconds, 2 min

3M™ SmartClip™ SL3 Self-Ligating Brackets

Details

Proprietary self-ligating clips which can save valuable chair time during wire engagement and disengagement

"Intelligent" clips that release the archwire if forces are excessive

Provide selective ligation capability, allowing clinicians to engage the archwire in a single clip

Precisely finished nickel titanium clips and smooth radiused surfaces for enhanced patient comfort

Offered in the MBT[™] System prescription, as well as Ricketts (High Torque) and Roth (Low Torque) for desired treatment philosophy options

3M[™] APC[™] Flash-Free and 3M[™] APC[™] II Adhesive coating options for bonding efficiency

Molar brackets with small, low profile tie-wings for treatment flexibility

Notched tie-wings for additional bicuspid identification

Low profile for patient comfort

Torque-in-base for level slot line up

Micro-etched, three point contact base anatomical contouring for better base to tooth fit

Vertical scribe line for accurate bracket placement

The traditional tie-wings allow for application of auxiliaries

The traditional tie-wings allow for application of auxiliaries

Typical properties

Adhesive Precoat: Non-coated

Angulation: +2°, +3°, +4°, +5°, +8°, +9°, 0°

Arch: Lower, Upper Bracket Style: Twin Brands: SmartClip™

Factory ISO Certification: No ISO Certification information available

Hook Included: Yes Hook Type: Distal Material: Stainless Steel

Offset Base: Yes

Prescription: High Torque, Low Torque, MBT™, Medium Torque, Ricketts*, Roth*

Side: Left, Left/Right, Right

Slot Size (Imperial): 0.018 in, 0.022 in Slot Size (Metric): 0.457 mm, 0.559 mm

Tooth: Anterior, Bicuspid, Central Incisor, Cuspid, First Bicuspid, First Molar,

Lateral Incisor, Second Bicuspid

Torque: +10°, +12°, +15°, +17°, +22°, +3°, +5°, +6°, +7°, -12°, -14°, -17°, -1°, -20°, -3°, -4°, -6°, -7°, -9°, 0°

Units per Pack: 4,5,20

3M™ Victory Series™ Low Profile Bracket System

3M ID B5005051042

Details

A unique rhomboid design with a flatter, compounded radius accommodates naturally occurring variation in tooth surfaces. This design assures three point contact for increased tooth adaptation

Mesial-distal dimension is not compromised even though the bracket has a low profile. You maintain the rotational control you need to get the results you desire Delivers uncompromising levels of control throughout the treatment process, assuring consistent patient outcomes with the use of the MBT™ System Positioning is made easier and more accurate because of their rhomboid shape that follows the natural vertical line of the tooth

Ample under tie-wing area is integral to the design so you do not have to sacrifice bracket size for easy ligation

Smooth radius surface for enhanced patient comfort

Torque in Base for level slot line-up

Uncompromising under tie-wing area for ease of ligation

Color-coded quadrant ID system

Vertical scribe line for more accurate bracket placement

Micro-etch base surface and anatomical contouring

Typical properties

Brands: Victory Series™

3M™ Sondhi™ Rapid-Set Indirect Bonding Adhesive, 712-070, 1/Kit

3M ID 7000004409

Details

Two-part resin system
No mixing of Resin A and Resin B
Chemical cure
Give your patient a faster, more comfortable, more confident bonding experience

Typical properties

Brands: Sondhi™

Catalog Number: 712-070 Cure Type: Chemical Cure Delivery System: Bottle

Factory ISO Certification: No ISO Certification information available

Product Type: Indirect Bonding Adhesive Kit

Units per Inner Pack: 1 /KT

3M[™] Transbond[™] IDB Pre-Mix Chemical Cure Adhesive Kit, 712-120, 1/Kit

3M ID 7000024142

Details

Can be used on labial or lingual indirect bonding systems
Strong enough so that the wire can be inserted as soon as three minutes after
the indirect bonding tray is removed
Ample work time and quick set time
Indirect bonding adhesive for the 3M™ Incognito™ Appliance System

Typical properties

Brands: Transbond[™]
Catalog Number: 712-120
Cure Type: Chemical Cure
Delivery System: Bottle

Factory ISO Certification: No ISO Certification information available

Product Type: Chemical Cure Adhesive Kit

Units per Inner Pack: 1 /KT

3M™ Peridex™ Chlorhexidine Gluconate 0.12% Oral Rinse

Details

Provides antimicrobial activity during oral rinsing

Reduces certain aerobic and anaerobic bacteria from a range of 54-97% through six months use

Shows no significant changes in bacterial resistance or adverse changes in the oral microbial ecosystem

Prescription only

Typical properties

Active Ingredient: Chlorhexidine

Brands: Peridex™ Condition: Gingivitis Delivery System: Bottle

Factory ISO Certification: No ISO Certification information available

Flavor: Mint

Material: Chlorhexidine

3M™ Protemp™ Plus Temporization Material Refill

Details

Record-breaking fracture resistance and toughness

Sophisticated filler technology provides strength and esthetics without polishing or glazing

Six natural-looking shades designed to match Filtek™ Supreme Ultra Flowable Restorative

Quick intraoral set

Total set time of only five minutes

Easy plaque removal to support gingival health

Typical properties

Brands: Protemp™

Catalog Number: 46955,46963,46964,46961,46973,46958

Delivery System: Garant™ Dispenser

Factory ISO Certification: No ISO Certification information available

Material: Bis-acrylic Composite

Number of Crowns: 12

Pre-formed: No

Shade: A1, A2, A3, A3.5, B3, Bleach eClass 14 Classification Group: 34110000

3M™ RelyX™ Luting Plus Cement Kit

Details

5-second tack light cure feature for clean-up on demand Higher bond strength than conventional glass ionomers Sustained fluoride release Virtually no post-op sensitivity Moisture tolerant

Interested in further simplifying your direct and indirect restorative workflows with a fully aligned system? Explore 3M™ RelyX™ Universal Resin Cement

Typical properties

Brands: RelyX™

Catalog Number: 3525TK Cure Type: Self-cure

Delivery System: Clicker™ Dispenser

Factory ISO Certification: No ISO Certification information available

Package Contents: 1 – 8.5g Syringe; 12 Mixing Tips,1 – 11g Clicker™ Dispenser

Units per Case: 5.0

Units per Inner Pack: 1 per Pack

3M™ Adper™ Prompt™ L-Pop™ Self-Etch Adhesive Refill, 41925, 40ct

3M ID 7000030679 UPC 04035077001915

Details

Self-etch adhesive Indicated for direct restorations Unit dose delivery Reduces risk of post-op sensitivity Two liquid components are mixed prior to application

Typical properties

Adhesive Type: Self-etch

Brands: Adper™

Catalog Number: 41925 Cure Type: Light Cure Delivery System: L-Pop™

Factory ISO Certification: No ISO Certification information available

Package Contents: 40 L-Pop™ applicators

Product Type: Dental Adhesive Units per Inner Pack: 40 per Pack

3M™ Imprint™ 4 Preliminary VPS Impression Material

Details

High surface quality and long-term dimensional stability ensure well-fitting temporaries

Multiple pours are possible without losing accuracy

Can be thoroughly disinfected while maintaining dimensional stability

Automated mixing makes your entire procedure faster and more predictable

Scannable for easy access to the digital workflow

For successful indirect procedures, capture precision impressions with 3M[™] Imprint[™] 4 VPS Impression Material or 3M[™] Impregum[™] Impression Material

Typical properties

Brands: Imprint™

Catalog Number: 71548,71523,71524,71521,71522,71540,71549,71539

Delivery System: Garant™ Dispenser, Pentamix™ Mixing Unit Factory ISO Certification: No ISO Certification information available

Impression Type: Preliminary

Material: VPS

Package Contents: 2 - base pastes, 300ml each; 2 - catalysts, 60ml each, 10 - 50ml cartridges, 30 - 50ml cartridges, 6 - base pastes, 300ml each; 6 - catalyst

pastes, 60ml each Product Color: Pink

Set Time: 1.5 Minutes, 2 min

Set Type: Regular Set, Super Quick Set

Setting Time: 1.5 Minute, 2 Minute

Viscosity: Medium

Work Time: 1 Minute 30 Seconds, 1 min eClass 14 Classification Group: 34110100

3M™ Pentamix™ Lite Automatic Mixing Unit - 120V/60Hz NA, 77907

3M ID 7100196264 UPC 04054596646339

Details

Offers a lightweight, portable design Efficient, automated device consistently delivers a void-free, homogeneous mix To mix 3M impression materials twice as fast, check out the 3M™ Pentamix™ 3 Automatic Mixing Unit

Typical properties

Brands: Pentamix™ Catalog Number: 77907

Factory ISO Certification: No ISO Certification information available

Package Contents: 1 - Pentamix[™] Lite 120V 60HZ Automatic Mixing Unit; 1 - 3M[™] Penta[™] Cartridge, 1 - Quick Start Guide, 1 - IFU, 1 - 3M[™] Pentamix[™] Lite

Automatic Mixing Unit Sticker

Product Type: Automatic Mixing Unit eClass 14 Classification Group: 34110100

3M™ Sof-Lex™ Extra-Thin Contouring and Polishing Discs Kit, 2380

3M ID 7000053446 UPC 30605861011793

Details

Easy-to-use system designed to create high-luster restorations that offer lifelike esthetics

Discs are color coded to allow for easy identification of grit sequence Pop-on disc technology designed to save you time Extra-thin (XT) discs allow for easy access to interproximal surfaces while offering a slightly stiffer disc

Typical properties

Brands: Sof-Lex™

Factory ISO Certification: No ISO Certification information available

Grade: Assorted

Model: XT

Package Contents: Kit includes 30 each Extra-thin discs in coarse, medium, fine,

and superfine grits in both 1/2" and 3/8" diameters, and 1 RA mandrel

Product Type: Extra-thin Contouring & Polishing Disc

3M™ Paradigm™ Nano Hybrid Universal Restorative, 1480 Series

Details

Easy to use for all classes of direct restorations

Easy to polish for predictable esthetics

Great handling without slumping or sticking to instruments

High strength, with good wear resistance

Available in 12 radiopaque, fluorescent shades

Complete your premium restoration with the gold standard of adhesion, 3M™

Scotchbond™ Universal Plus Adhesive

Typical properties

Brands: Paradigm™ Catalog Number:

1481A3.5,1481B2,1480A2,1481A1,1480A3,1481A2,1481B3,1480A4,1481A3,14 80A3.5,1481A4,1480OA2,1481OA3,1481OA2,1480D3,1480B1,1480C2,1480A1,

1480B2,1481D3,1481B1,1481C2

Composite Type: Universal Cure Type: Light Cure

Delivery System: Capsule, Syringe

Factory ISO Certification: No ISO Certification information available

Material Structure: Nanohybrid Composite

Net Weight (Metric): 0.2 g, 3 g

Opacity: Body, Dentin

Package Contents: Refill includes 1 - 3g syringe in Body Shade A3.5,Refill includes 1 - 3g syringe in Body Shade A3,Refill includes 20 - 02g capsules in Opaque A3 Shade,Refill includes 1 - 3g syringe in Body Shade A4,Refill includes 20 - 02g capsules in Opaque A2 Shade,Refill includes 20 - 02g capsules in Body shade A3,Refill includes 20 - 02g capsules in Body shade A2,Refill includes 20 - 02g capsules in Body shade B3,Refill includes 20 - 02g capsules in Body shade A4,Refill includes 1 - 3g syringe in Body Shade A1,Refill includes 1 - 3g syringe in Body Shade A2,Refill includes 20 - 02g capsules in Body shade B2,Refill includes 20 - 02g capsules in Body shade B1,Refill includes 1 - 3g syringe in Body Shade D3,Refill includes 1 - 3g syringe in Body Shade C2,Refill includes 20 - 02g capsules in Body shade B1,Refill includes 20 - 02g capsules in Body shade B1,Refill includes 20 - 02g capsules in Body shade C2,Refill includes 20 - 02g capsules in Body shade B1,Refill includes 20 - 02g capsules in Body shade C2,Refill includes 20 - 02g capsules in Body shade C2,Refill includes 20 - 02g capsules in Body shade C2,Refill includes 20 - 02g capsules in Body shade C2,Refill includes 20 - 02g capsules in Body shade C2,Refill includes 20 - 02g capsules in Body shade C2,Refill includes 20 - 02g capsules in Body shade C2,Refill includes 20 - 02g capsules in Body shade C2,Refill includes 20 - 02g capsules in Body shade C2,Refill includes 20 - 02g capsules in Body shade C2,Refill includes 20 - 02g capsules in Body shade C2,Refill includes 20 - 02g capsules in Body shade C2,Refill includes 20 - 02g capsules in Body shade C2,Refill includes 20 - 02g capsules in Body shade C2,Refill includes 20 - 02g capsules in Body shade C2,Refill includes 20 - 02g capsules in Body shade C2,Refill includes 20 - 02g capsules in Body shade C2,Refill includes 20 - 02g capsules in Body shade C2,Refill includes 20 - 02g capsules in Body shade C2,Refill includes 20 - 02g capsules in Body shade C2,Refill includes 20 - 02g capsules in Body shade C2,Ref

Shade: A1, A2, A3, A3.5, A4, B1, B2, B3, C2, D3, OA2, OA3

Units per Inner Pack: 1 per Pack, 20 per Pack

Units per Pack: 1, 20

3M™ Filtek™ P60 Posterior Restorative

Details

Excellent packability — packs well without sticking to instruments

Easy sculpting with exceptional resistance to slumping

Excellent marginal adaptation and interproximal contact

Low polymerization shrinkage

Excellent wear resistance and strength means you can be confident of long-

lasting posterior restorations

Fast placement technique takes only 20 seconds to cure an increment of 2.5mm

Want to provide your patients with both long-lasting performance and uncompromising esthetics? Explore 3M™ Filtek One Bulk Fill Restorative

Typical properties

Brands: Filtek™

Catalog Number: 8100C2,8100B2,8100A3,4720A3,4720C2,4720B2

Composite Type: Universal Cure Type: Light Cure Delivery System: Syringe

Factory ISO Certification: No ISO Certification information available

Inner Packs per Case: 10 Packs per Case Material Structure: Microhybrid Composite

Net Weight (Metric): 4 g

Package Contents: Refill includes 1 - 4g syringe in shade A3, Refill includes 1 -

4g syringe in shade C2, Refill includes 1 - 4g syringe in shade B2

Shade: A3, B2, C2

Total Units per Case: 10 per Case Units per Inner Pack: 1 per Pack

Units per Pack: 1

3M™ Lava™ Plus High Translucency Zirconia Disc

Details

High strength for tooth-preserving preparations and all crown and bridge indications

Scientific data proves Lava all-zirconia has antagonist preserving wear properties and is not abraded or roughened under occlusal load

Lava™ Plus Zirconia shading system enables highly esthetic all-zirconia or traditional layered restorations

Excellent aging stability

Easy adjustment

Easy to cement. Lava Plus zirconia restorations can be conventionally cemented, using a cement like 3M™ RelyX™ Unicem 2 Self-Adhesive Resin Cement

Typical properties

Block Size: S Brands: Lava™

CAD/CAM System: Lava

Catalog Number: 69273,69272,69271 Digital Materials: Lava Laboratory Materials

Factory ISO Certification: No ISO Certification information available

Material: Zirconia

Restoration Type: Bridge, Crown, Inlay, Onlay, Veneer

Shade: Multiple Shades

Dimensions and Classifications

Overall Diameter (Imperial): 3.858 in Overall Diameter (Metric): 98 mm

Overall Height (Imperial): 0.551 in, 0.709 in, 0.984 in, 1.181 in

Overall Height (Metric): 14 mm, 18 mm, 25 mm, 30 mm

3M™ OMNI Gel™ 0.4% Stannous Fluoride Brush-On Gel Refill, 12106

Details

Formulated with stannous fluoride to help your patients prevent tooth decay

Brush-on therapy designed for caries-prone patients

Provides antimicrobial activity for up to 8 hours

Promotes tooth remineralization

Helps inhibit demineralization for lasting oral health

Offered in Mint flavor

Typical properties

Active Ingredient: Stannous Fluoride

Brands: 3M™

Delivery System: Bottle, Tube

Factory ISO Certification: No ISO Certification information available

Features: Anti-cavity, Hypersensitivity

Flavor: Mint

Fluoride Content: 0.4 % Material: Stannous Fluoride Net Weight (Imperial): 4.3 oz Net Weight (Metric): 121.903 g

Product Form: Gel

Product Type: Brush-on Gel

Units per Inner Pack: 1 per Pack, 24 per Pack eClass 14 Classification Group: 25092006

3M™ Clinpro™ XT Varnish Durable Fluoride-releasing Coating

3M ID B00005785

Details

Delivers fluoride, calcium and phosphate
Creates a barrier against demineralization and acid erosion* *Based on in vitro
data
Helps repair demineralized enamel* *Based on in vitro data
Provides protection beyond the coating
Coating remains for 6 months or longer
Moisture-tolerant

Typical properties

Active Ingredient: Sodium Fluoride

Brands: Clinpro™

3M™ Express™ STD VPS Impression Material Putty Refill

Details

VPS impression material for 1-step and 2-step putty/wash techniques Provides accuracy and detail reproduction Multiple-pour capability Good dimensional stability

Typical properties

Brands: Express™ Catalog Number: 7312 Delivery System: Hand Mix

Factory ISO Certification: No ISO Certification information available

Impression Type: Final

Inner Packs per Case: 10 Packs per Case

Material: VPS

Package Contents: 1 - jar base; 1 - jar catalyst, 305ml each

Product Color: Orange

Set Time: 5 min

Set Type: Regular Set Setting Time: 5 Minute Work Time: 1 min

3M™ Pentamix™ 3 Automatic Mixing Unit, 77876

3M ID 7100215403 UPC 04054596735729

Details

Get reliable impressions and more working time with automated, void-free, homogeneous and consistent mixing of dental impression material — including putty materials

Automated mixing device provides reliable mixing for polyether, VPS and alginate replacement material consistencies

Enables fast, automated mixing for shorter impression procedures Convenient and ergonomic with buttons on both sides of the unit for easy operation whether you're right- or left-handed

Offers easy cartridge exchange and features easy-to-attach mixing tips

Typical properties

Brands: Pentamix™ Catalog Number: 77876

Factory ISO Certification: No ISO Certification information available

Package Contents: 1 - Pentamix™ 3 Automatic Mixing Unit

Product Type: Automatic Mixing Unit eClass 14 Classification Group: 34110100

3M™ Gold Anodized Temporary Crowns, 940000 Series

Details

Made from medium-hard aluminum for good durability
Gold anodization minimizes metallic taste and galvanic shock
Bicuspid and molar crown selection to cover adult posterior applications
Crown set is packaged in a unique, stackable tray that provides convenient, easy storage

Typical properties

Arch: Lower, Upper Brands: 3M™
Catalog Number:

 $940532,940334,940531,940333,940530,940332,940331,940526,940525,940646\\,940645,940524,940326,940644,940446,940523,940325,940643,940445,94052\\2,940324,940642,940444,940323,940521,940641,940443,940322,940520,9406\\40,940442,940321,940441,940516,940636,940515,940635,940514,940316,940\\634,940436,940513,940315,940633,940435,940314,940512,940434,940632,94\\0313,940511,940519,940433,940631,940510,940312,940630,940432,940311,9\\40431,940546,940545,940544,940346,940509,940543,940345,940542,940344,\\940541,940343,940342,940540,940341,940536,940535,940336,940534,940335,940533$

Dentition: Permanent

Factory ISO Certification: No ISO Certification information available

Material: Aluminum

Mesial/Distal Width (Imperial): 0.22 in, 0.224 in, 0.24 in, 0.244 in, 0.26 in, 0.28 in, 0.291 in, 0.299 in, 0.319 in, 0.339 in, 0.358 in, 0.37 in, 0.378 in, 0.39 in, 0.402 in, 0.409 in, 0.421 in, 0.429 in, 0.441 in, 0.449 in, 0.461 in, 0.469 in, 0.48 in, 0.488 in Mesial/Distal Width (Metric): 5.6 mm, 5.7 mm, 6.1 mm, 6.2 mm, 6.6 mm, 7.1 mm, 7.4 mm, 7.6 mm, 8.1 mm, 8.6 mm, 9.1 mm, 9.4 mm, 9.6 mm, 9.9 mm, 10.2 mm, 10.4 mm, 10.7 mm, 10.9 mm, 11.2 mm, 11.4 mm, 11.7 mm, 11.9 mm, 12.2 mm, 12.4 mm

Number of Crowns: 5 Side: Left, Right

Tooth: First Bicuspid, First Molar, Second Bicuspid, Second Molar

3M™ Lava™ Esthetic Fluorescent Full-Contour Zirconia Disc

Details

First zirconia with inherent, tooth-like fluorescence for high-esthetic full-contour crown and bridge cases.*

High translucency. Optimized for esthetic full-contour restorations

Shade match. Built-in shade gradient, designed for excellent shade match for all VITA classical shades including Bleach

Durability. Durable full-contour with 800 MPa through and through. Indicated for crowns, 3-unit bridges, inlays/onlays, veneers.**

Productivity. Benefit from a streamlined mill-sinter-glaze process. No shading liquids. No drying time

Typical properties

Block Size: S Brands: Lava™

CAD/CAM System: Lava

Catalog Number:

69369,69325,69368,69367,69323,69322,69366,69329,69328,69327,69624,69326,69361,69360,69365,69321,69364,69320,69342,69363,69341,69340,69362,69358,69336,69357,69356,69334,69378,69377,69933,69339,69338,69359,69337,69319,69371,69370,69332,69354,69376,69331,69353,69375,69330,69352,69374,69373

Digital Materials: Lava Laboratory Materials

Factory ISO Certification: No ISO Certification information available

Material: Zirconia

Restoration Type: Bridge, Crown, Inlay, Onlay, Veneer

Shade: A1, A2, A3, A3.5, A4, B1, B2, B3, B4, BL, C1, C2, C3, C4, D2, D3, D4

Dimensions and Classifications

Overall Diameter (Imperial): 3.858 in Overall Diameter (Metric): 98 mm

Overall Height (Imperial): 0.551 in, 0.709 in, 0.866 in Overall Height (Metric): 14 mm, 18 mm, 22 mm

3M™ Victory Series™ Buccal Tubes

Details

Trouble-free integral convertible cap

Solid, one piece construction

Easy archwire insertion

Smooth, patient-friendly hooks

Low-profile design

Side gripping notches for ease of handling and placement

Distal notch for ease of ligation

An alignment seating notch for more accurate placement

Includes a patented, trouble-free integral convertible cap that converts easily and

consistently

Designed with funneled and radiused entries on both the headgear tube and

archwire slot

Smooth, patient-friendly hooks improve comfort for your patients

Typical properties

Adhesive Precoat: APC™ II, APC™ PLUS, Non-coated

Arch: Lower, Upper, Upper/Lower Base Type: Bondable, Weldable

Brands: Victory Series™

Distal Offset: 0 Degree, 5 Degree, 6 Degree, 7 Degree, 8 Degree, 10 Degree, 14

Degree, 15 Degree Entry Style: Funneled

Factory ISO Certification: No ISO Certification information available

Hook Included: Yes Hook Type: Distal Material: Stainless Steel

Prescription: Low Torque, MBT™, Medium Torque, Sondhi

Side: Left, Left/Right, Right

Slot Size (Imperial): 0.018 in, 0.022 in Slot Size (Metric): 0.457 mm, 0.559 mm

Tooth: First Molar, Second Molar

Torque: -10°, -14°, -19°, -20°, -25°, -30°, 0°

Tube Style: Convertible, Double Convertible, Double Non-convertible, Single Convertible, Single Non-convertible, Triple Convertible, Triple Non-convertible

Units per Inner Pack: 1 per Pack, 4 per Pack, 5 per Pack

eClass 14 Classification Group: 34110000

Dimensions and Classifications

Overall Length (Metric): 3.6 mm, 3.9 mm, 4 mm, 4.3 mm, 4.6 mm, 5.1 mm

3M™ Incognito™ Appliance System

Details

The Incognito[™] Appliance System features fully customized brackets, archwires and bonding trays to deliver predictable, efficient and effective treatment outcomes

Designed to help orthodontists service the growing patient demand for invisible, comfortable and convenient care

Makes it easy for patients to smile with confidence during and after treatment Customized to the unique prescription of each patient, making it easy to treat a wide range of patients and orthodontic case types

The only 100% customized, invisible fixed orthodontic appliance system on the market today

Designed to make every treatment experience predictable and convenient Uses state-of-the-art technologies, including Computer-Aided Design, Computer-Aided Manufacturing, Rapid Prototyping and Robotics

Low profile bases are custom designed to optimize patient comfort and maximize the delivery of the customized prescription

Typical properties

Brands: Incognito™

Factory ISO Certification: No ISO Certification information available

Material: Gold, Stainless Steel

712-022 Unitek™ Light Sensitive Wells 4/Pk

3M ID 7100154486

Details

Contains 3 cavities

Designed to protect light sensitive materials when the cover slides over the wells

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

3M[™] Clinpro[™] 5000 1.1% Sodium Fluoride Anti-Cavity Toothpaste, 12000 Series

Details

Help your patients improve their everyday oral health with a prescription-strength toothpaste that strengthens teeth and repairs early lesions

Formulated with 1.1% NaF (5000 ppm fluoride ion) and an innovative calcium and phosphate ingredient

Proven remineralization throughout the lesion, not just on the surface Creates acid-resistant mineral in dentin tubules.* *3M Internal data Gentle formulation results in low abrasion, even for exposed root surfaces Patients can simply use this prescription-strength toothpaste in place of their normal brushing regimen, resulting in better compliance

Offered in three great flavors: Vanilla Mint, Spearmint and Bubble Gum

Typical properties

Active Ingredient: Sodium Fluoride

Comes in a convenient 4oz tube

Brands: Clinpro™ Delivery System: Tube

Factory ISO Certification: No ISO Certification information available

Features: Anti-cavity

Flavor: Bubble Gum, Spearmint, Vanilla Mint

Fluoride Content: 1.1 %
Material: Sodium Fluoride
Net Weight (Imperial): 4 oz
Net Weight (Metric): 113.398 g
Package Contents: 1 Tube

Product Form: Paste Product Type: Dentifrice

3M™ Vanish™ 5% Sodium Fluoride White Varnish, 12100 Series

Details

Contains 22,600 ppm fluoride and tri-calcium phosphate (TCP) Relieves hypersensitivity
High fluoride uptake, even after acid attack
Contains an exclusive calcium and phosphate ingredient
Advanced protection – long-lasting coating
Easy to apply, no drying needed
Can be applied to tooth surfaces where plaque is present
Fast, horizontal sweep application
Available in mint, melon and cherry

Typical properties

Active Ingredient: Sodium Fluoride

Brands: Vanish™

Factory ISO Certification: No ISO Certification information available

Flavor: Assorted, Cherry, Melon, Mint

Units per Inner Pack: 50 per Pack, 100 per Pack, 1000 per Pack

3M[™] Protemp[™] Plus Temporization Material, 46962, Shade A2, Introductory Kit

3M ID 7100294975 UPC 04064035090166

Details

Record-breaking fracture resistance and toughness

Sophisticated filler technology provides strength and esthetics without polishing or glazing

Six natural-looking shades designed to match Filtek™ Supreme Ultra Flowable Restorative

Typical properties

Brands: Protemp™ Catalog Number: 46962

Delivery System: Garant™ Dispenser

Factory ISO Certification: No ISO Certification information available

Material: Bis-acrylic Composite

Pre-formed: No Shade: A2

3M™ Durelon™ Carboxylate Luting Cement Intro Kit, 38019

3M ID 7000129019 UPC 04035077012737

Details

Simple technique requires no conditioning and no bonding agents for reduced chair time

Thin film thickness Fluoride release

Typical properties

Brands: Durelon™

Catalog Number: 38019 Cure Type: Self-cure Delivery System: Hand Mix

Factory ISO Certification: No ISO Certification information available

Material: Carboxylate

Package Contents: 1 – 20g Powder, Regular Set; 1 – 14ml Liquid; 1 Mixing Pad

(70144

3M™ Ketac™ Cem Plus Luting Cement Refills, 56930

3M ID B00007467

Details

5-second tack light cure feature for clean-up on demand
Higher bond strength than conventional glass ionomers
Sustained fluoride release
Virtually no post-op sensitivity
Moisture tolerant
Interested in further simplifying your direct and indirect restorative workflows with a fully aligned system? Explore 3M™ RelyX™ Universal Resin Cement

Typical properties

Brands: Ketac™

Catalog Number: 56930 Cure Type: Dual Cure

Delivery System: Clicker™ Dispenser Inner Packs per Case: 10 Packs per Case Material: Resin-modified Glass Ionomer Total Units per Case: 10 per Case Units per Inner Pack: 1 per Pack

3M™ Express™ VPS Impression Materials

Details

The first vinyl polysiloxane automix system used worldwide: for over twenty years, dentists who prefer a hand-mixed putty system have relied on 3M™ Express™ VPS Impression Material to deliver accurate, reliable results Great detail reproduction through superior hydrophilicity and flow behavior Offers technique flexibility, which gives dentists the freedom to use either the two-step or simultaneous putty-wash technique Good dimensional stability Excellent elastic recovery

Typical properties

Brands: Express™

Catalog Number: 7301,7322,7322T,7302,7302T,7301T,36799

Delivery System: Garant™ Dispenser

Factory ISO Certification: No ISO Certification information available

Impression Type: Final

Inner Packs per Case: 5 Packs per Case, 10 Packs per Case

Material: VPS

Package Contents: 4 - 50ml cartridges: 20 - 3M™ Garant™ Mixing Tips -

Yellow, 10 - 50ml cartridges

Product Color: Blue, Green, Purple

Set Time: 3.5 min, 5 min

Set Type: Quick Set, Regular Set Setting Time: 3.5 Minute, 5 Minute

Total Units per Case: 20 per Case, 40 per Case Units per Inner Pack: 2 per Pack, 4 per Pack

Viscosity: Light, Regular Work Time: 1 min, 30 sec

3M™ Penta™ Universal Cartridge for Pentamix™ 2, 76994

3M ID 7100215310 UPC 04054596734227

Details

For consistent, homogeneous and void-free mixing at the touch of a button Compatible with 3M™ Pentamix™ 2 Mixing Unit Reusable, metal-reinforced cartridge Helps get an accurate impression on the first take, reducing the need for costly adjustments, remakes and retakes

Typical properties

Brands: Penta™

Catalog Number: 76994

Factory ISO Certification: No ISO Certification information available

Package Contents: 1 - Penta™ Cartridge

Product Type: Cartridge

3M™ Pentamix™ Mounting Kit Twist Grip, 950916

3M ID 7100001564 UPC 04035077014212

Details

Includes 1 mounting kit twist grip

Typical properties

Brands: Pentamix™

Factory ISO Certification: No ISO Certification information available

Product Type: Case

3M™ Lava™ Milling Bur

Details

For Lava™ milling machines

Typical properties

Brands: Lava™

Bur Type: 10, 4, 5, 6, 7, 8

Factory ISO Certification: No ISO Certification information available

3M™ Disposable Wells

Details

Disposable 96/pack

Typical properties

Brands: 3M™

Factory ISO Certification: No ISO Certification information available

Number of Wells: 1

Package Contents: Refill includes 96 mixing wells eClass 14 Classification Group: 34110000

3M™ Paradigm™ 10 mm Light Guide, 76968

Part Number 76968 3M ID 7000129171 UPC 04035077006897

Details

Wide 10mm light guide illuminates up to 80% more surface compared to leading devices tested

Allows you to cure the entire cavity surface of even larger restorations in one shot Light guide rotates 360° for easy alignment of light tip to restoration

Typical properties

Brands: Paradigm™

Factory ISO Certification: No ISO Certification information available

Package Contents: Includes one 10mm light guide

Product Type: Light Guide

3M™ Paradigm™ DeepCure LED Curing Light, 76974, 120 V NA

3M ID 7100196359 UPC 04054596648111

Details

Economical, lightweight plastic model

A predictable, reliable cure, even at the bottom of the proximal box Better depth of cure when perfect light positioning is difficult

120 minutes of cure time from a full charge

Updated light guide geometry allows easy access to all tooth surfaces Simple two-button operation with one cure mode for all indications Built-in battery with charging plug

Cordless portability with 120 minutes of cure time from a full charge Lightweight, vent-free plastic housing for easy disinfection

Typical properties

Brands: Paradigm™

Factory ISO Certification: No ISO Certification information available

Package Contents: Includes one curing light handpiece, charger, 10mm light

guide, eye shield and three curing discs

Product Type: Curing Light

3M[™] Vitrebond[™] Light Cure Glass Ionomer Liner/Base Refill, 7510 Series

3M ID B00034690

Details

Forms a strong bond to dentin, sealing the tooth to reduce the incidence of microleakage, which helps reduce post-operative sensitivity

Produces zones of inhibition against bacteria commonly found in the oral cavity in-vitro

Sustained fluoride release to help protect the tooth long after the restoration is complete

Fast technique and excellent handling characteristics allow for quick and easy placement

Provides excellent adhesion to 3M dental adhesive systems, whether using composite or amalgam

Typical properties

Brands: Vitrebond™ Cure Type: Light Cure

3M™ Giz Accesory Pack Aplicap, 37160

3M ID 7100215417

Details

Easy activation and application of the Aplicap capsule 3M[™] Aplicap[™] applier provides easy delivery of 3M[™] Aplicap[™] capsule material to the tooth preparation or to the restoration

Typical properties

Brands: Aplicap™

Factory ISO Certification: No ISO Certification information available

Package Contents: Refill includes one 3M[™] Aplicap[™] Activator and one 3M[™]

Aplicap™ Applier

Product Type: Activator/Applier Set

3M™ Bendable Applicator Brushes, 3402, 120ct

3M ID 7100214002 UPC 30605861068827

Details

Disposable for improved infection control
Designed for adhesive application
Disposable applicator brushes packaged in a flip-top box
Refills contain 120 applicator brushes: 30 each of four colors (blue, red, yellow, green)

Typical properties

Brands: 3M™

Factory ISO Certification: No ISO Certification information available

Package Contents: Includes 120 applicator brushes: 30 each of four colors (blue,

red, yellow, green)

Product Type: Tip Applicators

3M™ Garant™ Dispenser 4:1/10:1 Plunger Refill, 77591

3M ID 7100215264 UPC 04054596733756

Typical properties

Brands: Garant™

Catalog Number: 77591

Factory ISO Certification: No ISO Certification information available

Product Type: Plunger Refill

3M™ RelyX™ Luting Plus Cement

Details

5-second tack light cure feature for clean-up on demand Higher bond strength than conventional glass ionomers Sustained fluoride release Virtually no post-op sensitivity

Moisture tolerant

Interested in further simplifying your direct and indirect restorative workflows with a fully aligned system? Explore 3M™ RelyX™ Universal Resin Cement

Typical properties

Brands: RelyX™

Catalog Number: 3525,3525A

Cure Type: Self-cure

Delivery System: Clicker™ Dispenser

Factory ISO Certification: No ISO Certification information available

Material: Resin-modified Glass Ionomer

Package Contents: 3 – 8.5g Syringe; 36 Mixing Tips,2 – 11g Clicker™

Dispensers

Product Type: Luting Cement

Shade: White Opaque

Units per Inner Pack: 2 per Pack

3M™ RelyX™ U200 Self-Adhesive Resin Cement Refill, 56878, A2

3M ID 7000003320 UPC 04035077005302

Details

Cement delivery system provides easy handling and flexible dose dispensing Saves time by eliminating etching*, priming and bonding steps *Except for 2 or 3-unit Maryland bridges and 3-unit inlay/onlay bridges. Selective etching is required to increase surface area for bonding

Virtually no post-operative sensitivity

Easy removal of excess cement

High intrinsic color stability for natural-looking esthetics

Cement resists discoloration

Interested in further simplifying your direct and indirect restorative workflows with a fully aligned system? Explore 3M™ RelyX™ Universal Resin Cement

Typical properties

Brands: RelyX™

Catalog Number: 56878 Cure Type: Dual Cure

Delivery System: Clicker™ Dispenser, Syringe

Factory ISO Certification: No ISO Certification information available

Material: Self-adhesive Resin

Product Type: Self-adhesive Resin Cement

Shade: A2, TR

3M[™] Scotchbond[™] Universal DCA Dual Cure Activator - Vial, 41264, 5ml

3M ID 7000055182 UPC 04035077006637

Details

Combined Total-Etch, Self-Etch and Selective-Etch adhesive Consistent bond strength to both moist and dry etched dentin

Virtually no post-op sensitivity

Effective on all surfaces

High bond strength to all indirect surfaces (zirconia, alumina, glass ceramics and metals) without a separate primer

Excellent marginal integrity for highly esthetic restorations

One-step, one-coat adhesive application

Self- and dual-cure compatible with 3M™ Scotchbond™ Universal DCA Dual

Cure Activator

No refrigeration needed

Store at room temperature for two years

Flip-top vial and unit-dose delivery Improved, easy-to-place etchant

Typical properties

Adhesive Type: Selective Etch, Self-etch

Brands: Scotchbond™ Catalog Number: 41264

Container Size (Imperial): 0.169 fl oz

Container Size (Metric): 5 mL

Cure Type: Dual Cure Delivery System: Vial

Factory ISO Certification: No ISO Certification information available

Package Contents: Refill includes 1 - 5ml vial of activator

Product Type: Dental Adhesive

3M™ Impregum™ Polyether Impression Materials

Details

Unique moisture tolerance results in excellent detail reproduction and void-free impressions, even in wet environments

Constant flowability throughout the whole working time makes it suitable for large and complex cases

Snap-set behavior minimizes risk of distortion for precise-fitting restorations Post-set material rigidity helps ensure accurate transfer of implant copings Available in a range of viscosities, set times and delivery systems for monophase and 1-step heavy/light impression techniques

Control moisture at every step and gently retract gingival tissue for more accurate impressions with 3M[™] Astringent Retraction Paste

Typical properties

Brands: Impregum™ Catalog Number:

P31722,P31779,31786,P31735,P31739,P31738,P31771,P31685
Delivery System: Garant™ Dispenser, Pentamix™ Mixing Unit
Factory ISO Certification: No ISO Certification information available

Impression Type: Final Material: Polyether

Package Contents: 1 -base paste 300ml, 1 catalyst paste 60ml; 1 - 3M™ Penta™ Cartridge; 10 - 3M[™] Penta[™] Mixing Tips, red; 2 - 50ml cartridges Light Body; 10 - 3M™ Garant™ Mixing Tips - White; 10 - 3M™ Garant™ intraoral Tips - White; 1 - 3M[™] Polyether Adhesive, 17ml,1 - base paste, 300ml; 1 - catalyst paste, 60ml; 1 - 3M™ Penta™ Cartridge; 1 - 3M™ Penta™ Elastomer Syringe; 10 - 3M™ Penta[™] Mixing Tips, red; 1 -3M[™] Polyether Adhesive, 17ml,1 -base paste, 300ml; 1 - catalyst paste, 60ml; 1 - 3M[™] Penta[™] Cartridge; 2 - 50ml Light Body cartidges; 1 - 3M™ Elastomer Syringe, 10 - 3M™ Penta™ Mixing Tips, red; 5 -3M™ Garant™ Mixing Tips - White; 5 - 3M™ Garant™ Intraoral Tips -White; 1 -3M[™] Polyether Adhesive, 17ml,1 - base paste, 300ml; 1 - catalyst paste, 60ml; 1 - 3M[™] Penta[™] Cartridge: 1 - 3M[™] Elastomer Syringe: 10 - 3M[™] Penta[™] Mixing Tips, red; 1 - 3M[™] Polyether Adhesive, 17ml,3 - 50ml cartridges, Medium Body; 1 - 50ml cartridge Light Body; 10 - 3M™ Garant™ Mixing Tips - Purple; 5 -3M™ Garant™ Mixing Tips - White; 5 - 3M™ Garant™ Intraoral Tips - White; 5 -3M™ Intra-oral Syringes - Purple, 1 - 3M™ Polyether Adhesive - 17ml,1 - base paste, 300ml; 1 - catalyst paste, 60ml; 1 - 3M™ Penta™ Cartridge; 1 - 3M™ Penta™ Elastomer Syringe; 10 - 3M™ Penta™ Mixing Tips, red; 1 - 3M™

Polyether Adhesive 17ml Product Color: Medium Purple Set Time: 4 min, 6 Minutes Set Type: Quick Set, Regular Set

Setting Time: 4 Minute, 6 Minute

Viscosity: Heavy, Medium, Multiple Viscosities

Work Time: 1 Minute 45 Seconds, 1 min, 2 min 30 sec, 2 min 45 sec

3M™ Sof-Lex™ Finishing and Polishing Discs, Kit, 1980

3M Product Number 1980 3M ID 7000053381 UPC 30605861011779

Details

Easy-to-use system designed to create high-luster restorations that offer lifelike esthetics

Discs are color coded to allow for easy identification of grit sequence Pop-on disc technology designed to save you time

Typical properties

Brands: Sof-Lex™

Factory ISO Certification: No ISO Certification information available

Grade: Assorted Model: Standard

Package Contents: Kit includes 30 each discs in coarse, medium, fine, and

superfine grits in both 1/2" and 3/8" diameters, and 1 RA mandrel

Product Type: Contouring & Polishing Disc eClass 14 Classification Group: 34111407

3M[™] Filtek[™] Supreme Ultra Universal Restorative Master's Kit, 6029M, Capsule

Part Number 6029M 3M ID 7000128914 UPC 30605861012011

Details

Exceptional esthetics
Excellent adaptation
Excellent polish retention
Exceptional handling
Outstanding strength for anterior and posterior
Excellent wear resistance
Natural-looking fluorescence
True nanotechnology
36 shades, 4 opacities

Typical properties

Brands: Filtek™

Catalog Number: 6029M Composite Type: Universal Delivery System: Capsule

Environmental claim being made: Yes Environmental claim relates to: Product

Factory ISO Certification: No ISO Certification information available

Material Structure: Nanocomposite

Opacity: Multiple Opacities

Package Contents: Kit includes: 132 – 02g capsules: 3 each of A1D, A2D, A3D, A4D, B3D, C4D, WD; 6 each of A1B, A2B, A3B, A35B, A4B, B1B, B2B, C2B; 3 each of B3B, B5B, C1B, C3B, D2B, D3B, WB, XWB; 6 of A2E; 3 each of A1E, A3E, B1E, B2E, D2E, WE, XWE; 3 each of Translucent shades AT, BT, CT, GT All in plastic reusable case with 1 Shade Selection Wheel; Technique Guides;

Instructions for Use Shade: Multiple Shades

3M™ Unitek™ Miniature Twin V-Slot Metal Brackets

Details

Permanent raised I.D. markings provide ready identification

A permanent vertical scribe line facilitates more accurate bracket placement on the long axis for bonding

Symmetrical tie-wings add balance to enhance tie-wing strength and improve patient esthetics

An integral base perimeter ridge facilitates easy flash cleanup

Torque-in-base design with compound base contour

Bracket design provides additional reference points to ensure proper bracket placement regardless of tooth position

A chamfered slot entrance promotes easier and more rapid archwire placement

Typical properties

Angulation: +11°, +5°, +9°, 0° Arch: Lower, Upper, Upper/Lower

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Hook Included: Yes Hook Type: Distal

Prescription: Roth*, Standard Edgewise

Side: Left, Left/Right, Right

Slot Size (Imperial): 0.018 in, 0.022 in Slot Size (Metric): 0.457 mm, 0.559 mm

Tooth: Anterior, Bicuspid, Central Incisor, Cuspid, Cuspid/Bicuspid, First

Bicuspid, Lateral Incisor, Second Bicuspid Torque: +12°, +8°, -11°, -17°, -1°, -22°, -7°, 0°

Units per Pack: 5

3M™ ESPE™ Stainless Steel Permanent Molar Crowns, PO-96, Kit

3M ID 7000053377 UPC 30605861011199

Details

Pre-trimmed, belled and crimped for fast and easy placement Superior longevity compared to multi-surface amalgam fillings Wide array of shapes and sizes Made in the U.S.A. of high quality stainless steel Crown set is packaged in a unique, stackable tray that provides convenient, easy storage of 3M[™] ESPE[™] prefabricated crowns

Typical properties

Arch: Upper/Lower Brands: 3M™

Catalog Number: PO-96 Dentition: Permanent

Factory ISO Certification: No ISO Certification information available

Number of Crowns: 96
Patient: Pediatric
Size: Assorted

Tooth: Assorted

3M™ Paradigm™ DeepCure Light Guide, 76982, 1 - 10mm Light Guide

Part Number 76982 3M ID 7100057733 UPC 04035077010061

Details

Friction-fit, 360 degree rotating light guide mount The angle of the light guide is designed to make it more comfortable for patients even if they have limited mouth-opening abilities

Typical properties

Brands: Paradigm™

Factory ISO Certification: No ISO Certification information available

Package Contents: Includes one 10mm light guide

Product Type: Light Guide

3M™ Ketac™ Nano Glass Ionomer Primer, 3304P, 6.5ml

3M ID 7000128813 UPC 30605861036437

Details

Primer prior to restoration to maximize adhesion

Typical properties

Brands: Ketac™

Catalog Number: 3304P Cure Type: Light Cure Delivery System: Vial

Factory ISO Certification: No ISO Certification information available

Package Contents: Refill includes 1 - 65ml primer

Total Units per Case: 10 per Case

3M™ Ketac™ Cem EasyMix Triple Kit EE, 56906

3M ID 7100277467

Details

Granular formulation for easy powder dispensation without residual dust Releases fluoride

Excellent liquid absorption by the granules for great mixing High marginal integrity and low solubility for a long-lasting restoration High compressive strength and bond strength for lasting durability Thin film thickness for precise work without raised bite

Typical properties

Brands: Ketac™

Catalog Number: 56906 Delivery System: Bottle

Factory ISO Certification: No ISO Certification information available

3M™ Giz Accessory Pack Maxicap, 37180

3M ID 7100288681 UPC 04064035074098

Details

3M[™] Maxicap[™] activator: Easy activation of the powder and liquid to prepare the 3M[™] Maxicap[™] capsule for mixing

3M™ Maxicap™ applier: Easy delivery of 3M™ Maxicap™ capsule material to the tooth preparation or to the restoration

Easy delivery of 3M[™] Maxicap[™] capsule material to the tooth preparation or to the restoration

Typical properties

Brands: Maxicap™

Factory ISO Certification: No ISO Certification information available

Package Contents: Refill includes one 3M™ Maxicap™ Activator and one 3M™

Maxicap™ Applier

Product Type: Activator/Applier Set

3M™ Elipar™ DeepCure-S Light Guide, 76981, 1 - 10mm Guide

Part Number 76981 3M ID 7100057117 UPC 04035077010054

Details

Designed so that light intensity, despite focusing, is comparable to that of a standard light guide with a diameter of 10mm

One-shot cure with up to 50% more optical power and up to 80% more curing surface compared to leading devices tested

Faster, deeper cure

Simple two-button operation with one cure mode for all indications

Tack cure function allows for predictable curing with a reproducible short light output

Silent operation - innovative design eliminates noisy fan, and includes option to turn off beep signals

Two light guides, both offering consistently high and focused light output

Typical properties

Brands: Elipar™

Catalog Number: 76981

Factory ISO Certification: No ISO Certification information available

Package Contents: Includes on 10mm light guide

Product Type: Light Guide

3M™ Elipar™ S10 Eye Shield, 76984, 1 - Eye Shield

Part Number 76984 3M ID 7100057118 UPC 04035077010238

Details

Protects patient and user from reflection and intensive scattered light

Typical properties

Brands: Elipar™

Catalog Number: 76984

Factory ISO Certification: No ISO Certification information available

Package Contents: Includes one eye shield

Product Type: Eye Shield

3M[™] Vitrebond[™] Plus Light Cure Glass Ionomer Liner/Base Refill, 7570 Series

Details

Fast and easy dispensing and the right mix every time with the award-winning 3M™ Clicker™ dispenser

Conventional paste-liquid formula offers the protection of a conventional glass ionomer liner without the mess of measuring

Award-winning Clicker dispenser assures a consistent mix

Fluoride release Radiopaque

Fluoride release of a true glass ionomer

Excellent handling characteristics for easy placement

Cures in 20 seconds - a faster cure time than powder-liquid formulas

Delivers excellent adhesion with any dental adhesive system

Radiopaque

No additional hardware required

Typical properties

Brands: Vitrebond™

Catalog Number: 7572,7571

Cure Type: Light Cure

Delivery System: Clicker™ Dispenser

Factory ISO Certification: No ISO Certification information available

Inner Packs per Case: 10 Packs per Case

Package Contents: Refill includes 1 - 10g 3M[™] Clicker[™] Clicker Dispenser, 1 Large Mixing Pad and 1 Technique Card,Refill includes 2 - 10g 3M[™] Clicker[™]

Clicker Dispenser, 1 Large Mixing Pad and 1 Technique Card

3M™ Penta™ Plunger Discs, 77780

3M ID 7000145294 UPC 04035077012904

Details

Includes one pair of plastic plunger disks Plastic plunger discs Easy to use

Typical properties

Brands: Penta™

Catalog Number: 77780

Factory ISO Certification: No ISO Certification information available

Package Contents: 2 - Plunger Discs

Product Type: Plunger Disc

3M™ RelyX™ Unicem Elongation Tips for 3M™ Aplicap™ Capsules, 77550

3M ID 7100250304 UPC 04054596875708

Details

Use with 3M™ RelyX™ Universal or 3M™ RelyX™ Unicem/Unicem 2 resin cements for post cementation

Unique design allows for elongation tip to reach into the root canal from bottom to top

Virtually eliminates voids which can compromise bond strength Tip makes the entire procedure faster, easier and more reliable Strong and moisture-tolerant

**Except for 2- or 3-unit Maryland bridges and 3-unit inlay/onlay bridges. Selective etching is required to increase surface area for bonding

Typical properties

Brands: RelyX™

Catalog Number: 77550

Factory ISO Certification: No ISO Certification information available

Package Contents: 10 - Elongation Tips

Product Type: Tip

3M™ Filtek™ Universal Restorative Essentials Kit

Details

Four Designer Shades A1, A2, A3, A3.5, Plus a PO
NaturalMatch Technology (A blend of nanofillers and proprietary monomers)
Creates beautiful, natural restorations with less shades
Blends to surrounding dentition - Chameleon effect
Designer shades match the VITA classical and bleach shades*
Excellent, creamy handling
AFM Monomer provides stress relief, Low-shrinkage stress
Long-lasting polish
Designer shades cover most patients
Capsules can be warmed (See IFU)
2mm depth of cure, PO 1mm

3M™ Mini Uni-Twin™ Brackets

Details

Single piece, low profile ball hook increases patient comfort

Permanent I.D. markings make for faster, easier bracket identification

The tooth-shaped angulation makes for easier, more accurate positioning Sculptured labial surfaces reduce soft tissue irritation and enhance patient

comfort

The mesial/distal pockets increase wire deflectability while improving bracket engagement and reducing forces on the teeth

Twin tie-wings provide rotational control and allow selective ligation for rotation ties

Typical properties

Adhesive Precoat: APC™ II, APC™ PLUS, Non-coated

Angulation: +10°, +11°, +4°, +5°, +6°, +9°, 0°

Arch: Lower, Upper, Upper Left/Lower Right, Upper Right/Lower Left,

Upper/Lower

Base Type: Weldable Brands: Mini Uni-Twin™

Factory ISO Certification: No ISO Certification information available

Hook Included: Yes Hook Type: Distal

In/Out (Imperial): 0.03 in, 0.032 in, 0.033 in, 0.034 in, 0.035 in, 0.04 in, 0.044 in,

0.045 in. 0.046 in

In/Out (Metric): 0.76 mm, 0.81 mm, 0.84 mm, 0.86 mm, 0.89 mm, 1.02 mm, 1.12

mm, 1.14 mm, 1.17 mm Material: Stainless Steel

Mesial/Distal Width (Metric): 2.46 mm, 2.9 mm, 3.05 mm, 3.3 mm, 3.45 mm

Offset Base: Yes

Prescription: Alexander*, Andrews*, Dischinger*, Hilgers*, Level Anchorage™,

Roth*, Samson*, Sondhi, Standard Edgewise, ZeroBase™

Side: Left, Left/Right, Right

Slot Size (Imperial): 0.018 in, 0.022 in Slot Size (Metric): 0.457 mm, 0.559 mm

Tooth: Anterior, Bicuspid, Central Incisor, Cuspid, First Bicuspid, Lateral Incisor,

Second Bicuspid

Torque: +10°, +12°, +14°, +17°, +22°, +3°, +7°, +8°, -10°, -11°, -17°, -1°, -22°, -

5°, -7°, 0°

Units per Pack: 5

3M™ Typodonts

Details

Available in a variety of models

Typical properties

Bracket Style: Low Profile, Self-ligating, Twin, V-slot

Brands: 3M™

Factory ISO Certification: No ISO Certification information available eClass 14 Classification Group: 34110000

3M™ Unitek™ Stainless Steel Permanent Molar Crowns, 902350, Kit

3M ID 7000054516 UPC 30605861012516

Details

Permanent molar crown applications

Typical properties

Arch: Upper/Lower Brands: Unitek™

Catalog Number: 902350 Dentition: Permanent

Factory ISO Certification: No ISO Certification information available

Number of Crowns: 84

Patient: Pediatric Size: Assorted Tooth: Assorted

3M™ Clarity™ Advanced Ceramic Brackets Kit

Details

Small, low profile bracket enhances patient comfort without compromising strength, while minimizing occlusal interference

Typical properties

Adhesive Precoat: APC™ Flash-Free, Non-coated

Arch: Upper, Upper/Lower

Bracket Style: Twin Brands: Clarity™

Factory ISO Certification: No ISO Certification information available

Material: Ceramic

Prescription: High Torque, MBT[™], Roth* Slot Size (Imperial): 0.018 in, 0.022 in Slot Size (Metric): 0.457 mm, 0.559 mm Units per Pack: 1, 6, 10, 12, 20

3M™ Unitek™ Nitinol Closed Coil Spring

Details

Length 9 mm and 12 mm lengths B Dimension, inner diameter (arbor) 0.030" (0.76 mm) Packaged 10 per box

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Force Rating: Heavy, Light, Medium

Material: Nitinol

Units per Inner Pack: 10 per Pack

eClass 14 Classification Group: 34110690

Dimensions and Classifications

Overall Length (Imperial): 0.354 in, 0.472 in Overall Length (Metric): 9 mm, 12 mm

3M[™] Transbond[™] MIP Moisture Insensitive Primer, 712-025, 6 ml, 1/Btl

3M ID 7000004378 UPC 00652221108883

Details

Moisture tolerant without compromising bond strength Compatible with all 3M light cure adhesives Quick and easy application with just one coat needed Helps to reduce chair time and increase office efficiency

Typical properties

Brands: Transbond™ Catalog Number: 712-025

Container Size (Imperial): 0.203 fl oz

Container Size (Metric): 6 mL

Delivery System: Bottle

Factory ISO Certification: No ISO Certification information available

Product Type: Moisture Insensitive Primer Kit

Units per Inner Pack: 1 /BT

3M™ Unitek™ Ball Clasp

Details

Stainless steel round wire with integrated round ball at one end Ball diameter approximately twice wire diameter

Typical properties

Archwire Size (Imperial): 0.024 in, 0.028 in, 0.032 in, 0.04 in Archwire Size (Metric): 0.61 mm, 0.711 mm, 0.813 mm, 1.016 mm

Brands: Unitek ™

Factory ISO Certification: No ISO Certification information available

Units per Inner Pack: 10 per Pack

3M™ AlastiK™ Bobbin Chain Dispenser, 406-927, 1/Ea

3M ID 7000025418

Details

Not made with natural rubber latex Holds 16 bobbins Removable handle Replaceable blade

Typical properties

Brands: AlastiK™

Factory ISO Certification: No ISO Certification information available

Units per Inner Pack: 1 per Pack

3M™ AlastiK™ Lingual Ligatures, 406-993, Navy Blue, 12/Stick, 12/Pk

3M ID 7100154423

Details

3M[™] AlastiK[™] Lingual Ligatures, when used as plain ties, apply similar tip control when compared to overties used with powerchain (Source: 3M Lab) The tab designed on the 3M[™] AlastiK[™] Lingual Ligatures ensures a safe grip and simplifies placement and handling during the ligation procedure The ligatures have been shown to break less than commonly used chains during installation and between appointments and they provide patient comfort at the same time

Optimized cross section delivers ideal forces

Easy to see on gold brackets

Innovative design for easy ligating

Proprietary design for easy application

Replaces overties

Recommended by 90% of doctors

Convenient tab enables secure grip

Typical properties

Brands: AlastiK™

Factory ISO Certification: No ISO Certification information available

Product Color: Navy Blue

Units per Inner Pack: 12 per Pack

3M™ Unitek™ Nitinol Super-Elastic Archwire

Details

Provides light continuous forces Easy engagement characteristics Good early torque control

Typical properties

Arch: Lower, Upper

Arch Shape: LA, Ovoid, Reverse Curve, Square

Archwire Size (Imperial): 0.014 in, 0.014 x 0.025 in, 0.016 in, 0.016 x 0.016 in, 0.016 x 0.022 in, 0.016 x 0.025 in, 0.017 x 0.025 in, 0.018 in, 0.018 x 0.018 in,

0.018 x 0.025 in, 0.019 x 0.025 in, 0.02 in, 0.021 x 0.025 in

Archwire Size (Metric): 0.356 mm, 0.356 x 0.635 mm, 0.406 mm, 0.406 x 0.406 mm, 0.406 x 0.559 mm, 0.406 x 0.635 mm, 0.432 x 0.635 mm, 0.457 mm, 0.457 x 0.635 mm, 0.457 x 0.635 mm, 0.508 mm, 0.533 x 0.635 mm

Brands: Unitek™

Cross Section Shape: Hybrid

Dimpled: Yes

Factory ISO Certification: No ISO Certification information available

Material: Nitinol

Units per Inner Pack: 10 per Pack

3M[™] APC[™] Inventory Dispensing System, 2600-022, Drawer Dust Cover Set, 1/Ea

3M ID 7100286981

Details

Helps protect APC blisters in drawer Clear protective cover Includes 5X5 upper/lower Palmer Notation labels to identify the contents of the drawers

Typical properties

Brands: 3M™

Factory ISO Certification: No ISO Certification information available

3M™ Unitek™ Prestige How Utility Pliers

Details

Matching serrated tips Grip surfaces securely hold wire or auxiliary at any convenient working angle 3M™ Unitek™ Prestige Series

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

3M™ Unitek™ Prestige How Utility Pliers, Os Tip, 900-702, 1/Ea

3M ID 7000004392

Details

For rectangular and square wires
Gripping surfaces are precision machined to hold wire at right angle without slipping
3M™ Unitek™ Prestige Series

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

3M™ Unitek™ Raised Arm Lingual Cleats

Details

Accommodates up to three heavy elastics on one cleat

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

3M™ OrthoForm™ Operatory Arch Form Templates, White, 701-724, Ea

3M ID 7100241940

Details

Designed to aid in arch form symmetry determination 3/pack, 1 each OrthoForm I, OrthoForm II, OrthoForm III

Typical properties

Brands: OrthoForm™

Factory ISO Certification: No ISO Certification information available

3M™ Unitek™ Gurin Loop Spring

Details

Anchor without welding or soldering Adjustable

Typical properties

Brands: Unitek™

Compatible Archwire Size (Imperial): 0.018 - 0.022 in

Factory ISO Certification: No ISO Certification information available

Side: Left, Right

Units per Inner Pack: 1 per Pack

3M™ Unitek™ Gurin Wrench

Details

Machined from stainless steel Has a knurled grasp

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

3M™ Unitek™ Clarity™ Ultra SL Open-Close Instrument, 804-155, 1/Ea

3M ID 7000001409 UPC 00652221110121

Details

1/each

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

3M™ Unitek™ Nitinol Classic Archwire

Details

Long working range 40% the forces of stainless steel

Typical properties

Arch: Lower, Upper

Arch Shape: LA, Ovoid, Square, Standard

Archwire Size (Imperial): 0.012 in, 0.014 in, 0.014 x 0.025 in, 0.016 in, 0.016 x 0.016 in, 0.016 x 0.022 in, 0.017 x 0.022 in, 0.017 x 0.025 in, 0.018 in, 0.018 x 0.018 x 0.018 x 0.018 x 0.018 x 0.018 x 0.018 in, 0.018 x 0.025 in, 0.018 x 0.025 in

Archwire Size (Metric): 0.305 mm, 0.356 mm, 0.356 x 0.635 mm, 0.406 mm, 0.406 x 0.406 mm, 0.406 x 0.559 mm, 0.432 x 0.559 mm, 0.432 x 0.635 mm, 0.457 mm, 0.457 x 0.457 mm, 0.457 x 0.635 mm, 0.457 x 0.635 mm, 0.538 mm, 0.538 mm, 0.538 mm, 0.538 mm, 0.538 mm, 0.538 mm

Brands: Unitek™

Cross Section Shape: Hybrid, Rectangular

Dimpled: Yes

Factory ISO Certification: No ISO Certification information available

Material: Nitinol

Units per Inner Pack: 10 per Pack

3M™ Victory Series™ Low Profile Brackets Kit

Details

Uncompromised mesial-distal dimension

Micro-etch mesh base surface with three point contact design and anatomical contouring

Vertical scribe line for accurate bracket placement

Uncompromising under tie-wing area for ease of ligation

Smooth radius surface

Color-coded quadrant ID system for bracket identification

Typical properties

Adhesive Precoat: APC™ Flash-Free, APC™ II, Non-coated

Arch: Upper/Lower Bracket Style: Twin Brands: Victory Series™

Factory ISO Certification: No ISO Certification information available

Hook Included: Yes Hook Type: Distal

Prescription: MBT™, Roth*

Slot Size (Imperial): 0.018 in, 0.022 in Slot Size (Metric): 0.457 mm, 0.559 mm

Units per Pack: 20, 24

3M™ Unitek™ Prestige Arch Bending Pliers, 900-721, 1/Ea

3M ID 7000025676

Details

"Sure-grip" precision-matched serrated pads
Outer surface of beaks are carefully rounded for safety and comfort to lips and
intraoral tissue
Smaller tips for better intraoral accessibility

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Pliers Type: Weingart

3M™ RelyX™ Ultimate Adhesive Resin Cement Refill

Details

High bond strength

Easy-to-use, two-component system

Versatile for total-etch, selective etch or self-etch techniques

Excellent color stability with natural fluorescence

Offers wear resistance

If you like 3M™ RelyX™ Ultimate Cement, you'll love 3M™ RelyX™ Universal

Resin Cement

Typical properties

Brands: RelyX™

Catalog Number: 56887,56920,56888,56921,56890,56889,56923

Cure Type: Dual Cure

Delivery System: Automix Syringe

Factory ISO Certification: No ISO Certification information available

Material: Adhesive Resin

Package Contents: 1 – 8.5g Syringe; 10 Mixing Tips (regular); 5 Mixing Tips

(wide); 5 Intra-oral Tips Shade: A1, A3, B0.5, TR

3M[™] Ortholux[™] Luminous Curing Light Rechargeable Battery, 704-455, 1/Ea

3M ID 7000148601 UPC 00652221149367

Details

Lithium ion battery Twist insert into handpiece

Typical properties

Brands: Ortholux™

Factory ISO Certification: No ISO Certification information available

Product Type: Battery

3M™ Filtek™ Bulk Fill Flowable Restorative, Capsule

Details

Flowable composite resin filling with 4mm depth of cure that reduces the need for incremental placement and risk of contamination

Flowable viscosity results in easy handling and adaptation

Blend of fillers and monomers maximizes strength, wear resistance and radiopacity

Complete your premium restoration using the gold standard of adhesion with 3M™ Scotchbond™ Universal Plus Adhesive

Typical properties

Brands: Filtek™

Catalog Number: 4861U,4861A1,4861A2,4861A3

Composite Type: Bulk Fill Flowable

Cure Type: Light Cure

Delivery System: Capsule, Syringe Environmental claim being made: Yes Environmental claim relates to: Product

Factory ISO Certification: No ISO Certification information available

Material Chemistry: Zirconia & Silica Filler Material Structure: Microhybrid Composite

Net Weight (Metric): 0.2 g Shade: A1, A2, A3, Universal

Units per Inner Pack: 15 per Pack, 20 per Pack

Units per Pack: 15, 20

3M™ Polyether Contact Tray Adhesive Refill, 69408

3M ID 7100196387 UPC 04054596648791

Details

Easy to apply Good adhesion after short one minute drying time Easy to remove

Typical properties

Brands: 3M™

Catalog Number: 69408

Container Size (Imperial): 0.57 oz Container Size (Metric): 17 mL

Cure Type: Air Dry Delivery System: Bottle

Factory ISO Certification: No ISO Certification information available

Material: Polyether

Package Contents: 1 refill - 17ml Product Type: Tray Adhesive

3M™ Chairside Zirconia for CEREC® Kit, 69956

3M ID 7100189194 UPC 04054596631267

Details

3M[™] Chairside Zirconia Intro Kit includes: 1- crown block A1, 2- crown blocks A2, 1- crown block A3, 1- bridge block A2, 1- Instructions for Use New Zirconia block for single-unit crown and three-unit bridge restorations* Optimal blend of high strength and esthetics, with high flexural strength of 1000 MPa (fast fired)

The esthetics you desire from a chairside zirconia, designed to match the VITA classical shade guide

Fast sintering time of approximately 20 minutes for a thin walled crown in CEREC® SpeedFire furnace**,***

Makes it possible to go from scan to seat in one hour Available in two pre-shaded block sizes and eight shades

Typical properties

Brands: 3M™

CAD/CAM System: CEREC® Catalog Number: 69956

Digital Materials: Digital Chairside Restorative Materials

Factory ISO Certification: No ISO Certification information available

Material: Zirconia

Shade: Multiple Shades

Dimensions and Classifications

Overall Diameter (Imperial): 0.63 in Overall Diameter (Metric): 16 mm

3M™ Single Bond Universal Adhesive Vial Refill

Details

Bonds dental restorations with uncompromising strength with a one-step, one-coat application

Compatible with all etching techniques

Compatible with all resin cements

Self- and dual-cure compatible with 3M™ Single Bond Universal DCA Dual-Cure Activator

Confidently bonds to both moist and dry dentin

Offers a strong bond to all indirect surfaces without an additional primer for more streamlined workflows

Results in virtually no post-operative sensitivity,* creating a better patient experience. *3M Internal data

Available in a flip-top vial for one-handed operation or unit dose delivery

Typical properties

Adhesive Type: Self-etch, Total Etch

Brands: 3M™

Catalog Number: 41266,41269,41282

Container Size (Imperial): 0.101 fl oz, 0.169 fl oz

Container Size (Metric): 3 mL, 5 mL

Cure Type: Light Cure

Delivery System: Bottle, Vial

Factory ISO Certification: No ISO Certification information available

Product Type: Dental Adhesive Units per Inner Pack: 1 per Pack

3M™ Victory Series™ Mini Tubes

Details

Available with 3M[™] APC[™] Adhesive Small base designed for placement on partially erupted second molars Seating notch designed for easy placement

Typical properties

Adhesive Precoat: APC™ Flash-Free, APC™ II, APC™ PLUS, Non-coated

Arch: Lower, Upper Base Type: Bondable Brands: Victory Series™

Catalog Number: 4066-4082,4066-4071,4066-5033,4066-5044,4066-4033,4066-

4044,4066-5082,4066-5071 Distal Offset: 0 Degree

Factory ISO Certification: No ISO Certification information available

Hook Included: Yes Hook Type: Distal

Material: Stainless Steel

Prescription: Low Torque, MBT™

Side: Left, Right

Slot Size (Imperial): 0.018 in, 0.022 in Slot Size (Metric): 0.457 mm, 0.559 mm

Tooth: Second Molar Torque: -10°, -14°, -25°, 0°

Tube Style: Single Non-convertible

Units per Inner Pack: 4 per Pack, 5 per Pack eClass 14 Classification Group: 34110000

Dimensions and Classifications

Overall Length (Imperial): 0.098 in Overall Length (Metric): 2.5 mm

3M™ AlastiK™ QuiK-StiK™ Ligatures

Details

Not made with natural rubber latex A-1 or A-2 modules 60 modules per stick Central thumb ring

Typical properties

Brands: AlastiK™

Factory ISO Certification: No ISO Certification information available

Ligature Size: A-1, A-2

Product Color: Assorted, Blue Metallic, Deep Black, Dusky Blue, Extreme Green,

Gold, Gray, Green Metallic, Hot Red, Kelly Green, Lavender, Light Blue,

Mediterranean Blue, Midnight Glow, Navy Blue, Obscure, Orange Blast, Plum Crazy, Purple Metallic, Radical Red, Screaming Pink, Silver, Smoke, Teal,

Transparent, White, Yellow

Units per Inner Pack: 17 per Pack

eClass 14 Classification Group: 34110603

Dimensions and Classifications

Overall Diameter (Imperial): 0.094 in, 0.125 in Overall Diameter (Metric): 2.4 mm, 3.2 mm

3M™ APC™ Inventory Dispensing System, Back-Up Storage Units, 7 x 7, 2600-032, 1/Ea

3M ID 70006629433

Details

Allows first-in/first-out inventory dispensing
It can accommodate a 7x7 appliance system
Contains complete set of labeling for drawer content identification

Typical properties

Brands: 3M™

Factory ISO Certification: No ISO Certification information available

3M[™] Concise[™] Orthodontic Chemical Cure Adhesive Enamel Bond Resin B, 196-010, 6 ml, 1/Btl

3M ID 7000024652 UPC 00652221047892

Details

One part of two-part resin and paste system Cures chemically Reliable bond strength

Typical properties

Adhesive Type: Direct Brands: Concise™

Catalog Number: 196-010

Container Size (Imperial): 0.203 fl oz

Container Size (Metric): 6 mL Cure Type: Chemical Cure

Factory ISO Certification: No ISO Certification information available

Product Type: Orthodontic Primer

Units per Inner Pack: 1 /BT

3M™ Adper™ Single Bond Plus Adhesive Refill, 51102

3M ID 7000128797 UPC 30605861036093

Details

Offers good bond strength and performance, providing your patients with lasting dental restorations

Convenient bottle design offers easier opening and dispensing

No shaking required thanks to stable nano-filler designed not to settle within the bottle

Reduces the risk of post-operative sensitivity, creating a better experience for your patients

Compatible with light-cured composite restoratives

Typical properties

Adhesive Type: Total Etch

Brands: Adper™

Catalog Number: 51102 Cure Type: Light Cure Delivery System: Vial

Factory ISO Certification: No ISO Certification information available

Package Contents: Provide what is included in the refill package contents

Product Type: Dental Adhesive

3M™ Unitek™ Preformed Vertical Lingual Archwire with Loop

Details

.036 in. (.91 mm) diameter medium-hard Permachrome wire Arches feature two-point contact with sheath by means of two distal posts – one double, one single

Typical properties

Arch Length (Imperial): 2.05 in, 2.12 in, 2.2 in, 2.28 in, 2.36 in, 2.44 in, 2.5 in, 2.6

in, 2.68 in, 2.75 in

Arch Length (Metric): 52.07 mm, 53.848 mm, 55.88 mm, 57.912 mm, 59.944 mm,

61.976 mm, 63.5 mm, 66.04 mm, 68.072 mm, 69.85 mm

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Material: Stainless Steel

Units per Inner Pack: 10 per Pack Wire Diameter (Imperial): 0.036 in Wire Diameter (Metric): 0.9144 mm

3M™ Filtek™ Supreme Flowable Restorative, Syringe

Details

Offers excellent adaptation and versatility for a wide range of indications and techniques, including injection molding

Upgraded syringe creates virtually no bubbles or material "run-on"

Bendable cannula improves access to deep restorations

Ergonomic design is easy to hold and requires less effort to inject

Made with 3M's TRUE nanotechnology for excellent polish retention and wear resistance

Available in a wide variety of shades to match the VITA Classical Shade Guide and body shades of 3M[™] Filtek[™] Supreme Ultra Universal Restorative for natural-looking esthetics

Can be warmed up to 70°C/158°F for up to 1 hour, up to 25 times

Complete your premium restoration using the gold standard of adhesion with 3M™ Scotchbond™ Universal Plus Adhesive

Typical properties

Brands: Filtek™ Catalog Number:

3930XW,6032A2,6032A3,6032C2,6032A1,6032A4,3930C2,6032W,3930A1,3930

A2,3930A3,7032A3,7032A1,7032A2,6032A2-S,7032A3.5,6032A3-

S,6032A3.5,6032D2,6032B1,6032B2,7032B1,3930A3.5,6032OA3,6032XW

Composite Type: Flowable Cure Type: Light Cure Delivery System: Syringe

Environmental claim being made: Yes Environmental claim relates to: Product Factory ISO Certification: ISO 13485 Inner Packs per Case: 10 Packs per Case Material Chemistry: Zirconia & Silica Filler

Material Structure: Nanocomposite

Net Weight (Metric): 1 g, 2 g

Opacity: Body, Extra White, Opaque

Package Contents: Refill includes 2 - 2g syringes in shade C2; 20 dispensing tips,Refill includes 2 - 2g syringes in shade A3; 20 dispensing tips,Refill includes 2 - 2g syringes in shade B1; 20 dispensing tips,Refill includes 2 - 2g syringes in shade XW; 20 dispensing tips,Refill includes 2 - 2g syringes in shade XW; 20 dispensing tips,Refill includes 2 - 2g syringes in shade A3.5; 20 dispensing tips,Refill includes 2 - 2g syringes in shade D2; 20 dispensing tips,Refill includes 2 - 2g syringes in shade A2; 20 dispensing tips,Refill includes 2 - 2g syringes in shade A1; 20 dispensing tips,Refill includes 2 - 2g syringes in shade A4; 20 dispensing tips,Refill includes 2 - 2g syringes in shade A4; 20 dispensing tips,Refill includes 2 - 2g syringes in shade B2; 20 dispensing tips

Product Form: Refill, Sample

Shade: A1, A2, A3, A3.5, A4, B1, B2, C2, D2, OA3, W, XW

Total Units per Case: 200 per Case

Units per Inner Pack: 1 per Pack, 2 per Pack

Units per Pack: 1, 2

3M™ AlastiK™ C Modules

Details

Not made with natural rubber latex Injection molded Three configurations

Typical properties

Brands: AlastiK™

Center to Center Distance (Imperial): 0.115 in, 0.15 in, 0.17 in, 0.175 in Factory ISO Certification: No ISO Certification information available

Force Rating: Heavy, Medium Ligature Size: C-1, C-2, CK Product Color: Gray, Transparent

Units per Inner Pack: 10 per Pack

3M™ AlastiK™ Separator Modules, 406-080 Series

Details

Not made with natural rubber latex Three sizes 100 modules/pack

Typical properties

Brands: AlastiK™

Factory ISO Certification: No ISO Certification information available

Product Color: Blue, Green, Transparent eClass 14 Classification Group: 34110603

Dimensions and Classifications

Overall Diameter (Imperial): 0.14 in, 0.175 in, 0.21 in Overall Diameter (Metric): 0.14 mm, 0.21 mm, 4.4 mm

3M™ Unitek™ Nitinol Open Coil Spring

Details

Length 7" (175 mm) Packaged 3 per tube

Typical properties

Arbor Diameter (Imperial): 0.03 in, 0.033 in

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Force Rating: Heavy, Light, Medium

Material: Nitinol

Units per Inner Pack: 3 per Pack

eClass 14 Classification Group: 34110690

Dimensions and Classifications

Overall Length (Imperial): 7 in Overall Length (Metric): 175 mm

3M™ Unitek™ Bracket Placement Instrument

Details

Stainless steel reverse action tweezer

Downward slope of the forearm improves operator visibility while placing the bracket onto the tooth

Constructed with a convex reverse end that can be used to seat brackets

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

3M[™] Unitek[™] Lingual Weingart Utility Pliers, 803-530, 60 Degree Bend, 1/Ea

3M ID 7100154550

Details

Used to make a 90° wire bend distal to the buccal tube Long handles assure convenient access to the area with patient comfort Serrated tips hold the wire firmly

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

3M™ Unitek™ Double-Ended Ligature Director, 803-530, 1/Ea

3M ID 7100266449

Details

Tips angled at 60° for easy access to the lingual arch Fine tapered tips help make working in the limited space easier Tip serrations provide a secure grip on the wire Tips are rounded for patient comfort and safety

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

3M™ Unitek™ Bondable Lingual Cleat, 709-050, 10/Pk

3M ID 7000004377

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

3M™ Wax-Pak™ Sticks, 603-063, 100/Pk

3M ID 7100154549

Details

Tasteless
Packaged conveniently for the patient to carry

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

3M™ Unitek™ Retraction Spring

Details

.188 in. (4.8 mm) closewound spring Intermaxillary retraction

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Force (Imperial): 8.0 ozf,4.0 ozf

Force Rating: Light

Units per Inner Pack: 10 per Pack

3M™ Scotchbond™ Universal Plus Adhesive Kits

Details

First radiopaque universal adhesive lets you provide reliable bonding while reducing the risk of X-ray misdiagnosis and overtreatment

Offers gold standard adhesion to all dental substrates without separate primers, silanes or activators

Bonds and seals caries-affected dentin to support minimally invasive dentistry Full dual and self-cure compatibility, without a separate activator BPA derivative-free formulation

Completes a fully aligned, two-component system with 3M™ RelyX™ Universal Resin Cement

Results in virtually no post-operative sensitivity,* helping provide a more positive experience for your patients. *3M Internal data

Suitable for all direct and indirect bonding indications and all etching techniques Available in a streamlined flip-top vial that provides up to 200 applications or unit dose for efficient hygiene management

Typical properties

Adhesive Type: Selective Etch

Brands: Scotchbond™

Catalog Number: 41297,41293

Cure Type: Light Cure

Delivery System: Unit Dose, Vial

Factory ISO Certification: No ISO Certification information available

Package Contents: 50 L-Pop,1x 3M[™] Scotchbond[™] Universal Plus Adhesive vial (5 ml), 50 disposable brushes, 1x 3M[™] Scotchbond[™] Universal Etchant

(3ml), 25 dispensing tips

Product Type: Dental Adhesive

3M™ Sof-Lex™ Mandrel Refill

Details

Low-profile mandrel design reduces risk of marring the restoration for easier use and improved results
Includes 3 replacement RA mandrels
Easier to use with better results

Typical properties

Brands: Sof-Lex™

Factory ISO Certification: No ISO Certification information available

For Use With: Contouring & Polishing Disc

Package Contents: Refill includes three RA mandrels, Refill includes three HP

mandrels, Refill includes three FG mandrels

3M™ XS Applicator Brushes, 41904, 50 ct

3M ID 7000054649 UPC 04035077012713

Details

Extra-small applicator works well with very small preparations and in hard-toreach areas of the mouth Single-use, disposable brush for hygiene management

Typical properties

Brands: 3M™

Factory ISO Certification: No ISO Certification information available

Package Contents: Refill includes 50 adhesive single-use XS applicator brushes

Product Type: Applicator Brush

3M™ Delivery Tips With Pistons, 3303DT, 50ct

3M ID 7100213641 UPC 30605861068797

Details

Allows for easy delivery of 3M™ Vitremer™ core buildup/restorative

Typical properties

Brands: 3M™

Factory ISO Certification: No ISO Certification information available

Package Contents: Refill includes 50 delivery tips

Product Type: Tip

3M™ Adper™ Scotchbond™ Multi-Purpose 3 Well Dispenser, 7544

3M ID 7000003064 UPC 30605861032378

Details

Reusable Protects adhesive from light exposure

Typical properties

Brands: Adper™

Factory ISO Certification: No ISO Certification information available

Number of Wells: 3

Package Contents: Refill includes 4 dispensers eClass 14 Classification Group: 34110000

3M™ Elastomer Syringe Knurled Hub Cap Refill, 71043

3M ID 7000054666 UPC 04035077012966

Details

Not compatible with the 3M[™] Penta[™] Elastomer Syringe

Typical properties

Brands: 3M™

Catalog Number: 71043

Factory ISO Certification: No ISO Certification information available

Package Contents: 1 - Syringe hub cap

Product Type: Knurled Hub Cap

3M™ Strip Crown Form - Pedo, 915100, Kit

3M ID 7000002613 UPC 30605861011090

Details

Automatically contours restorative material to match natural dentition
Strips off easily leaving a smooth surface
Thin interproximal walls allow for tight contacts
Crown set is packaged in a unique, stackable tray that provides convenient, easy storage

Typical properties

Arch: Upper/Lower Brands: 3M™

Catalog Number: 915100

Dentition: Primary

Factory ISO Certification: No ISO Certification information available

Number of Crowns: 120

Side: Left/Right Size: Assorted

Tooth: Central Incisor Units per Case: 120.0

3M™ Unitek™ Etching Gel System (Blue Tint) Kit

Details

Tips available in metal or plastic

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Package Contents: 20 Disposable Metal Dispensing Tips,25 Disposable Plastic

Dispensing Tips

3M[™] SmartClip[™] Appliance Wire Disengagement Hand Instrument Slide Reorder, 804-170, 5/Pk

3M ID 7000001448

Details

Deponding tool designed to work with many 3M brackets May be used with or without the archwire engaged in the bracket slot

Typical properties

Brands: SmartClip™

Factory ISO Certification: No ISO Certification information available

3M™ Unitek™ Stainless Steel Primary Anterior Crowns, 908100, Kit

3M ID 7000053486 UPC 30605861011175

Details

Primary central, laterals, cuspids

Typical properties

Arch: Upper/Lower Brands: Unitek™

Catalog Number: 908100

Dentition: Primary

Factory ISO Certification: No ISO Certification information available

Number of Crowns: 72 Patient: Pediatric Size: Assorted

Tooth: Assorted, Central Incisor

3M™ AlastiK™ Mini-StiK Patient Sticks

Details

Not made with natural rubber latex A-1 or A-3 modules 24 modules per stik Injection molded for smooth, comfortable surface Exciting colors

Typical properties

Brands: AlastiK™ Dispenser Type: Stick

Factory ISO Certification: No ISO Certification information available

Ligature Size: A-1, A-3

Product Color: Assorted, Blue Metallic, Deep Black, Dusky Blue, Extreme Green,

Gold, Gray, Green Metallic, Hot Red, Kelly Green, Lavender, Light Blue, Mediterranean Blue, Midnight Glow, Navy Blue, Obscure, Orange Blast, Plum Crazy, Purple Metallic, Radical Red, Screaming Pink, Silver, Smoke, Teal,

Transparent, White, Yellow

Units per Inner Pack: 42 per Pack, 417 per Pack eClass 14 Classification Group: 34110603

3M™ Transbond™ Plus Self Etching Primer, 712-093, Easy Roller, 1/Ea

3M ID 7100154494

Details

Facilitates activation and mixing of the two components in the Transbond Plus Self Etching Primer
Can be disinfected

Typical properties

Brands: Transbond™

Factory ISO Certification: No ISO Certification information available

3M™ AlastiK™ Rotation Wedge

Details

Not made with natural rubber latex Preformed tabs snap over tie-wings

Typical properties

Brands: AlastiK™

Factory ISO Certification: No ISO Certification information available

Product Color: Gray

Units per Inner Pack: 10 per Pack, 100 per Pack eClass 14 Classification Group: 34110690

3M™ AlastiK™ Loose Modules

Details

Not made with natural rubber latex Available in 2 colors Available in 2 force levels 1,000/pack

Typical properties

Brands: AlastiK™

Factory ISO Certification: No ISO Certification information available

Force Rating: Heavy, Medium

Ligature Size: A-1

Product Color: Gray, Transparent Units per Inner Pack: 100 per Pack

3M™ APC™ Reusable 7x7 Set-Up Tray, 2600-005, 1/Ea

3M ID 7100256146

Details

Allows visual identification of the contents prior to bonding Sterilizable Stackable

Typical properties

Brands: 3M™

Factory ISO Certification: No ISO Certification information available

3M™ Unitek™ Lingual Cinch Back Pliers, 803-531, 1/Ea

3M ID 7100266450

Details

With an enhanced gripping surface for security and a functional ergonomic design, the debonding pliers make the bracket debonding step easier Pliers accommodate the custom bases of the 3M™ Incognito™ Appliance System Brackets

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

3M™ Unitek™ Molar Band Seater, 811-001, 1/Ea

3M ID 7000025656

Details

"Sure-grip" precision-matched serrated pads
Outer surface of beaks are carefully rounded for safety and comfort to lips and
intraoral tissue
3M™ Unitek™ Prestige Series

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Pliers Type: Weingart

3M™ Unitek™ Low Profile Lingual Cleats

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Product Type: Lingual Cleat eClass 14 Classification Group: 34110603

3M™ Unitek™ Level Anchorage Palatal Bar

Details

Increases anchorage

Typical properties

Arch Size: 2, 3, 4 Brands: Unitek™

Factory ISO Certification: No ISO Certification information available Units per Inner Pack: 1 per Pack eClass 14 Classification Group: 34110602

3M™ Unitek™ Hi-T II Stainless Steel Coil Retraction Spring

Details

3 ft. (.91 mm) self-sealing bobbin Stainless steel spring

Typical properties

Arbor Diameter (Imperial): 0.02 in, 0.022 in, 0.03 in, 0.032 in, 0.036 in, 0.04 in,

0.3 in

Arbor Diameter (Metric): 0.559 mm, 0.813 mm, 0.914 mm, 1.016 mm, 7.62 mm

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Material: Stainless Steel

Units per Inner Pack: 1 per Pack

Wire Diameter (Imperial): 0.006 in, 0.007 in, 0.008 in, 0.009 in, 0.01 in

Wire Diameter (Metric): 0.152 mm, 0.178 mm, 0.203 mm, 0.229 mm, 0.254 mm

eClass 14 Classification Group: 34110690

Dimensions and Classifications

Overall Length (Imperial): 0.022 in, 0.03 in, 0.032 in, 0.04 in

Overall Length (Metric): 0.5588 mm, 0.762 mm, 0.8128 mm, 1.016 mm

3M[™] SmartClip[™] Appliance Wire Disengagement Hand Instrument, 804-162, 1/Ea

3M ID 7000024328

Details

1 kit

Typical properties

Brands: SmartClip™

Factory ISO Certification: No ISO Certification information available

3M™ Scotchbond™ Universal Plus Adhesive Unit Doses L-Pop

Details

First radiopaque universal adhesive lets you provide reliable bonding while reducing the risk of X-ray misdiagnosis and overtreatment

Offers gold standard adhesion to all dental substrates without separate primers, silanes or activators

Bonds and seals caries-affected dentin to support minimally invasive dentistry Full dual and self-cure compatibility, without a separate activator BPA derivative-free formulation

Completes a fully aligned, two-component system with 3M™ RelyX™ Universal Resin Cement

Results in virtually no post-operative sensitivity,* helping provide a more positive experience for your patients. *3M Internal data

Suitable for all direct and indirect bonding indications and all etching techniques Available in a streamlined flip-top vial that provides up to 200 applications or unit dose for efficient hygiene management

Typical properties

Adhesive Type: Universal Adhesive

Brands: Scotchbond™

Catalog Number: 41295,41294,41299,41298,41304

Container Size (Imperial): 0.169 fl oz

Container Size (Metric): 5 mL

Cure Type: Light Cure

Delivery System: Unit Dose, Vial

Factory ISO Certification: No ISO Certification information available Package Contents: 200 L-Pop,1 x 5ml,3 x 5ml,100 L-Pop,400 L-Pop

Product Type: Dental Adhesive

Units per Inner Pack: 1 per Pack, 3 per Pack, 100 per Pack, 200 per Pack, 400

per Pack

3M™ Unitek™ Multi-Cure Glass Ionomer Band Cement Liquid, 712-052, 25 g, 1/Btl

3M ID 7000004386

Details

Moisture and saliva tolerant
Fills voids between tooth surface and bands
Light cure or self cure
Blue tint for easy initial clean-up
Achieve adequate initial strength with light application
Band under moist conditions

Typical properties

Brands: Unitek™

Catalog Number: 712-052 Delivery System: Bottle

Factory ISO Certification: No ISO Certification information available

Units per Inner Pack: 1 /BT

3M[™] Just For Kids[™] 0.4% Stannous Fluoride Brush-On Gel Refill, 12107

Details

Formulated with stannous fluoride to help your patients prevent tooth decay

Brush-on therapy designed for young, caries-prone patients

Provides antimicrobial activity for up to 8 hours

Promotes tooth remineralization

Helps inhibit demineralization for lasting oral health

Offered in Grapey Grape and Bubble Gum flavors

Typical properties

Active Ingredient: Stannous Fluoride

Brands: Just For Kids™

Delivery System: Bottle, Tube

Factory ISO Certification: No ISO Certification information available

Features: Anti-cavity, Hypersensitivity

Flavor: Bubble Gum, Fruit Punch, Grapey Grape

Fluoride Content: 0.4 % Material: Stannous Fluoride Net Weight (Imperial): 4.3 oz Net Weight (Metric): 121.903 g

Product Form: Gel

Product Type: Brush-on Gel

Units per Inner Pack: 1 per Pack, 24 per Pack eClass 14 Classification Group: 25092006

3M[™] Sof-Lex[™] XT Pop-on Extra Thin Disc Refills - Orange Series, 30 x 12.7mm, Coarse, 4931C

3M ID 7000003147 UPC 30605861034853

Details

Easy to use system allows you to create high-luster restorations that offer lifelike esthetics

Discs are color coded to allow for easy identification of grit sequence Pop-on disc technology makes changing discs easier to save you time Pop-on discs conform to the contours of the teeth for easier finishing and polishing

Typical properties

Brands: Sof-Lex™

Factory ISO Certification: No ISO Certification information available

Grade: Coarse Model: XT

Product Color: Orange

eClass 14 Classification Group: 34111407

Dimensions and Classifications

Overall Diameter (Imperial): 0.5 in Overall Diameter (Metric): 12.7 mm

3M[™] Sof-Lex[™] XT Pop-on Extra Thin Disc Refills - Orange Series, 30 x 12.7mm, Superfine, 4931SF

3M ID 7100006357 UPC 30605861034914

Details

Easy to use system allows you to create high-luster restorations that offer lifelike esthetics

Discs are color coded to allow for easy identification of grit sequence Pop-on disc technology makes changing discs easier to save you time Pop-on discs conform to the contours of the teeth for easier finishing and polishing

Typical properties

Brands: Sof-Lex™

Factory ISO Certification: No ISO Certification information available

Grade: Super Fine

Model: XT

Product Color: Orange

eClass 14 Classification Group: 34111407

Dimensions and Classifications

Overall Diameter (Imperial): 0.5 in Overall Diameter (Metric): 12.7 mm

3M™ Unitek™ Miniature Twin V-Slot Hook Pins, 461-140, 100/Pk

3M ID 70202082544

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Units per Pack: 100

3M[™] Concise[™] Orthodontic Chemical Cure Adhesive Bonding System, 196-001, 1/Kit

3M ID 7000003066 UPC 00652221047823

Details

Two-part resin and paste system Chemical cure Reliable bond strength

Typical properties

Adhesive Type: Direct Brands: Concise™

Catalog Number: 196-001 Cure Type: Chemical Cure

Factory ISO Certification: No ISO Certification information available

Product Type: Chemical Cure Adhesive Kit

Units per Inner Pack: 1 /KT

3M™ Unitek™ Miniature Twin Metal Brackets

Details

A microetched base features strength comparable to mesh

Ample area under tie wings allows easy ligation

Torque-in-base design with compound base contour

A permanent vertical scribe line facilitates more accurate placement of brackets on the long axis for bonding

A bi-directional hook eases elastic engagement and presents a smooth labial surface for increased patient comfort

Symmetrical tie wings add balance to enhance tie wing strength and improve patient aesthetics

A chamfered slot entrance promotes easier and more rapid archwire placement An integral base perimeter ridge facilitates easy flash clean up Permanent recessed I.D. markings provide ready identification

Typical properties

Adhesive Precoat: Non-coated

Angulation: +11°, +2°, +3°, +4°, +5°, +8°, +9°, 0°

Arch: Lower, Upper Base Type: Weldable Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Hook Included: Yes Hook Type: Distal Material: Stainless Steel Prescription: MBT™, Roth* Side: Left, Left/Right, Right

Slot Size (Imperial): 0.018 in, 0.022 in Slot Size (Metric): 0.457 mm, 0.559 mm

Tooth: Anterior, Bicuspid, Central Incisor, Cuspid, First Bicuspid, Lateral Incisor,

Second Bicuspid

Torque: +12°, +17°, +8°, -11°, -12°, -17°, -1°, -22°, -6°, -7°, 0°

Units per Pack: 5

3M™ Clarity™ Aligners

Details

100% doctor satisfaction with 3M[™] Clarity[™] Aligners Flex in an initial evaluation Excellent initial and sustained force persistence

Straight trimline contributes to torque control, retentive forces, esthetics and patient comfort

Nearly invisible and more stain-resistant than Invisalign® SmartTrack™ Material Thin, flexible design

Confidence backed by 75 years of experience moving teeth 100% Customer Satisfaction with 3M Clinician Support dedicated to Clarity Aligners

Available on 3M[™] Oral Care Portal, an integrated software platform that also supports 3M[™] Filtek[™] Matrix restorative treatment and 3M[™] Clarity[™] Retainers Oral Care Portal is an open system that accepts any high quality STL file

Typical properties

Brands: Clarity™

Factory ISO Certification: No ISO Certification information available

3M™ Unitek™ 804-155 Open-Close Instrument

Part Number 804-155 3M Product Number 804-155 3M ID 70000915168

Details

The 3M™ Unitek™ Debonding Instrument for the Ultra SL System is required for a proper debonding experience

3M™ Clarity™ Ultra Self-Ligating Brackets should be debonded with the archwire and ligation removed

High quality stainless steel construction

Precision-manufactured

The close end of the instrument assists with archwire insertion into the slot of the Clarity™ Ultra Self-Ligating Bracket

The use of the Unitek Open-Close Instrument is highly recommended

High quality stainless steel construction with a ceramic tip for opening the bracket door

Precision-manufactured

Typical properties

Brands: SmartClip™

Factory ISO Certification: No ISO Certification information available

3M™ Filtek™ Universal Restorative, 6500 Series

Details

Simplified shade selection with 1 universal opacity and 8 designer shades to cover most restorations

Pink Opaquer easily masks dark spots caused by metal or stained dentition Made with 3M's TRUE nanotechnology for excellent polish retention and wear resistance

Innovative monomers protect against shrinkage and shrinkage stress for lasting fracture resistance

Fluorescent pigments give restorations a natural-looking appearance 3M filler technology increases radiopacity for easier identification at recall Capsules can be safely warmed up to 70°C/158°F for up to 1 hour Complete your premium restoration with the gold standard of adhesion, 3M[™] Scotchbond[™] Universal Plus Adhesive

3M™ PerioMed™ 0.63% Stannous Fluoride Oral Rinse

Details

Stannous fluoride provides antimicrobial activity for up to 8 hours

Non-alcohol formula

Costs about the same as the leading over-the-counter mouthwashes

Helps control gum disease, tooth decay, tooth sensitivity and plaque bacteria for a healthier smile

Second-generation anti-microbial

Non-alcohol formula

Costs about the same as the leading OTC mouthwashes

Typical properties

Active Ingredient: Stannous Fluoride

Brands: PerioMed™

Condition: Caries Prevention, Gingivitis, Hypersensitivity

Delivery System: Bottle

Factory ISO Certification: No ISO Certification information available

Flavor: Mint, Tropical Fruit Material: Stannous Fluoride

3M™ RelyX™ Unicem 2 Automix Self-Adhesive Resin Cement Value Pack

Details

Cement delivery system provides easy handling and flexible dose dispensing Saves time by eliminating etching*, priming and bonding steps *Except for 2 or 3 unit Maryland bridges and 3-unit inlay/onlay bridges. Selective etching is required to increase surface area for bonding

Virtually no post-operative sensitivity

Easy removal of excess cement

High intrinsic color stability for natural-looking esthetics

Cement resists discoloration

Interested in further simplifying your direct and indirect restorative workflows with a fully aligned system? Explore 3M™ RelyX™ Universal Resin Cement

Typical properties

Brands: RelyX™

Catalog Number: 56874,56875,56876,56849,56857,56846,56847,56858,56848

Cure Type: Dual Cure

Delivery System: Automix Syringe

Factory ISO Certification: No ISO Certification information available

Material: Self-adhesive Resin

Package Contents: 3 – 8.5g Syringes; 30 Mixing Tips (regular); 20 Mixing Tips (wide); 15 Endo Tips; 5 Intra-oral Tips,Refill includes one 8.5g syringe of shade A2 Universal; 10 mixing tips (regular); 5 mixing tips (wide); 5 endo tips.,1 – 11.0g Clicker™ Dispenser; 1 Mixing Pad (70144),1 – 8.5g Syringe; 15 Mixing Tips

(regular); 5 Mixing Tips (wide); 5 Endo Tips

Shade: A2, A3, TR

3M™ Scotchbond™ Universal Adhesive Unit Dose Intro Pack

Details

Bonds dental restorations with uncompromising strength with a one-step, one-coat application

Compatible with all etching techniques

Compatible with all resin cements

Self- and dual-cure compatible with 3M™ Scotchbond™ Universal DCA Dual-

Cure Activator

Confidently bonds to both moist and dry dentin

Offers a strong bond to all indirect surfaces without an additional primer for more streamlined workflows

Results in virtually no post-operative sensitivity,* creating a better patient experience. *3M Internal data

Available in a flip-top vial for one-handed operation or unit dose delivery

Typical properties

Adhesive Type: Selective Etch, Self-etch, Total Etch

Brands: Scotchbond™

Catalog Number: 41257,41256,41255,41254,41259,41258

Cure Type: Light Cure

Delivery System: Unit Dose, Vial

Factory ISO Certification: No ISO Certification information available

Package Contents: Bulk Pack includes: 200 - 3M[™] Scotchbond[™] Universal Adhesive - Unit Dose (011ml liquid each), 1 - 3ml etchant syringe, 25 dispensing tips, 7 technique cards, 2 IFUs, Giant Pack includes: 100 - 3M[™] Scotchbond[™] Universal Adhesive - Unit Dose (011ml liquid each), 1 - 3ml etchant syringe, 25 dispensing tips, 7 technique cards, 2 IFUs, Intro kit includes: 1 - 5ml 3M[™] Scotchbond[™] Universal Adhesive - Vial, 50 mix wells, 50 microbrushes, 1 - 3ml etchant syringe, 25 dispensing tips, 7 technique cards, 2 IFUs, Intro kit includes: 40 - 3M[™] Scotchbond[™] Universal Adhesive - Unit Dose (011ml liquid each), 1 - 3ml etchant syringe, 25 dispensing tips, 7 technique cards, 2 IFUs, Refill includes 1 - 5ml vial of adhesive

Product Type: Dental Adhesive

Units per Inner Pack: 40 per Pack, 100 per Pack

3M™ Impregum™ Soft Polyether Impression Material

Details

Unique moisture tolerance results in excellent detail reproduction and void-free impressions, even in wet environments

Constant flowability throughout the whole working time makes it suitable for large and complex cases

Snap-set behavior minimizes risk of distortion for precise-fitting restorations Post-set material rigidity helps ensure accurate transfer of implant copings Available in a range of viscosities, set times and delivery systems for monophase and 1-step heavy/light impression techniques

Control moisture at every step and gently retract gingival tissue for more accurate impressions with 3M™ Astringent Retraction Paste

Typical properties

Brands: Impregum™ Catalog Number:

P3773,31740,31731,31752,31741,31785,31642,31792,31769,31766,31788,3176

5,31746,31745,31767

Delivery System: Garant™ Dispenser, Pentamix™ Mixing Unit Factory ISO Certification: No ISO Certification information available

Impression Type: Final Material: Polyether

Package Contents: 4 - 50ml cartridges; 10 - 3M[™] Garant[™] Mixing Tips - Purple,4 - 50ml cartridges; 10 - 3M[™] Garant[™] Mixing Tips - White; 5 - 3M[™] Intra-oral Syringes - Purple,1 - base paste, 300ml; 1 catalyst, 60ml,4 - 50ml cartridges; 10 - 3M[™] Garant[™] Mixing Tips -White; 5 - 3M[™] Intra-oral Syringes -

Purple, 2 - base paste, 300ml each; 2 - catalyst, 60ml each Product Color: Dark Purple, Medium Purple, Pink, Purple

Set Time: 4 min, 6 Minutes

Set Type: Quick Set, Regular Set Setting Time: 4 Minute, 6 Minute Viscosity: Heavy, Light, Medium

Work Time: 1 Minute 45 Seconds, 1 min, 2 min, 2 min, 30 sec, 2 min 45 sec, 3

min 15 sec

3M™ Unitek™ Prestige Universal Cutter

Details

Holds loose distal ends Shear-cuts hard wire 3M™ Unitek™ Prestige Series

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

3M™ Iso-Form™ Temporary Crowns

Details

Pre-finished, contoured and belled to save placement time

Available for adult posterior applications

Crown set is packaged in a unique, stackable tray that provides convenient, easy storage of prefabricated crowns

Faster placement than acrylics and chemical resin for single units for a quicker procedure

Typical properties

Arch: Lower, Upper Brands: Iso-Form™ Catalog Number:

L61,L60,L63,L62,L65,L64,L67,L66,L69,L68,U40,U42,U41,U44,U43,U46,U45,U48,U47,U49,L70,L72,L71,L74,L73,L76,L75,L78,L77,L79,U50,U52,U54,U56,U58,L41,L40,L43,L42,L45,L44,L47,U60,L46,L49,U62,L48,U61,U64,U63,U66,U65,U68,U67,U69,L50,L52,L51,L54,L53,L56,L55,L58,U71,L57,U70,U73,L59,U72,U75,U74,U77,U76,U79,U78

Dentition: Permanent

Factory ISO Certification: No ISO Certification information available Mesial/Distal Width (Imperial): 0.236 in, 0.252 in, 0.26 in, 0.264 in, 0.268 in, 0.272 in, 0.276 in, 0.28 in, 0.283 in, 0.287 in, 0.291 in, 0.295 in, 0.299 in, 0.303 in, 0.307 in, 0.315 in, 0.323 in, 0.335 in, 0.354 in, 0.374 in, 0.386 in, 0.39 in, 0.402 in, 0.406 in, 0.409 in, 0.413 in, 0.421 in, 0.425 in, 0.437 in, 0.441 in, 0.449 in, 0.453 in, 0.457 in, 0.461 in, 0.472 in, 0.488 in

Mesial/Distal Width (Metric): 6 mm, 6.4 mm, 6.6 mm, 6.7 mm, 6.8 mm, 6.9 mm, 7 mm, 7.1 mm, 7.2 mm, 7.3 mm, 7.4 mm, 7.5 mm, 7.6 mm, 7.7 mm, 7.8 mm, 8 mm, 8.2 mm, 8.5 mm, 9 mm, 9.5 mm, 9.8 mm, 9.9 mm, 10.2 mm, 10.3 mm, 10.4 mm, 10.5 mm, 10.7 mm, 10.8 mm, 11.1 mm, 11.2 mm, 11.4 mm, 11.5 mm, 11.6 mm, 11.7 mm, 12 mm, 12.4 mm

Number of Crowns: 5 Side: Left, Right

Tooth: First Bicuspid, First Molar, Second Bicuspid, Second Molar

3M™ Unitek™ Lingual Curved Mathieu Pliers, 803-571, 1/Ea

Part Number 804-152 3M Product Number 804-152 3M ID 7000001407

Details

1/each

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Slot Size (Imperial): 0.022 in Slot Size (Metric): 0.559 mm

3M™ Unitek™ Stainless Steel Primary Molar Crowns, 902150, Kit

3M ID 7000053485 UPC 30605861011212

Details

Primary molars

Typical properties

Arch: Upper/Lower Brands: Unitek™

Catalog Number: 902150

Dentition: Primary

Factory ISO Certification: No ISO Certification information available

Number of Crowns: 112

Patient: Pediatric Size: Assorted

Tooth: Assorted, First Molar

3M™ Unitek™ Bicuspid Bands

Details

Universal anatomy fits rights and lefts Regular temper Permanent laser ID

Typical properties

Arch: Lower, Upper, Upper/Lower

Band Size: 0, 0+, 00+, 1, 1+, 10, 10+, 11, 11+, 12, 12+, 13, 13+, 14, 14+, 15, 15+, 16, 16+, 17, 17+, 18, 18+, 19, 19+, 2, 2+, 20, 21, 22, 23, 24, 25, 26, 27, 28, 29, 3, 3+, 30, 31, 32, 33, 34, 35, 36, 37, 38, 4, 4+, 5, 5+, 6, 6+, 7, 7+, 8, 8+, 9, 9+

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Material: Stainless Steel Product Type: Band Tooth: Bicuspid

3M™ Unitek™ MIA Quad Helix

Details

Palatal and lingual appliance
Easy removal and reinsertion
Allows greater flexibility and control during treatment
Available in various preformed sizes with straight or curved retention sheaths

Typical properties

Arch: Lower, Upper Arch Size: 1, 2, 3, 4 Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Retention: Curved

Units per Inner Pack: 5 per Pack

3M™ Transbond™ IDB Pre-Mix Chemical Cure Adhesive Part B, 712-122, 10 g, 1/Btl

3M ID 7100008361

Details

Good bond strength for reliable indirect bonding
Permits wire insertion upon tray removal and clean-up (within three minutes)
Can be used with custombase adhesives such as 3M™ Transbond™ XT Light
Cure Adhesive
Ample work time and quick set time
Mix Resin A & B before application

Typical properties

Brands: Transbond™
Catalog Number: 712-122
Cure Type: Chemical Cure
Delivery System: Bottle

Factory ISO Certification: No ISO Certification information available

Product Type: Chemical Cure Adhesive Kit

Units per Inner Pack: 1/BT

3M™ Unitek™ Modified Arrowhead Clasp

Details

Permachrome resilient Grips single teeth independently

Typical properties

Archwire Size (Imperial): 0.024 in, 0.028 in Archwire Size (Metric): 0.61 mm, 0.711 mm

Brands: Unitek ™

Factory ISO Certification: No ISO Certification information available

Units per Inner Pack: 10 per Pack

3M™ Lava™ Ultimate CAD/CAM Restorative for Lava™ Mills

3M ID B00008186

Details

Fast, easy to mill and a lasting polish, without firing Brilliant esthetics with long-lasting polish Less wear to opposing enamel than glass ceramics Resilient, not brittle

Resin nano ceramic material delivers exceptional durability, excellent polish retention

Available in eight shades and two translucencies

Typical properties

Brands: Lava™

3M™ Transbond™ XT Light Cure Adhesive Primer, 712-034, 6 mL, 1/BT

3M ID 7000004380 UPC 00652221108944

Details

Orthodontic primer for use in dry environments Comes in a bottle

Typical properties

Adhesive Type: Direct Brands: Transbond™ Catalog Number: 712-034

Container Size (Imperial): 0.203 fl oz

Container Size (Metric): 6 mL

Cure Type: Light Cure Delivery System: Bottle

Factory ISO Certification: No ISO Certification information available

Product Type: Orthodontic Primer

Units per Inner Pack: 1/BT

3M™ Mini Uni-Twin™ Brackets Kit

Typical properties

Adhesive Precoat: Non-coated

Bracket Style: Twin Brands: Mini Uni-Twin™

Factory ISO Certification: No ISO Certification information available

Material: Stainless Steel

Slot Size (Imperial): 0.018 in, 0.022 in Slot Size (Metric): 0.457 mm, 0.559 mm

Units per Pack: 20

3M™ AlastiK™ KX Modules

Details

Not made with natural rubber latex Greater force delivery than K modules Designed for intermaxillary application

Typical properties

Brands: AlastiK™

Factory ISO Certification: No ISO Certification information available

Force Rating: Medium

Ligature Size: KX-1, KX-2, KX-3 Product Color: Gray, Transparent Units per Inner Pack: 100 per Pack

3M™ APC™ Inventory Dispensing System, Drop In Organizer Tray, 2600-007, 1/Ea

3M ID 70006629474

Details

Contains 2 holding tabs for easy insertion/removal of tray in the drawer Each channel holds four 5-packs of APC blisters

Typical properties

Brands: 3M™

Factory ISO Certification: No ISO Certification information available

3M™ Unitek™ Universal Separating Pliers, 800-734, 1/Ea

3M ID 7100069181

Details

Facilitates the making of Omega, "Bull" and closing loops

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Pliers Type: Tweed Loop Forming

3M™ Unitek™ Lingual Buttons

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available eClass 14 Classification Group: 34110603

3M™ Unitek™ Lingual Straight Mosquito Forceps, 803-573, 1/Ea

3M ID 7100154553

Details

Double-ended instrument Doublesided head is designed to work interproximally, even where space is limited

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

3M™ Unitek™ Gurin Finger Spring, For Archwire .018/.022, 565-200, 1/Pk

3M ID 70202083617

Details

Anchor without welding or soldering Adjustable

Typical properties

Brands: Unitek™

Compatible Archwire Size (Imperial): 0.018 - 0.022 in

Factory ISO Certification: No ISO Certification information available

Units per Inner Pack: 1 per Pack

3M™ Unitek™ Gurin Nuts

Details

Anchor without welding or soldering Adjustable

Typical properties

Brands: Unitek™

Compatible Archwire Size (Imperial): 0.014 - 0.016 in, 0.018 - 0.022 in, 0.038 -

0.045 in

Factory ISO Certification: No ISO Certification information available

3M[™] Unitek[™] Prestige Double-Ended Dual Purpose Instrument, 900-815, 1/Ea

3M ID 7000025687

Details

Sterilizable

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

3M™ Filtek™ Z250 XT Nano Hybrid Universal Restorative

Details

Easy to use for all classes of direct restorations

Easy to polish for predictable esthetics

Great handling without slumping or sticking to instruments

High strength, with long-lasting wear resistance Available in 12 radiopaque, fluorescent shades

Complete your premium restoration with the gold standard of adhesion, 3M™

Scotchbond™ Universal Plus Adhesive

Typical properties

Brands: Filtek™ Catalog Number:

1470OA3,1470A3.5,1470C2,1470B1,1470A4,1471A3,1470A3,1471A2,1470B3,1

471A3.5,1471B2,1470A2,1471A1,1470A1,1471B1,1470B2

Composite Type: Universal Cure Type: Light Cure

Delivery System: Capsule, Syringe

Factory ISO Certification: No ISO Certification information available

Inner Packs per Case: 10 Packs per Case Material Structure: Nanohybrid Composite

Net Weight (Metric): 0.2 g, 3 g

Product Form: Refill

Shade: A1, A2, A3, A3.5, A4, B1, B2, B3, C2, OA3 Total Units per Case: 10 per Case, 200 per Case Units per Inner Pack: 1 per Pack, 20 per Pack

Units per Pack: 1, 20

3M™ RelyX™ Luting Cement Refill

Details

Easy clean-up
Sustained fluoride release
Virtually no post-operative sensitivity
If you like 3M™ RelyX™ Luting Cement, you'll love 3M™ RelyX™ Universal
Resin Cement

Typical properties

Brands: RelyX™

Catalog Number: 3505L,3505P

Cure Type: Self-cure Delivery System: Hand Mix

Factory ISO Certification: No ISO Certification information available

Material: Resin-modified Glass Ionomer

Package Contents: 1 – 9ml Liquid,1 – 16g Powder, Spoon

Shade: White Opaque

3M™ Scotchbond™ Universal Adhesive Vial Intro Kit

3M ID B00007530

Details

Bonds dental restorations with uncompromising strength with a one-step, one-coat application

Compatible with all etching techniques

Compatible with all resin cements

Self- and dual-cure compatible with 3M™ Scotchbond™ Universal DCA Dual-Cure Activator

Confidently bonds to both moist and dry dentin

Offers a strong bond to all indirect surfaces without an additional primer for more streamlined workflows

Results in virtually no post-operative sensitivity,* creating a better patient experience. *3M Internal data

Available in a flip-top vial for one-handed operation or unit dose delivery

Typical properties

Brands: Scotchbond™ Cure Type: Light Cure

Product Type: Dental Adhesive

3M™ Impregum™ F Polyether Impression Material

Details

Hydrophilic properties for moisture control and accuracy
Thixotropic material flows smoothly and does not slump
Ample working time and snap set behavior reduces potential for distortion
Polyether characteristics help ensure predictable, accurate impressions
For use with the monophase impression technique

Typical properties

Brands: Impregum™

Catalog Number: 31474,31504,31712

Delivery System: Hand Mix

Factory ISO Certification: No ISO Certification information available

Impression Type: Final Material: Polyether

Package Contents: 1 - tube catalyst, 15ml,1 - tube base, 120 ml,2 - tubes base,

120ml each; 2 - tube catalyst, 15ml each; 1 - mixing pad

Set Time: 6 Minutes Set Type: Regular Set Setting Time: 6 Minute Viscosity: Medium

Work Time: 2 min 45 sec

3M™ Imprint™ 4 Penta™ Putty VPS Impression Material Refill

Details

VPS impression material with innovative active self-warming to accelerate set time

Very short intra-oral setting time

Sufficient working time to reduce stress in your procedure

Hydrophilic to capture even the finest details in wet environments

Can be used for one-step or two-step techniques

Control moisture at every step and gently retract gingival tissue for more accurate impressions with 3M™ Astringent Retraction Paste

Typical properties

Brands: Imprint™ Catalog Number:

71488,71532,71489,71486,71487,71531,71484,71485,71493,71518,71492,7149

0

Delivery System: Garant™ Dispenser, Pentamix™ Mixing Unit Factory ISO Certification: No ISO Certification information available

Impression Type: Final, Preliminary

Material: VPS

Package Contents: 4 - 50ml cartridges; 5 - 3M[™] Garant[™] Mixing Tips - Yellow, 5 - 3M[™] Intra-oral Syringes - Green,4 - 50ml cartridges; 5 - 3M[™] Garant[™] Mixing Tips - Yellow; 5 - 3M[™] Intra-oral Syringes - Green,4 - 50ml cartridges; 10 - 3M[™] Garant[™] Mixing Tips - Green,10 - 50ml cartridges,2 - base paste, 300ml each; 2 - catalyst, 60ml each

Product Color: Lavender, Lime Green, Rose Pink, Royal Blue, Sky Blue, Yellow

Set Time: 1.25 min, 2 min

Set Type: Regular Set, Super Quick Set Setting Time: 1.25 Minute, 2 Minute

Units per Case: 10.0

Viscosity: Heavy, Light, Regular

Work Time: 1 Minute 15 Seconds, 1 min, 2 min, 35 Second

3M™ Paradigm™ DeepCure Eye Shield, 76986

Part Number 76986 3M ID 7100057735 UPC 04035077010252

Details

Protects patient and user from reflection and intensive scattered light

Typical properties

Brands: Paradigm™

Factory ISO Certification: No ISO Certification information available

Package Contents: 1 - Paradigm DeepCure Eye Sheild

Product Type: Eye Shield

3M[™] Vitremer[™] Core Buildup/Restorative Finishing Gloss Refill, 3303FG

3M ID 7000030407 UPC 30605861032774

Details

Light curable finishing gloss to maximize the esthtics of 3M[™] Vitremer[™] core buildup/restorative system

Unique dark-cure chemistry allows for bulk placement which eliminates the need for layering and saves time

Displays a true glass ionomer cure for long-term fluoride release and adhesion to tooth structure which provides better patient protection and increased doctor confidence

As a restorative material, it provides good handling properties and extended working time in the mouth for easier use

Available in a variety of shades to better match existing tooth structure for esthetic restorations. Pedo shade is available for pedodontic applications and Blue shade for core buildup applications

Typical properties

Brands: Vitremer™

Catalog Number: 3303FG Cure Type: Light Cure Delivery System: Bottle

Factory ISO Certification: No ISO Certification information available

Inner Packs per Case: 10 Packs per Case

Package Contents: Refill includes 65ml finishing gloss

Shade: Translucent

Total Units per Case: 10 per Case Units per Inner Pack: 1 per Pack

3M™ Composite Placement Instruments (10 of sizes 1, 2), 1948

Part Number 1948 3M ID 7000129232 UPC 30605861057319

Details

Allows for easy intraoral placement of dental composites Double-ended instrument Disposable for improved infection control

Typical properties

Brands: 3M™

Factory ISO Certification: No ISO Certification information available

3M™ Clinpro™ Sealant Mixing Well With Cover, 12621MW, 4 per box

3M ID 7000128791 UPC 30605861034938

Details

Mixing wells that work with Clinpro™ Sealant system

Typical properties

Brands: Clinpro™

Factory ISO Certification: No ISO Certification information available

Number of Wells: 3

3M[™] Penta[™] Elastomer Syringe Spare Syringe Barrels Refill, 71200, 10/Pack

3M ID 7100215422

Details

Spare barrels for 3M[™] Penta[™] Elastomer Syringe which allows for easy, intraoral delivery of Penta[™] impression materials

Typical properties

Brands: Penta™

Catalog Number: 71200

Factory ISO Certification: No ISO Certification information available

Package Contents: 10 - Syringe Barrels

Product Type: Barrel

3M™ Putty Spoons Refill, 7314

3M ID 7000053455 UPC 30605861026599

Details

For hand mixing putty materials

Typical properties

Brands: 3M™

Catalog Number: 7314

Factory ISO Certification: No ISO Certification information available

Package Contents: 2 - Putty spoons

Product Type: Putty Spoon

3M™ Directed Flow Impression Tray

Details

Revolutionary self-retentive fleece strip eliminates the need for a tray adhesive for fewer handling steps

Simplified process saves time--up to five minutes per procedure

Retentive fleece features striations that act as a barrier to avoid tooth/tray contact and occlusal flat-spots for better impression results

Innovative, directed-flow design minimizes flow defects and distal voids for more accurate impressions

Surrounding retention lip reduces material overflows

Built-in reservoir reduces gagging for better patient comfort

Single-use tray means there is no need to disinfect and clean, as with metal trays, for improved infection control

Delivers outstanding quality, accurate results and easy handling

Rigid material provides the stability that is ideal for both preliminary and final impressions

No unpleasant odor or mess from a tray adhesive

Tray can be customized, if needed

Typical properties

Arch: Lower, Upper Brands: 3M™

Catalog Number: 71620,71619,71617,71618,71615,71616

Factory ISO Certification: No ISO Certification information available

Package Contents: 10 - trays, size medium (M), lower, 10 - trays, size medium (M), upper, 10 - trays, size large (L), lower, 10 - trays, size small (S), lower, 10 -

trays, size small (S), upper,10 - trays, size large (L), upper

Tray Size: Large, Medium, Small

3M™ Unitek™ Stainless Steel Permanent Bicuspid Crowns

Details

Easy to customize with parallel walls to allow for clinician fit control Pre-trimmed to optimum length and contour to save placement time Made in the U.S.A. of high quality stainless steel Shallow occlusal anatomy means less tooth reduction Thick occlusal surface to prevent bite-through Available in molar styles to meet your adult patient needs Crown set is packaged in a unique, stackable tray that provides convenient, easy storage of 3M™ ESPE™ prefabricated crowns Wide array of shapes and sizes

Typical properties

Arch: Lower, Upper Brands: Unitek™ Catalog Number:

900536,900514,900635,900513,900535,900634,900534,900512,900633,900511,900533,900632,900532,900510,900631,900531,900630,900530,900509,900526,900646,900525,900645,900546,900524,900644,900545,900523,900643,900544,900522,900642,900543,900521,900641,900542,900520,900541,900640,900540,900519,900516,900636,900515

Dentition: Permanent

Factory ISO Certification: No ISO Certification information available

Mesial/Distal Width (Imperial): 0.22 in, 0.24 in, 0.244 in, 0.26 in, 0.28 in, 0.299 in,

0.319 in, 0.339 in, 0.358 in

Mesial/Distal Width (Metric): 5.6 mm, 6.1 mm, 6.2 mm, 6.6 mm, 7.1 mm, 7.6 mm,

8.1 mm, 8.6 mm, 9.1 mm Number of Crowns: 5 Patient: Pediatric Side: Left, Right

Tooth: First Bicuspid, Second Bicuspid eClass 14 Classification Group: 34111100

3M™ Sof-Lex™ Finishing and Polishing System Kit, 2385P

3M ID 7100066535 UPC 30605861059092

Details

Easy-to-use system designed to create high-luster restorations that offer lifelike esthetics

Discs are color coded to allow for easy identification of grit sequence

Pop-on disc technology designed to save you time

Finishing and polishing strips are designed to produce a smooth restoration to help prevent interproximal plaque buildup for superior finishing and polishing results

Typical properties

Brands: Sof-Lex™

Factory ISO Certification: No ISO Certification information available

Grades Included: Assorted

Grit Size: Assorted

Kit Type: Finishing & Polishing System Kit

Package Contents: System kit includes 240 contouring and polishing discs, 240 extra-thin contouring and polishing discs, 45 finishing and polishing strips, and 1

RA mandrel

3M™ Filtek™ Supreme Ultra Universal Restorative

Details

Combines exceptional esthetics, uncompromising fracture toughness and outstanding strength for any restoration

Excellent polish retention for lasting shine

Superior wear resistance

Versatile application for anterior and posterior

Exceptional adaptation and handling

Easy-to-use shading system color-coded by opacity

Natural-looking results with a selection of 36 radiopaque shades and 4 opacities Complete your premium restoration using the gold standard of adhesion with 3M™ Scotchbond™ Universal Plus Adhesive

Typical properties

Brands: Filtek™ Catalog Number:

6028A2E,6029B1E,6029A3D,6028A2B,6029C3B,6029B1B,6029A3E,6028C2B,6 028A2D,6028GT,6028CT,6029AT,6028A6B,6029B5B,6029A3B,6029B2E,6029D 2E,6028B1E,6029A4D,6028A3B,6029B2B,6029C4D,6028C3B,6028A3D,6029D2 B,6029A3.5B,6028B1B,6028A3E,6028WB,6028WD,6028WE,6028BT,6029A4B,6028B5B,6029A1E,6028D2B,6028B2E,6028D2E,6029B3B,6028A4D,6029A1B,6 029D3B,6029B3D,6028B2B,6029C1B,6029A1D,6028C4D,6029GT,6029XWE,60 28AT,6028A3.5B,6029CT,6029XWB,6028C1B,6028A1D,6028A1E,6029A2B,602 8B3B,6029C2B,6029A2D,6028A1B,6028D3B,6029A2E,3028A4B,6029 WB,6029WE,6029WD,6029BT,6028XWE,6028XWB,6029A6B,6029A2B-S

Composite Type: Universal

Delivery System: Capsule, Syringe Environmental claim being made: Yes Environmental claim relates to: Product Factory ISO Certification: ISO 13485 Material Structure: Nanocomposite Net Weight (Metric): 0.2 g, 4 g

Opacity: Body, Dentin, Enamel, Translucent

Package Contents: Refill includes 10 - 02g capsules in Translucent shade Gray,Refill includes 10 - 02g capsules in Dentin shade A4D,Refill includes 1 - 4g syringe in Translucent shade Gray,Refill includes 20 - 02g capsules in Body shade B3B,Refill includes 20 - 02g capsules in Body shade B3B,Refill includes 20 - 02g capsules in Body shade D3B,Refill includes 1 - 4g syringe in Enamel shade White,Refill includes 1 - 4g syringe in Dentin shade B3D,Refill includes 10 - 02g capsules in Enamel shade A2E,Refill includes 1 - 4g syringe in Dentin shade A1D,Refill includes 1 - 4g syringe in Body shade C3B,Refill includes 1 - 4g syringe in Body shade White,Refill includes 1 - 4g syringe in Body shade A3B,Refill includes 1 - 4g syringe in Body shade B5B,Refill includes 1 - 4g syringe in Body shade B5B,Refill includes 1 - 4g syringe in Body shade B5B,Refill includes 1 - 02g capsules in Dentin shade C4D,Refill includes 1 - 4g syringe in Body shade B1B,Refill includes 10 - 02g capsules in Translucent shade Amber,Refill includes 20 - 02g capsules in Body shade A6B,Refill includes 1 - 4g syringe in Translucent

shade Clear, Refill includes 20 - 02g capsules in Body shade A2B, Refill includes 20 - 02g capsules in Body shade White, Refill includes 10 - 02g capsules in Enamel shade B1E, Refill includes 20 - 02g capsules in Body shade C2B, Refill includes 1 - 4g syringe in Dentin shade A2D, Refill includes 10 - 02g capsules in Translucent shade Clear, Refill includes 10 - 02g capsules in Enamel shade A3E, Refill includes 1 - 4g syringe in Body shade B2B, Refill includes 10 - 02g capsules in Enamel shade White, Refill includes 1 - 4g syringe in Body shade A4B, Refill includes 10 - 02g capsules in Enamel shade Extra White, Refill includes 1 - 4g syringe in Enamel shade A2E, Refill includes 10 - 02g capsules in Dentin shade A3D, Refill includes 1 - 4g syringe in Body shade D2B, Refill includes 1 - 4g syringe in Dentin shade White, Refill includes 20 - 02g capsules in Body shade C3B.Refill includes 20 - 02g capsules in Body shade A3B.Refill includes 20 - 02g capsules in Body shade B5B, Refill includes 10 - 02g capsules in Enamel shade D2E, Refill includes 20 - 02g capsules in Body shade B1B, Refill includes 10 - 02g capsules in Enamel shade B2E, Refill includes 10 - 02g capsules in Dentin shade White, Refill includes 1 - 4g syringe in Dentin shade A3D, Refill includes 1 - 4g syringe in Enamel shade A1E, Refill includes 1 - 4g syringe in Translucent shade Amber, Refill includes 1 - 4g syringe in Body shade A1B, Refill includes 1 - 4g syringe in Body shade B3B, Refill includes 1 - 4g syringe in Body shade C1B, Refill includes 1 - 4g syringe in Body shade D3B, Refill includes 10 - 02g capsules in Dentin shade A2D, Refill includes 20 -02g capsules in Body shade B2B, Refill includes 20 - 02g capsules in Body shade Extra White.Refill includes 20 - 02g capsules in Body shade A4B,Refill includes 20 - 02g capsules in Body shade D2B, Refill includes 1 - 4g syringe in Dentin shade C4D, Refill includes 10 - 02g capsules in Enamel shade A1E, Refill includes 1 - 4g syringe in Dentin shade A4D, Refill includes 1 - 4g syringe in Enamel shade D2E, Refill includes 1 - 4g syringe in Enamel shade B2E, Refill includes 1 - 4g syringe in Body shade A6B, Refill includes 1 - 4g syringe in Body shade A2B, Refill includes 10 - 02g capsules in Translucent shade Blue, Refill includes 10 - 02g capsules in Dentin shade A1D, Refill includes 1 - 4g syringe in Translucent shade Blue, Refill includes 20 - 02g capsules in Body shade A35B, Refill includes 1 - 4q syringe in Enamel shade Extra White, Refill includes 1 - 4g syringe in Body shade C2B, Refill includes 10 - 02g capsules in Dentin shade B3D, Refill includes 1 - 4g syringe in Body shade Extra White Product Form: Sample

Shade: A1B, A1D, A1E, A2B, A2D, A2E, A3.5B, A3B, A3D, A3E, A4B, A4D, A6B, AT, B1B, B1E, B2B, B2E, B3B, B3D, B5B, BT, C1B, C2B, C3B, C4D, CT, D2B, D2E, D3B, GT, WB, WD, WE, XWB, XWE

Units per Inner Pack: 1 per Pack, 10 per Pack, 20 per Pack

Units per Pack: 1, 10, 20

3M™ AlastiK™ Continuous Chains Bobbin

Details

Not made with natural rubber latex 15 ft. (457 cm) per bobbin Die-cut chain Color coordinated bobbin for easy identification Medium force

Typical properties

Brands: AlastiK™

Center to Center Distance (Imperial): 0.11 in, 0.14 in

Dispenser Type: Spool

Factory ISO Certification: No ISO Certification information available

Ligature Size: C-1, CK

Product Color: Assorted, Blue Metallic, Deep Black, Dusky Blue, Extreme Green,

Gold, Gray, Green Metallic, Hot Red, Kelly Green, Lavender, Light Blue,

Mediterranean Blue, Midnight Glow, Navy Blue, Obscure, Orange Blast, Plum Crazy, Purple Metallic, Radical Red, Screaming Pink, Silver, Smoke, Teal,

Transparent, White, Yellow

Units per Inner Pack: 1 per Pack

3M™ Unitek™ Prestige Band Seating Pliers

3M ID B00045652

Details

Serrated tips grip prewelded bracket wings securely while fitting bands $3M^{\text{TM}}$ Unitek $^{\text{TM}}$ Prestige Series

Typical properties

Brands: Unitek™

3M™ AlastiK™ Easy-To-Tie Ligatures

Details

Not made with natural rubber latex
Choose from 26 exciting colors
All colors are fade resistant for a minimum of one month in artificial saliva, 8% ethanol in water, cottonseed oil and carbonated citrated sodas
One size fits most twin brackets
Available as ligatures or chains

Typical properties

Brands: AlastiK™ Dispenser Type: Stick

Factory ISO Certification: No ISO Certification information available

Product Color: Assorted, Blue Metallic, Deep Black, Dusky Blue, Extreme Green,

Gold, Gray, Green Metallic, Hot Red, Kelly Green, Lavender, Light Blue, Mediterranean Blue, Midnight Glow, Navy Blue, Obscure, Orange Blast, Plum Crazy, Purple Metallic, Radical Red, Screaming Pink, Silver, Smoke, Teal,

Transparent, White, Yellow

Units per Inner Pack: 100 per Pack

3M[™] Unitek[™] Miniature Twin V-Slot Metal Brackets Kit, Roth*, .022, U/L, 5x5, Cuspid Hk, 017-134, 20/Kit

3M ID 70202011980

Details

Permanent raised I.D. markings provide ready identification

A permanent vertical scribe line facilitates more accurate bracket placement on the long axis for bonding

Symmetrical tie-wings add balance to enhance tie-wing strength and improve patient esthetics

An integral base perimeter ridge facilitates easy flash cleanup

Torque-in-base design with compound base contour

Bracket design provides additional reference points to ensure proper bracket placement regardless of tooth position

A chamfered slot entrance promotes easier and more rapid archwire placement

Typical properties

Arch: Upper/Lower Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Hook Included: Yes Hook Type: Distal Prescription: Roth*

Slot Size (Imperial): 0.022 in Slot Size (Metric): 0.559 mm

Units per Pack: 20

3M™ Unitek™ Debonding Instrument

Details

3M[™] Clarity[™] Ultra Self-Ligating Brackets have one instrument used for both opening and closing the doors

3M™ Unitek™ Open-Close Instrument is designed with a ceramic tip for opening the bracket door

The open end of the instrument allows you to open the Clarity™ Ultra Self-Ligating Bracket without overstressing the door

Closing the Clarity™ Ultra Self-Ligating Bracket doors with the 3M™ Unitek™ Open-Close Instrument is easy

The close end of the instrument assists with archwire insertion into the slot of the Clarity™ Ultra Self-Ligating Bracket

The use of the Unitek Open-Close Instrument is highly recommended

High quality stainless steel construction with a ceramic tip for opening the bracket door

Precision-manufactured

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

3M™ Unitek™ Prestige Pin and Ligature Cutter, 900-753, 1/Ea

3M ID 7000025681

Details

Tips for finite cutting of ligature wires 3M™ Unitek™ Prestige Series

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

3M™ Unitek™ Tweed Loop Forming Pliers, 801-720, 1/Ea

3M ID 7000024327

Details

5/pack

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

3M[™] Transbond[™] Plus Self Etching Primer, 712-092, Applicator 100/Pk

3M ID 7100154493

Details

100/pack

Typical properties

Brands: Transbond™

Factory ISO Certification: No ISO Certification information available

Product Type: Sponge Applicator

3M™ Adper™ Single Bond 2 Vial Refill

Details

Offers good bond strength and performance, providing your patients with lasting dental restorations

Convenient bottle design offers easier opening and dispensing

No shaking required thanks to stable nano-filler designed not to settle within the bottle

Reduces the risk of post-operative sensitivity, creating a better experience for your patients

Compatible with light-cured composite restoratives

Typical properties

Adhesive Type: Total Etch

Brands: Adper™

Catalog Number: 51202 Cure Type: Light Cure Delivery System: Bottle, Vial

Delivery System. Bottle, viai

Factory ISO Certification: No ISO Certification information available

Inner Packs per Case: 10 Packs per Case

Product Type: Dental Adhesive Total Units per Case: 10 per Case Units per Inner Pack: 1 per Pack, 1 /BT eClass 14 Classification Group: 34110205

3M™ Transbond™ IDB Pre-Mix Chemical Cure Adhesive Part A, 712-121, 10 g, 1/Btl

3M ID 7100008360

Details

Good bond strength for reliable indirect bonding
Permits wire insertion upon tray removal and clean-up (within three minutes)
Can be used with custombase adhesives such as 3M™ Transbond™ XT Light
Cure Adhesive
Ample work time and quick set time
Mix Resin A & B before application

Typical properties

Brands: Transbond™
Catalog Number: 712-121
Cure Type: Chemical Cure
Delivery System: Bottle

Factory ISO Certification: No ISO Certification information available

Product Type: Chemical Cure Adhesive Kit

Units per Inner Pack: 1/BT

3M™ Clinpro™ Glycine Prophy Powder, 67008, 6.5oz, 2 bottles

3M ID 7100036776 UPC 04035077009386

Details

For supra- and subgingival plaque biofilm removal

Proven to remove biofilm in periodontal pockets more effectively and gentler on dental tissues than conventional procedures using currettes and scalers Safe removal of plaque biofilm from soft root surfaces, dental restorations, implants and around orthodontic brackets

Can be used in commercial air-polishing equipment

Pleasant tasting

Effective for areas difficult to reach with paste

Typical properties

Brands: Clinpro™ Cup Flexibility: Firm

Factory ISO Certification: No ISO Certification information available

Latex-free: Yes

3M™ Adper™ Scotchbond™ Multi-Purpose Adhesive Refill, 7543, 8ml

3M ID B00007557 UPC 30605861035638

Details

Improve patient outcomes with a total-etch system that offers excellent bond strength to a variety of dental surfaces

Creates a complete, versatile dental bonding system when used with 3M[™] Adper[™] Scotchbond[™] Multi-Purpose Primer

Obtain strong and reliable bonds both to wet and dry surfaces, even when complete isolation from moisture is difficult

Reduces microleakage that can contribute to post-operative sensitivity, helping to ensure a positive patient experience

Typical properties

Adhesive Type: Total Etch

Brands: Adper™

Catalog Number: 7543

Container Size (Imperial): 0.271 fl oz

Container Size (Metric): 8 mL

Cure Type: Light Cure Delivery System: Vial

Factory ISO Certification: No ISO Certification information available

Package Contents: Refill includes 8 ml bottle

Product Type: Dental Adhesive

3M™ Imprint™ 4 VPS Impression Materials

Details

VPS impression material with innovative active self-warming to accelerate set time

Very short intra-oral setting time

Sufficient working time to reduce stress in your procedure

Hydrophilic to capture even the finest details in wet environments

Can be used for one-step or two-step techniques

Control moisture at every step and gently retract gingival tissue for more accurate impressions with 3M™ Astringent Retraction Paste

Typical properties

Brands: Imprint™ Catalog Number:

71536,71534,71535,71499,71554,71497,71498,71553,71495,71496,71494 Delivery System: Garant™ Dispenser, Multiple Deliveries, Pentamix™ Mixing

Unit

Factory ISO Certification: No ISO Certification information available

Impression Type: Final

Material: VPS

Package Contents: 1 - base paste, 300ml each; 1 - catalyst, 60ml each; 1 - 50ml cartridge Super Quick Light; 1 - 50ml cartridge Super Quick Regular; 1 - 3M™ Penta[™] cartridge; 10 -3M[™] Penta Mixing Tips, red; 5 - 3M[™]Intra-oral Syringes -Green; 5 - 3M™ Garant™ Mixing Tips - Yellow; 5 - 3M™ Garant™ intraoral Tips - Yellow; 1 - 3M[™] VPS tray adhesive, 17ml,2 - 50ml cartridges Heavy; 1 - 50ml cartridge Light; 1 - 50ml cartridge Regular; 10 - 3M™ Intra-oral Syringes -Green; 10 - 3M™ Garant™ Mixing Tips - Green; 1 - 3M™ VPS tray adhesive, 17ml,1 -Imprint™ 4 Penta Heavy Refill; 1 - Light cartridge, 50ml; 1 - Regular cartidge, 50ml; 1 - 3M[™] Penta[™] Cartridge; 10 - 3M[™] Penta[™] Mixing Tips, red; 5 - 3M[™] Intra-oral Syringes - Green; 5 - 3M™ Garant™ Mixing Tips - Yellow; 5 - 3M™ Garant™Intra-oral Tips - Yellow; 1 - 3M™ VPS Tray Adhesive, 17ml.,1 -Imprint[™] 4 Penta[™] Super Quick Heavy base paste, 300ml; 1 - catalyst, 60ml; 1 cartridge Super Quick Light, 50ml; 1 - cartridge Super Quick Regular, 50 ml; 1-3M™ Penta™ Cartridge; 10 3M™ Penta™ Mixing Tips, red; 5 - 3M™ Intra-oral Syringes - Green; 5 - 3M™ Garant™ Mixing Tips - Yellow; 5 - 3M™ Garant™ Intraoral Tips - Yellow; 1 - 3M[™] VPS Tray Adhesive, 17ml,2 - 50ml cartridges Super Quick Heavy; 1 - 50 ml cartridge Super Quick Light; 1 - 50 ml cartridge Super Quick Regular; 10 - 3M™ Intra-oral Syringes - Green; 10 - 3M™ Garant™ Mixing Tips - Green; 1 - 3M™ VPS Tray Adhesive, 17ml

Set Time: 1.25 min, 2 min, 2.5 min Set Type: Regular Set, Super Quick Set

Setting Time: 2 Minute

Viscosity: Heavy, Multiple Viscosities

Work Time: 1 Minute 15 Seconds, 1 Minute 30 Seconds, 2 min

3M[™] Express[™] 2 Penta[™] H VPS Impression Material Cartridge for Pentamix[™] 3, P3810 Series

3M ID B00007763

Details

For consistent, homogeneous and void-free mixing at the touch of a button Compatible with 3M™ Pentamix™ 3 Mixing Unit Reusable, metal-reinforced cartridge Helps get an accurate impression on the first take, reducing the need for costly adjustments, remakes and retakes

Typical properties

Brands: Express™

3M™ AlastiK™ K Modules

Details

Not made with natural rubber latex For intermaxillary Class II and intramaxillary Class III applications

Typical properties

Brands: AlastiK™

Factory ISO Certification: No ISO Certification information available

Force Rating: Heavy, Medium Ligature Size: K-0, K-1, K-2, K-3 Product Color: Gray, Transparent Units per Inner Pack: 100 per Pack

3M™ Unitek™ Prestige Universal Cutter, Long Handle, 900-757, 1/Ea

3M ID 7000025684

Details

Facilitates 3M™ AlastiK module placement in hard to reach places

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Units per Inner Pack: 1 per Pack

3M™ APC™ Inventory Dispensing System, 2600-023, Label Set, 1/Ea

3M ID 7100286982

Details

Helps identify the content of the drawers and back-up storage unit Drawer labels identify bracket type and maxillary or mandibular brackets

Typical properties

Brands: 3M™

Factory ISO Certification: No ISO Certification information available

3M™ Unitek™ Prestige Jaraback Pliers, 900-704, 1/Ea

3M ID 7000025659

Details

Consistent bends are assured with precision-aligned tips 3M™ Unitek™ Prestige Series

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Pliers Type: Three-jaw Wire Bending

3M™ Unitek™ Combination Lingual Hook/Seating Lugs

Details

Hold elastics securely from mesial, distal, or occlusal Straight occlusal edge for band seatin Smooth edges to protect against tissue irritation Strip scored for easy breakaway

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

3M™ MBT™ Appliance System Bracket Positioning Guide, 701-725, Plastic, 1/Ea

3M ID 7100240269

Details

Provide a range of heights in mm to guide in placement for all teeth 4/pack

Typical properties

Brands: MBT™

Factory ISO Certification: No ISO Certification information available

3M™ Unitek™ Action Line Tubing

Details

Not made with natural rubber latex Natural color elastomeric tubing Self-sealing bobbin

Typical properties

Arbor Diameter (Imperial): 0.02 in, 0.025 in, 0.03 in

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Units per Inner Pack: 1 Each

eClass 14 Classification Group: 34110690

Dimensions and Classifications

Overall Length (Imperial): 25 ft

3M™ Unitek™ Gurin Lock

Details

Anchor without welding or soldering Adjustable

Typical properties

Brands: Unitek™

Compatible Archwire Size (Imperial): 0.014 - 0.016 in, 0.018 - 0.022 in, 0.028 -

0.036 in, 0.038 - 0.045 in

Factory ISO Certification: No ISO Certification information available

3M™ Sof-Lex™ Square Eyelet Discs, 1958

Details

Easy-to-use system allows you to create high-luster restorations that offer lifelike esthetics

Discs are color coded to allow for easy identification of grit sequence Pop-on disc technology makes changing discs easier to save you time Pop-on discs conform to the contours of the teeth for easier finishing and polishing

Typical properties

Brands: Sof-Lex™

Factory ISO Certification: No ISO Certification information available

Grade: Coarse, Fine, Medium, Super Fine

Model: Standard

Package Contents: Refill includes 100 discs of 3M[™] Sof-Lex[™] Square Eyelet Discs, 5/8" Diameter, Medium grit,Refill includes 100 discs of 3M[™] Sof-Lex[™] Square Eyelet Discs, 5/8" Diameter, Superfine grit,Refill includes 100 discs of 3M[™] Sof-Lex[™] Square Eyelet Discs, 5/8" Diameter, Fine grit,Refill includes 100

discs of 3M™ Sof-Lex™ Square Eyelet Discs, 5/8" Diameter, Coarse grit

Product Color: Black, Dark Blue, Light Blue, Medium Blue

Product Type: Contouring & Polishing Disc eClass 14 Classification Group: 34111407

Dimensions and Classifications

Overall Diameter (Imperial): 0.625 in Overall Diameter (Metric): 15.8 mm

3M™ Sof-Lex™ Diamond Polishing System, 5092

Details

Imparts paste-like gloss in the convenience of a rubberized system Unique, flexible shape adapts to all tooth surfaces Multi-use, can be sterilized and reused High, long-lasting gloss when used with Filtek™ Supreme Ultra Universal Restorative

Typical properties

Brands: Sof-Lex™

Catalog Number: 5092-IM

Factory ISO Certification: No ISO Certification information available

Kit Type: Diamond Polishing System Kit

Package Contents: Refil includes 5 pre-polishing spirals beige, 5 diamond

polishing spirals pink, 3 RA mandrels eClass 14 Classification Group: 34111407

3M™ Filtek™ Supreme Ultra Universal Restorative Single Shade Intro Kit

3M ID B00007966

Details

Exceptional esthetics
Excellent adaptation
Excellent polish retention
Exceptional handling
Outstanding strength for anterior and posterior
Excellent wear resistance
Natural-looking fluorescence
True nanotechnology
36 shades, 4 opacities

Typical properties

Brands: Filtek™

Composite Type: Universal Product Type: Introduction Kit

3M™ Transbond™ Supreme LV Adhesive

Details

Flowable orthodontic adhesive for indirect bonding Excellent strength, flow and wear properties Stays in place without slumping or running prior to light curing Metal syringe tips increase accuracy, reducing waste Compatible with all 3M primers

Typical properties

Brands: Transbond™
Catalog Number: 712-046
Cure Type: Light Cure
Delivery System: Syringe

Factory ISO Certification: No ISO Certification information available

Product Type: Light Cure Band Adhesive Kit

Units per Inner Pack: 1 /KT

3M™ Vitremer™ Core Buildup/Restorative Primer, 3303P, 6.5 mL

3M ID 7000054177 UPC 30605861032736

Details

Used to prime enamel and dentin surfaces prior to restorations and core-buildups

Typical properties

Brands: Vitremer™
Catalog Number: 3303P
Cure Type: Light Cure
Delivery System: Bottle

Factory ISO Certification: No ISO Certification information available

Inner Packs per Case: 10 Packs per Case
Package Contents: Refill includes 65ml primer

Total Units per Case: 10 per Case Units per Inner Pack: 1 per Pack

3M™ Disposable Mini-Sponge Applicators

Details

Designed for adhesive application

Typical properties

Brands: 3M™

Factory ISO Certification: No ISO Certification information available

Package Contents: Refill includes approximately 200

Product Type: Sponge Applicator eClass 14 Classification Group: 34110000

3M™ Elastomer Syringe and Filling Device WE/EE, 71209

3M ID 7100287344 UPC 04064035048396

Details

The included filling device makes for easy filling of the syringe

Typical properties

Brands: 3M™

Catalog Number: 71209

Factory ISO Certification: No ISO Certification information available

Package Contents: 1 - Elastomer Syringe; 1 - Filling Device; 4 - Intraoral Tips

Product Type: Filling Device

3M™ Brush Holders

Details

Designed to hold brush tips 4/pack

Typical properties

Brands: 3M™

Factory ISO Certification: No ISO Certification information available Product Type: Brush Holder eClass 14 Classification Group: 34110000

3M™ Elastomer Syringe Screw Refill, 71076

3M ID 7000054667 UPC 04035077012973

Details

Replacement screw for the $3M^{\,\text{\tiny TM}}$ Elastomer Syringe

Typical properties

Brands: 3M™

Catalog Number: 71076

Factory ISO Certification: No ISO Certification information available

Package Contents: 1 - Syringe Screw

Product Type: Screw

3M™ Elastomer Standard Spatula Refill, 71233

3M ID 7000054665 UPC 04035077012959

Details

For hand mixing dental materials

Typical properties

Brands: 3M™

Catalog Number: 71233

Factory ISO Certification: No ISO Certification information available

Package Contents: 1 - Spatula

Product Type: Spatula

3M™ Victory Series™ First Molar Bands

Details

Permanent laser ID marking

Straight interproximal walls for reduced separation

Reduced gingival/occlusal dimensions for less occlusal and gingival interferences

Available with micro-etched inner surface for increased retention

Curved festooned occlusal interproximal

Increased occlusal "edge" strength for reduced deformation

Typical properties

Arch: Lower, Upper

Band Size: 131, 131+, 132, 132+, 133, 133+, 134, 134+, 135, 135+, 136, 136+, 137, 137+, 138, 138+, 139, 139+, 140, 140+, 141, 141+, 142, 142+, 143, 143+, 144, 31, 31+, 32, 32+, 33, 33+, 34, 34+, 35, 35+, 36, 36+, 37, 37+, 38, 38+, 39,

39+, 40, 40+, 41, 41+, 42, 42+, 43, 43+, 44

Brands: Victory Series™

Factory ISO Certification: No ISO Certification information available

Material: Stainless Steel Product Type: Band Side: Left, Right Tooth: First Molar

3M[™] TheraMints[™] 100% Xylitol Sweetened Mints, 12119M, Mint Flavor, 520 Pieces Per Bottle

3M ID 7100302144 UPC 10605861081488

Details

Great-tasting mints sweetened with 100% xylitol natural sugar substitute
Do not promote tooth decay
Stimulate saliva flow to fight the harmful effects of dry mouth
Available in mint
Buy TheraMints™ 100% Xylitol Sweetened Mints on Amazon.com/3MOralCare

Typical properties

Active Ingredient: Xylitol Brands: TheraMints™

Factory ISO Certification: No ISO Certification information available

Flavor: Mint Material: Xylitol

3M™ RelyX™ Temp Zinc Oxide Non-Eugenol Temporary Cement Refill

Details

Low film thickness helps ensure an optimal fit Strong adhesion for high retention to the tooth, yet removes easily for final cementation Easy removal of excess cement

Typical properties

Brands: RelyX™

Catalog Number: 56660 Cure Type: Self-cure

Delivery System: Hand Mix

Factory ISO Certification: No ISO Certification information available

Material: Zinc Oxide Non-eugenol

Package Contents: 1 – 36g Base Paste; 1 – 16g Catalyst Paste; 1 Mixing Pad

(70144)

Shade: White Opaque

3M™ RelyX™ Unicem Self-Adhesive Universal Resin Cement Aplicap™ Capsule Kit

3M ID B00007432

Details

Saves time by eliminating etching, priming and bonding steps
Virtually no post-operative sensitivity
Less discoloration
Strong and moisture-tolerant
If you like 3M™ RelyX™ Unicem Self-Adhesive Resin Cement, you'll love 3M™
RelyX™ Universal Resin Cement

Typical properties

Brands: RelyX™

Cure Type: Dual Cure

Delivery System: Aplicap™ Capsule

Material: Self-adhesive Resin

Product Type: Self-adhesive Resin Cement

3M™ Unitek™ Stainless Steel Primary Anterior Crowns

Details

Easy to customize with parallel walls to allow for clinician fit control Pre-trimmed to optimum length and contour to save placement time Made in the U.S.A. of high quality stainless steel Shallow occlusal anatomy means less tooth reduction Thick occlusal surface to prevent bite-through Crown set is packaged in a unique, stackable tray that provides convenient, easy storage of 3M[™] ESPE[™] prefabricated crowns Wide array of shapes and sizes

Typical properties

Arch: Lower, Upper Brands: Unitek™ Catalog Number:

 $907031,907053,907054,907032,907051,907052,907013,907035,907014,907036,\\907011,907033,907055,907034,907056,907012,907015,907016,907042,907043,907021,907041,907002,907046,907024,907003,907025,907022,907044,907001,907045,907023,907006,907004,907026,907005$

Dentition: Primary

Factory ISO Certification: No ISO Certification information available Mesial/Distal Width (Imperial): 0.165 in, 0.181 in, 0.189 in, 0.197 in, 0.205 in, 0.213 in, 0.22 in, 0.228 in, 0.236 in, 0.244 in, 0.252 in, 0.26 in, 0.268 in, 0.276 in, 0.283 in, 0.291 in, 0.299 in, 0.307 in, 0.315 in, 0.323 in Mesial/Distal Width (Metric): 4.2 mm, 4.6 mm, 4.8 mm, 5 mm, 5.2 mm, 5.4 mm,

5.6 mm, 5.8 mm, 6 mm, 6.2 mm, 6.4 mm, 6.6 mm, 6.8 mm, 7 mm, 7.2 mm, 7.4 mm, 7.6 mm, 7.8 mm, 8 mm, 8.2 mm

Number of Crowns: 5 Patient: Pediatric Side: Left, Right Size: 1, 2, 3, 4, 5, 6

Tooth: Central, Central Incisor, Cuspid, Lateral Incisor

3M™ Lava™ Zirconia Blocks Refill, 68583, 20mm

3M ID 7000003303 UPC 04035077004541

Details

High strength for tooth-preserving preparations and all crown and bridge indications

Brilliant esthetics with long-lasting polish

Less wear to opposing enamel than glass ceramics

Lava™ Zirconia substructure is available in 8 shades

Excellent marginal fit

Conventional cementation

Superior strength (1440 MPa, 3-point bend strength)

10 years clinical experience

Typical properties

Block Size: 20 Brands: Lava™

CAD/CAM System: Lava

Digital Materials: Lava Laboratory Materials

Factory ISO Certification: No ISO Certification information available

Material: Zirconia

Restoration Type: Bridge, Crown, Inlay, Onlay, Veneer

Dimensions and Classifications

Overall Height (Imperial): 0.787 in Overall Height (Metric): 20 mm

3M™ Unitek™ MIA Mobile Intraoral Arch, 080-810, U/L, Curved, 1/Ea

3M ID 7000003706

Details

Available in straight or curved sheaths

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Retention: Curved

Units per Inner Pack: 1 per Pack

3M™ RotoMix™ Capsule Mixing Unit, 76310, 120 V/60 Hz

Part Number 76310 3M ID 7100196270 UPC 04054596646629

Details

Easy loading and unloading of capsules

Touchpad features both pre-programmed and programmable setting options for simple operation

Typical mixing speed: 2850 rotations per minute Typical centrifuging speed: 2950 rotations per minute

Extremely quiet operation

Rotating platform with simultaneously counter-rotating capsule holder produces the most complete mix

Automatic centrifugal cycle produces a denser mix and virtually void-free material

Typical properties

Brands: RotoMix™

Factory ISO Certification: No ISO Certification information available

Net Weight (Metric): 4000 g Operating Voltage: 120 V

Package Contents: Includes one capsule mixing unit (120V/60Hz)

Product Type: Capsule Mixing Device eClass 14 Classification Group: 34110000

Dimensions and Classifications

Overall Depth (Metric): 250 mm Overall Height (Metric): 210 mm Overall Width (Metric): 180 mm

3M™ Penta™ Elastomer Nozzle, 71225

3M ID 7100288606 UPC 04064035065485

Details

Allows for easy, intraoral delivery of $3M^{TM}$ Penta impression materials or hand-mixed light and medium bodied materials Compatible with the $3M^{TM}$ Elastomer Syringe and $3M^{TM}$ Penta Elastomer Syringe

Typical properties

Brands: Penta™

Catalog Number: 71225

Factory ISO Certification: No ISO Certification information available

Package Contents: 12 - intraoral tips

Product Type: Tip

3M™ Unitek™ Permachrome Posted Archwire

Details

High tensile strength Stainless steel wires with brass alloy posts Available in folder packs only

Typical properties

Arch: Lower, Upper

Arch Shape: Ovoid, Square

Archwire Size (Imperial): 0.016 x 0.022 in, 0.019 x 0.025 in Archwire Size (Metric): 0.406 x 0.559 mm, 0.483 x 0.635 mm

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Inter-post Spacing (Metric): 3 mm, 22 mm, 24 mm, 26 mm, 28 mm, 30 mm, 32

mm, 34 mm, 36 mm, 38 mm, 40 mm, 42 mm, 44 mm

Material: Stainless Steel

Units per Inner Pack: 10 per Pack

3M™ Transbond™ Plus Light Cure Band Adhesive Kit

Details

Light cure, no-mix adhesive simplifies banding process Blue color for easier initial application and final clean-up Pre-loaded syringe for easy dispensing Efficient, low-waste delivery system Less chair time

Typical properties

Brands: Transbond™

Catalog Number: 712-080,712-081

Cure Type: Light Cure Delivery System: Syringe

Factory ISO Certification: No ISO Certification information available Package Contents: 5 - 4g syringes, 50 disposable angulated tips

Product Type: Light Cure Band Adhesive Kit

Units per Inner Pack: 1 /KT

3M™ ESPE™ Ketac™ Cem Radiopaque Glass Ionomer Luting Cement Intro Kit, 37201

3M ID 7100215212 UPC 04054596733220

Details

Glass ionomer luting cement offers chemical bond to enamel and dentin Releases fluoride High biocompatibility Extremely low film thickness optimizes fit and marginal integrity Radiopaque

Typical properties

Brands: Ketac™

Catalog Number: 37201 Delivery System: Hand Mix

Factory ISO Certification: No ISO Certification information available

Package Contents: 1 – 33g Powder; 1 – 12ml Liquid; 1 Mixing Pad (70144)

Shade: White

3M™ Unitek™ Prestige Hollow Chop Pliers, 900-716, 1/Ea

3M ID 7000025671

Details

Replaceable occlusal pad Replaceable double-ended blade 3M™ Unitek™ Prestige Series

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

3M™ Fiber Tip Applicators, 1919F, 100ct

3M ID 7100213642 UPC 30605861068810

Details

To be used with 3M[™] Applicator Handles

Typical properties

Brands: 3M™

Factory ISO Certification: No ISO Certification information available

Package Contents: Refill includes 100 fiber tip applicators

Product Type: Tip Applicators

3M™ Mixing Pad Refill

Details

Individual mixing sheets are simply discarded after use, for fast and easy clean up

Tear-off, disposable coated sheets allow for easy hand-mixing of dental materials

Typical properties

Brands: 3M™

Catalog Number: 1994P,1961P,70190,70128,70144,70180

Factory ISO Certification: No ISO Certification information available

Package Contents: 10 - Mixing Pads,4 - Mixing Pads

Units per Case: 10.0

eClass 14 Classification Group: 34110110

Dimensions and Classifications

Overall Length (Imperial): 2 in, 3 in, 3.5 in, 6 in, 7.5 in, 9.25 in

Overall Length (Metric): 50.8 mm, 76.2 mm, 88.9 mm, 152.4 mm, 190.5 mm,

234.95 mm

Overall Width (Imperial): 1.5 in, 2.5 in, 3.5 in, 4.75 in, 5.75 in

Overall Width (Metric): 38.1 mm, 6

3M™ Elastomer Syringe Washer Refill, 71050

3M ID 7000054668 UPC 04035077012980

Details

This refill includes one replacement white Teflon ring

Typical properties

Brands: 3M™

Catalog Number: 71050

Factory ISO Certification: No ISO Certification information available

Package Contents: 1 - Syringe Washer

Product Type: Washer

3M™ Victory Series™ Active Self-Ligating Brackets Kit

Details

Enhanced patient comfort
Easy to open and close door
Accurate bracket placement and positioning
Design advantages for treatment flexibility

Typical properties

Adhesive Precoat: Non-coated

Arch: Upper/Lower

Brands: Victory Series™

Factory ISO Certification: No ISO Certification information available

Material: Stainless Steel

Prescription: Roth*

Slot Size (Imperial): 0.022 in Slot Size (Metric): 0.559 mm

Units per Pack: 20

3M™ Silane Coupling Agent, 68300, 1 per Pack

3M ID 7000054803 UPC 04035077001199

Details

Creates a silane layer on the silicate-coated surface of an indirect restoration to produce a chemical bond between an opaquer material like Sinfony™ Opaquer Material, and a veneering material such as Sinfony™ Indirect Lab Composite Proven in tests to be as strong or stronger than other popular brands and to possess inherent advantages in flexibility, impact strength and color stability Particle components are melted on the surface due to their high impact energy and are retentively anchored

Typical properties

Brands: 3M™

Catalog Number: 68300 Delivery System: Bottle

Factory ISO Certification: No ISO Certification information available

Product Type: Dental Primer Units per Inner Pack: 1 per Pack

3M™ Fast Release 5% Sodium Fluoride Varnish, 12150 Series

Details

Contains 22,600 ppm fluoride

Easy all-in-one tray and brush delivery

Virtually clear

Smooth coating

Does not change appearance of ceramic or metal orthodontic brackets

Relieves hypersensitivity

Contains 5% sodium fluoride and xylitol

Three fun and tasty flavors: Caramel, Strawberry and Mint

Typical properties

Active Ingredient: Sodium Fluoride

Brands: 3M™

Flavor: Caramel, Mint, Strawberry

Inner Packs per Case: 5 Packs per Case, 6 Packs per Case

Total Units per Case: 50 per Pack, 200 per Pack Units per Inner Pack: 50 per Pack, 200 per Pack

3M[™] Impregum[™] Penta[™] Super Quick Medium Body/Light Body Polyether Impression Material Intro Kit

Details

Superfast polyether with 45-second working time and quick 2-minute setting time Unique moisture tolerance provides excellent detail reproduction for void-free impressions, even in wet environments

Very high level of accuracy for precise-fitting restorations

Improved patient comfort thanks to fresh, minty taste and shorter setting time Scannable for more efficient workflow

Available in a range of viscosities and delivery systems for monophase and 1step impression techniques

Control moisture at every step and gently retract gingival tissue for more accurate impressions with 3M™ Astringent Retraction Paste

Typical properties

Brands: Impregum™

Catalog Number: 69411,69421,69404,69380

Delivery System: Garant™ Dispenser, Pentamix™ Mixing Unit Factory ISO Certification: No ISO Certification information available

Impression Type: Final Material: Polyether

Package Contents: 2 cartridge Impregum Super Quick MB - 50ml, 1 cartridge Impregum Super Quick LB - 50ml, 1 cartridge Impregum Super Quick ULB - 50ml, 5 3M[™] Intra-oral Syringes purple, 5 3M[™] Garant Mixing Tips purple, 5 3M[™] Garant Intra-oral Tips white, 1 3M[™] Polyether Tray Adhesive - 17ml, 1 base paste - 300ml, 1 catalyst - 60 ml, 10 3M[™] Penta[™] Mixing Tips, Red, 1 base paste - 300ml, 1 catalyst, 60 ml, 1 cartridge Impregum Super Quick LB - 50ml, 10 3M[™] Penta[™] Mixing Tips, Red, 5 3M[™] Intra-oral Syringes purple, 5 3M[™] Garant Mixing Tips purple, 5 3M[™] Garant Intra-oral Tips white

Product Color: Purple

Set Time: 2 min

Set Type: Super Quick Set Setting Time: 2 Minute

Viscosity: Light, Medium, Multiple Viscosities

Work Time: 45 seconds

3M™ Clarity™ Ultra Self-Ligating Brackets Kit

Details

Translucent and stain-resistant, Clarity Ultra brackets are an excellent fullyesthetic choice for you and your patients

Innovative, smooth ceramic bracket surfaces and edges are designed for patient comfort so your patients can live confidently

Clarity™ Ultra Self-Ligating Brackets are passive, but can be activated with the proprietary oval ligation design

A wide ceramic door assists rotational control with a broad mesial-distal span Trust-worthy and reliable door mechanism designed to resist unplanned opening and closing

Unique stress concentrator for predictable squeeze debonding

Clarity Ultra brackets do not require elastomeric ligatures that may stain during treatment

Color-coded with identification dots, vertical and horizontal reference marks for easy bracket selection, and positioning

Compatible with many 3M appliances to give orthodontists custom control and achieve results based on treatment preferences and the unique needs of the patient

Typical properties

Adhesive Precoat: APC™ Flash-Free, Non-coated

Arch: Upper, Upper/Lower

Bracket Style: Twin Brands: Clarity™

Factory ISO Certification: No ISO Certification information available

Hook Included: Yes Hook Type: Distal Material: Ceramic Prescription: MBT[™]

Slot Size (Imperial): 0.022 in Slot Size (Metric): 0.559 mm Units per Pack: 10, 20

3M™ Victory Series™ Brackets Kit

Details

The low profile hook and small ball are milled from a single piece of metal for increased patient comfort

3M offers an innovative color coding system that assigns a different color to each quadrant of the mouth for immediate identification

Gingivally Off-Set Bicuspid Bracket features a three point contact design along with rounded corners for a great base-to-tooth fit

Typical properties

Adhesive Precoat: APC™ II, Non-coated

Arch: Upper/Lower Bracket Style: Twin Brands: Victory Series™

Factory ISO Certification: No ISO Certification information available

Hook Included: Yes Hook Type: Distal

Prescription: Alexander*, MBT™, Roth* Slot Size (Imperial): 0.018 in, 0.022 in Slot Size (Metric): 0.457 mm, 0.559 mm

Units per Pack: 5, 20

3M™ Unitek™ Multi-Cure Glass Ionomer Band Cement, 712-050, 1/Kit

3M ID 7000004384

Details

Moisture and saliva tolerant
Fills voids between tooth surface and bands
Light cure or self cure
Blue tint for easy initial clean-up
Achieve adequate initial strength with light application
Band under moist conditions

Typical properties

Brands: Unitek™

Catalog Number: 712-050

Factory ISO Certification: No ISO Certification information available

Units per Inner Pack: 1 /KT

3M™ Victory Series™ Brackets

Details

The low profile hook and small ball are milled from a single piece of metal for increased patient comfort

3M offers an innovative color coding system that assigns a different color to each quadrant of the mouth for immediate identification

Gingivally Off-Set Bicuspid Bracket features a three point contact design along with rounded corners for a great base-to-tooth fit

Typical properties

Adhesive Precoat: APC™ II, Non-coated

Angulation: +10°, +2°, +3°, +4°, +5°, +6°, +8°, +9°, -4°, -6°, 0°

Arch: Lower, Upper

Bracket Style: Single, Twin Brands: Victory Series™

Factory ISO Certification: No ISO Certification information available

Hook Included: Yes

Hook Type: Distal, Mesial Material: Stainless Steel

Offset Base: Yes

Prescription: Alexander*, Andrews*, Bench, High Torque, Hilgers*, Level

Anchorage™, Low Torque, MBT™, Medium Torque, Ricketts*, Roth*, Samson*,

Sondhi, ZeroBase™

Side: Left, Left/Right, Right

Slot Size (Imperial): 0.018 in, 0.022 in Slot Size (Metric): 0.457 mm, 0.559 mm

Tooth: Anterior, Bicuspid, Central Incisor, Cuspid, First Bicuspid, Lateral Incisor,

Second Bicuspid

Torque: +10°, +14°, +17°, +22°, +3°, +7°, -10°, -11°, -12°, -17°, -1°, -22°, -3°, -5°,

-6°, -7°, 0°

Units per Pack: 5

3M™ Rocatec™ Junior Bonding System Refill

Details

Provides outstanding bond strengths without marginal gap

Compact and economical unit is mobile and can be used wherevercompressed air is available

Does not require a power supply

Easy-to-view pressure gauge allows you to monitor and regulate the unit's performance

Pressure-reduction valve with integral water collector supplies compressed air free of water or oil to ensure reliable operation

Intended for use with the Rocatec Bonding System, including Rocatec Pre Bonding Material, Rocatec Plus Bonding Material and Rocatec Soft Bonding Material

Typical properties

Brands: Rocatec™

Catalog Number: 68360,68340,68350

Digital Materials: Other Laboratory Products

Factory ISO Certification: No ISO Certification information available

3M™ Ortholux™ Luminous Curing Light Eye Shield, 704-457, 1/Ea

3M ID 7100018476 UPC 00652221149381

Details

Orange
Can be disinfected

Typical properties

Brands: Ortholux™

Factory ISO Certification: No ISO Certification information available

Product Type: Eye Shield

3M™ Dual Arch Impression Tray

Details

Used with Polyether or VPS materials with medium or heavy body consistency Not for use with Putty

Captures crown/bridge arch impression, bite registration and opposing arch impression in one

Typical properties

Brands: 3M™

Catalog Number: 71622,71623

Factory ISO Certification: No ISO Certification information available

Package Contents: 1 refill - 35 trays,1 refill - 48 trays

Tray Size: Small

3M™ Chairside Zirconia for CEREC®

Details

Zirconia block for single-unit crown and three-unit bridge restorations*

Optimal blend of high strength and esthetics, with high flexural strength of 1000 MPa (fast fired)

The esthetics you desire from a chairside zirconia, designed to match the VITA classical shade guide

Fast sintering time of approximately 20 minutes for a thin walled crown in CEREC® SpeedFire furnace**,***

Makes it possible to go from scan to seat in one hour Available in two pre-shaded block sizes and eight shades

Typical properties

Block Size: 20, 39 Brands: 3M™

CAD/CAM System: CEREC®

Catalog Number:

69941,69952,69951,69940,69950,69945,69955,69944,69943,69954,69953,6994

2,69949,69948,69947,69946

Digital Materials: Digital Chairside Restorative Materials

Factory ISO Certification: No ISO Certification information available

Material: Zirconia

Restoration Type: Bridge, Crown

Shade: A1, A2, A3, A3.5, B1, BL, C1, D2

Dimensions and Classifications

Overall Diameter (Imperial): 0.787402 in, 1.53543 in

Overall Diameter (Metric): 20 mm, 39 mm Overall Height (Imperial): 0.787 in, 1.535 in Overall Height (Metric): 20 mm, 39 mm

3M™ Unitek™ Prestige Wire Cutter, 900-753, 1/Ea

3M ID 7000025680

Details

Designed to cut soft wires
Sharp pointed diamond-honed tips to cut ligature wire at very point of cutter blades
3M™ Unitek™ Prestige Series

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Tip Shape: Straight

3M™ Unitek™ Prestige Hard Wire Cutter, 900-751, 1/Ea

3M ID 7000025679

Details

Cutting capacity: Maximum archwires up to .021 x .025 in. Hi-T and Nitinol

Cutting capacity: minimum .010 in. diameter Hi-T

3M™ Unitek™ Prestige Series

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Tip Shape: Straight

3M™ Unitek™ Lingual Sheaths

Details

Lingual sheaths not preassembled on bands

Typical properties

Arch: Upper Left/Lower Right, Upper Right/Lower Left, Upper/Lower

Brands: Unitek™

Compatible Wire Size (Imperial): 0.03 in, 0.036 in

Factory ISO Certification: No ISO Certification information available

Hook Included: Yes

Hook Type: Distal, Distal Ball, Mesial Ball, Round Product Type: Lingual Cleat, Lingual Sheaths eClass 14 Classification Group: 34110603

3M™ Transbond™ Plus Disposable Tips, 712-082, 25/Pack

3M ID 7100154492

Details

Disposable 25/pack

Typical properties

Brands: Transbond™ Catalog Number: 712-082 Cure Type: Light Cure

Factory ISO Certification: No ISO Certification information available

Units per Inner Pack: 25 per Pack

3M™ Transbond™ Plus Self Etching Primer

Details

One-step etch and prime system

High performance in moist or dry environments

Each single use foil pack contains enough primer for one complete arch Added patient comfort because of the elimination of the etchant rinsing process Exclusive delivery system is easily mixed by hand or used with the 3M[™] Transbond[™] Easy Roller

Use with moisture-tolerant light cure adhesives for good bond strength A shorter treatment window saves you valuable time and builds greater productivity for your practice (1)

Patient comfort is increased because of the elimination of the etchant rinsing process

Performs both etching and penetration of monomers into the demineralized enamel, in one step, resulting in less enamel loss when compared to traditional etching

Typical properties

Application: Primer Brands: Transbond™

Factory ISO Certification: No ISO Certification information available

Package Contents: 100/Pk,20/Pk

3M[™] Unitek[™] Lingual Heavy Wire Flush Cut Distal End Cutter, 803-551, 1/Ea

3M ID 7100267823

Details

Used for a wide range of wire sizes
Long handles and a narrow design
Cutting tips made of Tungsten Carbide for durability
Clean square end cut reduces deburring
Safety hold feature secures the loose wire after the archwire is cut

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

3M™ Unitek™ Gurin Intraoral Hook

Details

Anchor without welding or soldering Adjustable

Typical properties

Brands: Unitek™

Compatible Archwire Size (Imperial): 0.014 - 0.016 in, 0.018 - 0.022 in, 0.038 -

0.045 in

Factory ISO Certification: No ISO Certification information available

3M™ Sondhi™ Rapid-Set Indirect Bonding Adhesive Resin B, 712-073, 10 ml, 1/Btl

3M ID 7000025773

Details

Less bracket repositioning More control and accurate bracket placement second molar to second molar No mixing of Resin A & B required

Typical properties

Brands: Sondhi™

Catalog Number: 712-073

Container Size (Imperial): 0.338 fl oz

Container Size (Metric): 10 mL

Delivery System: Bottle

Factory ISO Certification: No ISO Certification information available

Product Type: Indirect Bonding Adhesive Kit

Units per Inner Pack: 1 /BT

3M™ Unitek™ Permachrome Standard Archwire

Details

Better malleability than Permachrome Resilient and Hi-T™ II Lowest tensile strength

Typical properties

Arch: Lower, Upper Arch Shape: Standard

Archwire Size (Imperial): 0.016 x 0.022 in, 0.017 x 0.025 in, 0.018 in, 0.018 x

0.025 in, 0.019 x 0.025 in

Archwire Size (Metric): 0.406 x 0.559 mm, 0.432 x 0.635 mm, 0.457 mm, 0.457 x

0.635 mm, 0.483 x 0.635 mm

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Material: Stainless Steel

Units per Inner Pack: 10 per Pack, 50 per Pack eClass 14 Classification Group: 34110604

3M™ Transbond™ XT Dispensing Gun, 712-032, 1/Ea

3M ID 7100154489

Details

Fits all 3M Adhesive capsules Sterilizable/reusable Easy loading and dispensing of capsules

Typical properties

Brands: 3M™

Catalog Number: 712-032

Factory ISO Certification: No ISO Certification information available

3M™ Forsus™ Fatigue Resistant Device Kits

Details

Esthetic option for patients with no headger

Increased versatility with three installation options

One appointment installation, chairside, with no lab work required for an out-of-the-box solution

A simple one-step installation into the headgear tube

A telescoping coaxial spring designed to allow for enhanced patient jaw opening

A comfortable stainless steel mechanism

Efficient treatment that works 24/7

Open spring design to facilitate good oral hygiene

Independent left and right components; use unilaterally if desired

L-pin Module for flexible, more mobile installation

EZ2 Module for simple, secure installation

Use in the headgear tube of buccal tubes welded to bands or with the Wire

Mount option for buccal tubes bonded to molars

Typical properties

Brands: Forsus™

Factory ISO Certification: No ISO Certification information available

Module Style: EZ2, L-Pin

Product Type: EZ2 Module Kit, L-pin Module Kit

Side: Left

Units per Inner Pack: 1 /KT

eClass 14 Classification Group: 34110690

Dimensions and Classifications

Overall Length (Metric): 22 mm, 25 mm, 29 mm, 32 mm, 35 mm, 38 mm

3M[™] Sondhi[™] Rapid-Set Indirect Bonding Adhesive Resin A, 712-072, 10 ml, 1/Btl

3M ID 7000025772

Details

Less bracket repositioning More control and accurate bracket placement second molar to second molar No mixing of Resin A & B required

Typical properties

Brands: Sondhi™

Catalog Number: 712-072

Container Size (Imperial): 0.338 fl oz

Container Size (Metric): 10 mL

Delivery System: Bottle

Factory ISO Certification: No ISO Certification information available

Product Type: Indirect Bonding Adhesive Kit

Units per Inner Pack: 1 /BT

3M™ Unitek™ Beta Titanium Lateral Development Archwire

Details

Available in three radiuses

Typical properties

Archwire Size (Imperial): 0.012 in, 0.013 in, 0.014 in, 0.014 x 0.025 in, 0.016 in, 0.016 x 0.022 in, 0.016 x 0.025 in, 0.017 x 0.025 in, 0.018 in, 0.018 x 0.025 in, 0.019 x 0.025 in, 0.02 in

Archwire Size (Metric): 0.305 mm, 0.33 mm, 0.356 mm, 0.356 x 0.635 mm, 0.406 mm, 0.406 x 0.559 mm, 0.406 x 0.635 mm, 0.432 x 0.635 mm, 0.457 mm, 0.483 x 0.635 mm, 0.508 mm

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Material: Beta III Titanium, Nitinol, Stainless Steel Radius (Imperial): 0.945 in, 1.024 in, 1.102 in Radius (Metric): 24 mm, 26 mm, 28 mm

Stops Included: Yes

Units per Inner Pack: 10 per Pack

3M™ RelyX™ Luting Plus Automix Tips, 3540

3M ID 7100111950 UPC 04035077014359

Details

For easy, automatic mixing of 3M[™] RelyX[™] Luting Plus Automix cement Automix delivery allows for direct delivery into the crown (no hand mixing) Virtually no post-operative sensitivity Resin Modified Glass Ionomers (RMGI)

Typical properties

Brands: RelyX™

Catalog Number: 3540

Factory ISO Certification: No ISO Certification information available

Package Contents: 12 Mixing Tips

Product Type: Tip

3M™ Unitek™ Lingual Ball Hooks

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available eClass 14 Classification Group: 34110603

3M™ AlastiK™ EZ3 Holder, 406-044, 10/Pk

3M ID 7100266647

Details

Not made with natural rubber latex
Pre-bonding lingual wire support
Eliminates need for wax, ligature wire, floss threader or manual assistance

Typical properties

Brands: AlastiK™

Factory ISO Certification: No ISO Certification information available

Units per Inner Pack: 10 per Pack

3M™ Unitek™ Preformed Ligature Wire

Details

Dead soft permachrome wire in color coded bundles of 200

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Units per Inner Pack: 100 per Pack

eClass 14 Classification Group: 34110603

Dimensions and Classifications

Overall Diameter (Imperial): 0.008 in, 0.009 in, 0.01 in, 0.011 in, 0.012 in, 0.014

in

Overall Diameter (Metric): 0.01 mm

2600-021 APC™ PLUS Adhesive Inventory Dispensing System Drawer 2/Pk

3M ID 7100154383

Details

Allows first-in/first-out inventory dispensing Each drawer insert has 10 channels Contains label set for APC Adhesive identification

Typical properties

Brands: 3M™

Factory ISO Certification: No ISO Certification information available

3M™ Unitek™ Convertible Cap Removing Pliers, 805-035, 1/Ea

3M ID 7000004387

Details

Replacement tips available

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

3M™ Unitek™ Prestige Three-Jaw Wire Bending Pliers, 900-705, 1/Ea

3M ID 7000025660

Details

Forms terminal double-back bends 3M™ Unitek™ Prestige Series

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

3M™ Unitek™ Molar Hooks 480-001, Ea

3M ID 7000025566

Details

Easy to weld on bands Strong enough for heavy cross elastics No sharp edges to irritate tissue or damage elastics

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

3M™ Mathieu Needle Holder, 500-083, 1/Ea

3M ID 7100272490

Details

Precision stainless steel ligating instrument with fine, serrated tips Quick-release positive locking mechanism on the rear of the handles Freesliding leaf spring on the inner surface opens the beaks when the lock is released Useful for placing 3M™AlastiK Force Modules Can also be used with metal ligatures

Typical properties

Brands: 3M™

Factory ISO Certification: No ISO Certification information available

3M™ Unitek™ Action Line Thread

Details

Not made with natural rubber latex Clear polyurethane elastomeric thread Self-sealing bobbin

Typical properties

Arbor Diameter (Imperial): 0.02 in, 0.025 in, 0.03 in

Brands: Forsus™

Factory ISO Certification: No ISO Certification information available

Units per Inner Pack: 1 Each

eClass 14 Classification Group: 34110690

Dimensions and Classifications

Overall Length (Imperial): 25 ft

3M™ Unitek™ Gurin Retraction Loop

Details

Anchor without welding or soldering Adjustable

Typical properties

Brands: Unitek™

Compatible Archwire Size (Imperial): 0.018 - 0.022 in Factory ISO Certification: No ISO Certification information available

3M™ Friction Stop Adjustable Facebows

Details

1/pack

Typical properties

Brands: 3M™

Factory ISO Certification: No ISO Certification information available

Inner Bow Size: Universal Outer Bow Angulation: 0°

Outer Bow Size: Long, Short, Super Short

Dimensions and Classifications

Overall Diameter (Imperial): 0.045 in, 0.051 in

3M™ Unitek™ Prestige Slimnose Weingart Utility Pliers, 900-735, 1/Ea

3M Product Number 800-423-4588 3M ID 7000025685

Details

Double-ended instrument One curved, tapered tip with sharpened inside radii One tip flat serrated to push bracket

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

3M[™] Filtek[™] Bulk Fill Posterior Restorative, 4863A3, 1 - 4g Syringe, A3

3M ID 7100054399 UPC 30605861000957

Details

One increment up to 5 mm

Excellent handling and sculptability

No expensive dispensing devices or time-consuming layers

Posterior Restorations made simple

Easy adaptation

Excellent wear resistance

Want to provide your patients with both long-lasting performance and uncompromising esthetics? Explore 3M™ Filtek One Bulk Fill Restorative Want to provide your patients with both long-lasting performance and uncompromising esthetics? Explore 3M™ Filtek One Bulk Fill Restorative

Typical properties

Brands: Filtek™

Catalog Number: 4863A3 Composite Type: Bulk Fill Cure Type: Light Cure Delivery System: Syringe

Factory ISO Certification: No ISO Certification information available

Inner Packs per Case: 10 Bags per Case

Net Weight (Metric): 4 g

Shade: A3

Total Units per Case: 10 per Case Units per Inner Pack: 1 per Pack

Units per Pack: 1

3M™ Filtek™ Bulk Fill Posterior Restorative, 4863A2, 1 - 4g Syringe, A2

3M ID 7100067363 UPC 30605861000933

Details

One increment up to 5 mm

Excellent handling and sculptability

No expensive dispensing devices or time-consuming layers

Posterior Restorations made simple

Easy adaptation

Excellent wear resistance

Want to provide your patients with both long-lasting performance and uncompromising esthetics? Explore 3M™ Filtek One Bulk Fill Restorative Want to provide your patients with both long-lasting performance and uncompromising esthetics? Explore 3M™ Filtek One Bulk Fill Restorative

Typical properties

Brands: Filtek™

Catalog Number: 4863A2 Composite Type: Bulk Fill Cure Type: Light Cure Delivery System: Syringe

Factory ISO Certification: No ISO Certification information available

Inner Packs per Case: 10 Bags per Case

Net Weight (Metric): 4 g

Shade: A2

Total Units per Case: 10 per Case Units per Inner Pack: 1 per Pack

Units per Pack: 1

3M™ Disposable Dispensing Tips

Details

Tips available in metal or plastic

Typical properties

Brands: 3M™

Factory ISO Certification: No ISO Certification information available Material: Metal, Plastic

3M™ Unitek™ Prestige End Cutter, Distal, 900-750, 1/Ea

3M ID 7000025678

Details

15° cutting angle
Diamond-honed to cut hard wire
Precision-machined joint
3M™ Unitek™ Prestige Series

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Tip Shape: Straight

3M™ AlastiK™ Elastomeric Ligature with Guard

Typical properties

Brands: AlastiK™

Factory ISO Certification: No ISO Certification information available

Product Color: Silver, Transparent Units per Inner Pack: 100 per Pack

3M[™] Unitek[™] Ligature Ties, 390-170, Preformed Coated, Tooth Colored, .012 in., 100/Pk

3M ID 7100154543

Details

Esthetic, tooth colored ligature tie Stainless Steel tie with a Polymer Coating (Ethylene/Tetrafluoroethylene Copolymer) .012 diameter

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Units per Inner Pack: 100 per Pack

eClass 14 Classification Group: 34110603

Dimensions and Classifications

Overall Diameter (Imperial): 0.012 in

3M[™] APC[™] Inventory Dispensing System Drawer, One Drawer Set, 2600-014, 1/Ea

3M ID 7100240300

Details

Allows first-in/first-out inventory dispensing
Drawer insert has 10 channels
Contains complete set of labeling for drawer content identification
Can be used to store U/L, left/right APC molar blisters

Typical properties

Brands: 3M™

Factory ISO Certification: No ISO Certification information available

3M[™] Unitek[™] Prestige Tweed Rectangular Arch-Forming Pliers, 900-703, 1/Ea

3M ID 7000025658

Details

For precise bending of lighter wires
Precision-aligned grooves with serrated tips in both beaks assure accurate bending of closing loops
3M™ Unitek™ Prestige Series

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

3M™ Unitek™ Lingual Button, 20400-02, Curved Bondable, 10/Pk

3M ID 7100154382

Details

10/pack

Typical properties

Arch: Upper/Lower Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

3M™ Unitek™ Bracket Positioning Gauges, 900-841,One - Four, 1/Kit

3M ID 7100272489

Details

Hinged instrument Ultra-fine tips

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

3M™ Forsus™ Push Rod

Details

Individual components available for re-order For 5-patient and 20-patient kits

Typical properties

Brands: Forsus™

Factory ISO Certification: No ISO Certification information available

Force Rating: Medium

Product Type: Fatigue Resistant Device

Side: Left, Right

Units per Inner Pack: 5 per Pack

eClass 14 Classification Group: 34110690

Dimensions and Classifications

Overall Length (Imperial): 0.984 in, 1.142 in, 1.26 in, 1.378 in

Overall Length (Metric): 22 mm, 25 mm, 29 mm, 32 mm, 35 mm, 38 mm

3M™ Unitek™ Gurin Extraoral Hook

Details

Anchor without welding or soldering Adjustable

Typical properties

Brands: Unitek™

Compatible Archwire Size (Imperial): 0.018 - 0.022 in Factory ISO Certification: No ISO Certification information available

Side: Left, Right

3M™ Friction Fit Stops, .045 in, 519-562, 100/Pk

3M ID 70202083120

Details

Springy for easy fit Solder in place Inside diameter .045 in. (1.14 mm)

Typical properties

Brands: 3M™

Factory ISO Certification: No ISO Certification information available

Dimensions and Classifications

Overall Diameter (Imperial): 0.045 in

#N/A

Details

Easy to incorporate Class Two Correction mid-treatment Molar bands no longer required T-hook design adds stability Universal for left and right use

Typical properties

Brands: Forsus™

Compatible Archwire Size (Imperial): 0.018 in, 0.022 in

Factory ISO Certification: No ISO Certification information available

Product Type: Fatigue Resistant Device Units per Inner Pack: 10 per Pack

3M[™] Ortho Laboratory Adhesive for Indirect Bonding, 704-050, 3 Syringes/Pk

3M ID 7000003355

Details

Water-soluble

Used for the temporary bonding of orthodontic bases to dry plaster or stone model

Typical properties

Brands: 3M™

Catalog Number: 704-050 Delivery System: Syringe

Factory ISO Certification: No ISO Certification information available

Product Type: Orthodontic Adhesive Units per Inner Pack: 3 per Pack

3M™ Unitek™ Etching Gel, 704-060, Blue Tint, 9 ml, 1/Btl

3M ID 7000026050 UPC 00652221108296

Details

9 ml bottle

Typical properties

Brands: Unitek™

Delivery System: Bottle

Factory ISO Certification: No ISO Certification information available

3M™ Unitek™ Nitinol Heat-Activated Lateral Development Archwire

Details

Can be cooled or chilled resulting in a softer, more pliable wire for easy engagement

Provides light continuous forces

Force activation occurs around body temperature

Good bend retention when bent in hot water

Available in three radiuses

Available with pre-loaded stops

Typical properties

Archwire Size (Imperial): 0.014 x 0.025 in, 0.016 x 0.025 in, 0.017 x 0.025 in,

0.018 x 0.025 in, 0.019 x 0.025 in

Archwire Size (Metric): 0.356 x 0.635 mm, 0.406 x 0.635 mm, 0.432 x 0.635 mm,

0.457 x 0.635 mm, 0.483 x 0.635 mm

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Material: Nitinol

Radius (Imperial): 0.945 in, 1.024 in, 1.102 in Radius (Metric): 24 mm, 26 mm, 28 mm

Stops Included: Yes

Units per Inner Pack: 10 per Pack

3M™ ESPE Concise Light Cure White Sealant Resin, 1930W, 1 Each/Pack

3M ID 7000054172 UPC 30605861032231

Details

Light-cure pit and fissure sealant Easy application due to white coloring Unfilled

Typical properties

Brands: Concise™

Factory ISO Certification: No ISO Certification information available

Inner Packs per Case: 10 Packs per Case Total Units per Case: 10 per Case Units per Inner Pack: 1 per Pack

3M™ Clarity™ Aligners-Flex + Force

Details

Proprietary multi-layer copolymer extruded in two thicknesses to express different mechanical force levels in a custom treatment design

Offered in one seamless digital workflow for 3M™ Clarity™ Aligners, 3M™

Clarity™ Retainers, 3M™ Digital Bonding and combination Tx designs

Excellent initial and sustained force persistence

Straight trimline contributes to torque control and retentive forces to facilitate movement

Highly durable and scratch-resistant

Nearly invisible and more stain-resistant than Invisalign® aligners

Typical properties

Brands: Clarity™

Factory ISO Certification: No ISO Certification information available

3M™ AlastiK™ Easy-To-Tie Ligatures Starter Kit

Typical properties

Brands: AlastiK™

Factory ISO Certification: No ISO Certification information available

Product Color: Assorted

Units per Inner Pack: 200 per Pack

3M™ Unitek™ Prestige Miniature Pin and Ligature Cutter, 900-754, 1/Ea

3M ID 7000025682

Details

Sharp pointed diamond-honed tips to cut ligature wire at very point

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Units per Inner Pack: 1 per Pack

3M™ Unitek™ Prestige Bird Beak Pliers, 900-709, 1/Ea

3M ID 7000025664

Details

Beaks are gradually tapered to bend small diameter loops 3M™ Unitek™ Prestige Series

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

3M[™] Unitek[™] Prestige Posterior Band Removing Pliers, Replacement Pads, 1/4 in., 900-715, 4/Pk

3M ID 7000025670

Details

Ideal for reshaping and contouring Nitinol wire and arches 3M™ Unitek™ Prestige Series

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

3M™ Unitek™ Seating Lugs, 55264, 1/Ea

3M ID 7000000085

Details

Standard or custom lingual position

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

3M™ Unitek™ Pedodontic Molar Bands Kit

Details

Available prewelded with 3M[™] Victory Series[™] Buccal Tubes Universal anatomy fits rights and lefts Regular temper Narrow occlusal-gingival height fits primary molars

Typical properties

Arch: Lower, Upper, Upper/Lower

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Material: Stainless Steel

Product Type: Introduction Kit, Spectrum Kit, Trial Kit

Tooth: Molar

3M™ Facebow Deluxe Carrying Case, 4x9 in, Blue, 706-004, 10/Pk

3M ID 70202089002

Details

Vinyl Zipper with pull tab

Typical properties

Brands: 3M™

Factory ISO Certification: No ISO Certification information available

3M™ Unitek™ Modular Bond Kit Organizer

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

3M™ Elipar™ DeepCure-S LED Curing Light

Details

A predictable, reliable cure, even at the bottom of the proximal box Significantly better depth of cure even when light positioning is difficult Better light uniformity and intensity distribution at clinically relevant distances Updated light guide geometry allows easy access to all tooth surfaces The angle of the light guide has been redesigned to make it more comfortable for patients even if they have limited mouth-opening abilities One-piece, high-quality stainless steel casing is robust and easy to disinfect Simple two-button operation with one cure mode for all indications Charging base with built-in light intensity meter and battery charge indicator

Typical properties

Brands: Elipar™

Catalog Number: 76973 Curing Light Type: LED

Factory ISO Certification: No ISO Certification information available

Light Intensity: 1.47 mW/cm²

Maximum Light Wavelength: 480 nm Minimum Light Wavelength: 430 nm

Net Weight (Imperial): 250 g Number of Curing Modes: 6

Package Contents: 1 - Elipar DeepCure-S LED Curing Light

Product Type: Adapter, Curing Light, Light Guide

Voltage: 110 V - 240 V, 120 V

3M™ RelyX™ Veneer Cement Kit

Details

Light-cure formula
Available in customer-preferred delivery and shading system
Excellent color stability
Easy handling and clean-up
Easy to dispense, apply and seat
Ensure precise and convenient color-matching with 3M™ RelyX™ Try-In Paste

Typical properties

Brands: RelyX™

Catalog Number: 8716TK,8718,8716

Cure Type: Light Cure Delivery System: Syringe

Factory ISO Certification: No ISO Certification information available

Inner Packs per Case: 10 Packs per Case

Material: Adhesive Resin

Package Contents: 1 – 3g Syringe 3M[™] RelyX[™] Veneer Cement in Shades: A1/Light Yellow, A3 Opaque/Yellow Opaque, White Opaque, Translucent, B0.5/White; 1 – 2g Syringe 3M[™] RelyX[™] Try-In Paste in Shades: A1/Light Yellow, A3 Opaque/Yellow Opaque, White Opaque, Translucent, B0.5/White; 1 – 5ml Vial 3M[™] Scotchbond[™] Universal Adhesive; 1 – 3ml Syringe 3M[™] Scotchbond[™] Universal Etchant; Disposable Mixing Wells; Microbrushes; Etchant Dispensing Tips,1 – 3g Syringe 3M[™] RelyX[™] Veneer Cement Shade Translucent; 1 – 2g Syringe 3M[™] RelyX[™] Veneer Try-In Paste Shade Translucent; 20 Unit Dose 3M[™] Scotchbond[™] Universal Adhesive; 1 – 3ml Syringe 3M[™] Scotchbond[™] Etchant; 25 – 3M[™] Scotchbond[™] Etchant Dispensing Tips; Technique Cards; 3M[™] RelyX[™] Veneer Shade Guide

Shade: Multiple

Total Units per Case: 10 per Case Units per Inner Pack: 1 per Pack

3M™ VPS Tray Adhesive 17 ml Refill, 7307

3M ID 7100196366 UPC 04054596648609

Details

Ensures adhesion between the VPS impression material and the impression tray

Typical properties

Adhesive Type: Self-etch

Brands: 3M™

Catalog Number: 7307 Cure Type: Air Dry

Factory ISO Certification: No ISO Certification information available

Material: VPS

Product Type: Tray Adhesive

3M™ Express™ VPS Impression Material Intro Kit, 7300

3M ID 7000054277 UPC 30605861036055

Details

Hydrophilicity and flow behavior helps result in void-free impressions
Elastic recovery helps prevent distortion for a more accurate impression
Resists tearing upon removal thanks to material toughness
Recovery from deformation helps result in fewer retakes
Available in a full range of viscosities, deliveries and setting times
Control moisture at every step and gently retract gingival tissue for more accurate impressions with 3M™ Astringent Retraction Paste

Typical properties

Brands: Express™ Catalog Number: 7300

Delivery System: Multiple Deliveries

Factory ISO Certification: No ISO Certification information available

Impression Type: Final

Material: VPS

Package Contents: 1 - jar of base, 305ml; 1 - jar of catalyst, 305 ml; 2 - 50ml cartridges Light Body regular set; 20 - 3M™ Garant™ Mixing Tips - Yellow; 10-3M™ Garant™ Intraoral Tips - Yellow; 1 - 3M™ VPS tray adhesive 17ml; 2 -

putty spoons; 1 - 3M™ Garant™ Dispenser 1:1/2:1

Set Time: 5 min Set Type: Regular Set

Setting Time: 5 Minute

Viscosity: Multiple Viscosities

Work Time: 1 min

3M[™] Express[™] XT Putty Quick VPS Impression Material, 36973, Putty, Refill

3M ID 7100011745 UPC 04035077003933

Details

Hydrophilicity and flow behavior helps result in void-free impressions
Elastic recovery helps prevent distortion for a more accurate impression
Resists tearing upon removal thanks to material toughness
Recovery from deformation helps result in fewer retakes
Available in a full range of viscosities, deliveries and setting times
Control moisture at every step and gently retract gingival tissue for more accurate impressions with 3M™ Astringent Retraction Paste

Typical properties

Brands: Express™

Catalog Number: 36973

Factory ISO Certification: No ISO Certification information available

Impression Type: Final

Material: VPS

Set Type: Quick Set

Viscosity: Multiple Viscosities

3M™ Impregum™ Penta™ Cartridge for Pentamix™ 3

Details

For consistent, homogeneous and void-free mixing at the touch of a button Compatible with 3M™ Pentamix™ 3 Mixing Unit

Reusable, metal-reinforced cartridge

Helps get an accurate impression on the first take, reducing the need for costly adjustments, remakes and retakes

Typical properties

Brands: Impregum™

Catalog Number: P3794,69401,P3792,P3788,69387,P3776,P3796,P3774

Delivery System: Pentamix™ Mixing Unit

Factory ISO Certification: No ISO Certification information available

Package Contents: 1 - Penta™ Cartridge,1 3M™ Penta™ Cartridge for the

Pentamix™ 3

Product Type: Cartridge

3M™ Unitek™ Stainless Steel Permanent Molar Crowns

Details

Easy to customize with parallel walls to allow for clinician fit control Pre-trimmed to optimum length and contour to save placement time Made in the U.S.A. of high quality stainless steel Shallow occlusal anatomy means less tooth reduction Thick occlusal surface to prevent bite-through Available in molar styles to meet your adult patient needs Crown set is packaged in a unique, stackable tray that provides convenient, easy storage of 3M™ ESPE™ prefabricated crowns Wide array of shapes and sizes

Typical properties

Arch: Lower, Upper Brands: Unitek™ Catalog Number:

900316,900436,900315,900337,900435,900336,900314,900335,900313,900434,900433,900334,900312,900333,900432,900311,900332,900431,900331,900327,900326,900325,900446,900347,900346,900445,900324,900345,900444,900323,900344,900443,900322,900442,900343,900321,900342,900441,900341,900317

Dentition: Permanent

Factory ISO Certification: No ISO Certification information available Mesial/Distal Width (Imperial): 0.37 in, 0.378 in, 0.39 in, 0.398 in, 0.409 in, 0.421 in, 0.429 in, 0.449 in, 0.461 in, 0.472 in, 0.476 in, 0.48 in, 0.492 in, 0.512 in Mesial/Distal Width (Metric): 9.4 mm, 9.6 mm, 9.9 mm, 10.1 mm, 10.4 mm, 10.7 mm, 10.9 mm, 11.4 mm, 11.7 mm, 12 mm, 12.1 mm, 12.2 mm, 12.5 mm, 13 mm Number of Crowns: 5

Patient: Pediatric Side: Left, Right Size: 2, 3, 4, 5, 6

Tooth: First Molar, Second Molar

3M™ Imprint™ 4 Penta™ Cartridge for Pentamix™ 3, 77780 Series

Details

For consistent, homogeneous and void-free mixing at the touch of a button Compatible with 3M™ Pentamix™ 3 Mixing Unit Reusable, metal-reinforced cartridge

Helps get an accurate impression on the first take, reducing the need for costly adjustments, remakes and retakes

Typical properties

Brands: Imprint™

Catalog Number: 77785,77786

Factory ISO Certification: No ISO Certification information available

Package Contents: 1 - Penta™ Cartridge

Product Type: Cartridge

3M™ Unitek™ Ligature Wire

Details

Preformed ligature slips easily over all brackets
Dead soft permachrome wire in color coded bundles of 200
Coated ligature wire reduces friction between archwire and bracket

Typical properties

Brands: Unitek™

Dispenser Type: Spool

Factory ISO Certification: No ISO Certification information available

Units per Inner Pack: 1 per Pack

eClass 14 Classification Group: 34110603

Dimensions and Classifications

Overall Diameter (Imperial): 0.008 in, 0.009 in, 0.01 in, 0.011 in, 0.012 in, 0.014

3M™ RelyX™ Fiber Post 3D Intro Kit, 56957

3M ID 7100288485 UPC 04064035073206

Details

Higher mechanical retention to the core build-up material thanks to coronal 3D macro retentions

Strong bond in the root canal due to its microporous surface without a silanating step

Better post position control thanks to improved radiopacity that is among the highest of available glass fiber posts

Low risk of root fracture due to dentin-like elasticity

More esthetic option than metal post

Excellent light translucency for efficient light cure deep into the root canal

Typical properties

Brands: RelyX™

Catalog Number: 56957

Factory ISO Certification: No ISO Certification information available

Material: Fiber Post

Package Contents: 5 Posts of each size: 0 (1.1mm dia.), 1 (1.3mm dia.), 2 (1.6mm dia.), 3 (1.9mm dia.); 1 Drill of each size: 0, 1, 2, 3, Universal

Product Type: Fiber Post 3D

3M™ Disposable Applicator Brush Tips, 1919B, 60ct

3M ID 7100214035 UPC 30605861068773

Details

To be used with 3M[™] Applicator Handles

Typical properties

Brands: 3M™

Factory ISO Certification: No ISO Certification information available

Package Contents: Refill includes 60 brush tips

Product Type: Brush Tip

3M™ Mixing Sticks (50), 1994S

3M ID 7000054184 UPC 30605861032873

Details

Allows for easy mixing of $3M^{\text{TM}}$ ESPETM dental materials Disposable for improved infection control

Typical properties

Brands: 3M™

Catalog Number: 1994S

Factory ISO Certification: No ISO Certification information available

Package Contents: Refills are available in packages of 50

Product Type: Stick

3M™ Transbond™ Plus Color Change Adhesive

Details

Moisture tolerance for added confidence in bonding when moisture control is difficult

Colored adhesive for enhanced bracket positioning and easy excess cleanup prior to curing

After light curing the color fades away

Available in easy-to-use syringes and capsules

Bond strength comparable to 3M[™] Transbond[™] XT Light Cure Adhesive Efficient, low-waste packaging – reorder bonding components only as needed

Typical properties

Adhesive Type: Direct Brands: Transbond™

Catalog Number: 712-101,712-103,712-102,712-105,712-104,712-107,712-

106,712-108

Cure Type: Light Cure

Delivery System: Capsule, Syringe

Factory ISO Certification: No ISO Certification information available Product Type: Color Change Adhesive Kit, Orthodontic Adhesive

Units per Inner Pack: 1 /KT, 1 per Pack, 2 per Pack, 4 per Pack, 25 per Pack

3M™ Transbond™ MIP Moisture Insensitive Primer Kit, 712-021, 1/Kit

3M ID 7000025650

Details

Moisture tolerant without compromising bond strength Compatible with all 3M light cure adhesives Quick and easy application with just one coat needed Helps to reduce chair time and increase office efficiency

Typical properties

Brands: Transbond™
Catalog Number: 712-021
Cure Type: Light Cure
Delivery System: Bottle

Factory ISO Certification: No ISO Certification information available

Product Type: Moisture Insensitive Primer Kit

Units per Inner Pack: 1 /KT

3M™ Victory Series™ V-Slot Brackets

Details

Full contoured mesh base provides precision placement and consistent bonds The low profile torque-in-base design achieves level slot line-up and easy bracket placement

Chamfered .018 x .018 V-slot allows easy insertion and removal of various auxiliaries to aid tooth alignment

Typical properties

Adhesive Precoat: APC™ II, APC™ PLUS, Non-coated

Angulation: +10°, +5°, +6°, +8°, 0°

Arch: Lower, Upper Base Type: Weldable Brands: Victory Series™

Factory ISO Certification: No ISO Certification information available

Hook Included: Yes Hook Type: Distal Offset Base: Yes Prescription: Roth*

Side: Left, Left/Right, Right

Slot Size (Imperial): 0.018 in, 0.022 in Slot Size (Metric): 0.457 mm, 0.559 mm

Tooth: Anterior, Bicuspid, Central Incisor, Cuspid, First Bicuspid, Lateral Incisor,

Second Bicuspid

Torque: +14°, +17°, +1°, +22°, +7°, -17°, -1°, -22°, -3°, -7°

Units per Pack: 5

3M™ Victory Series™ Active Self-Ligating Brackets

Details

Enhanced patient comfort

Easy to open and close door

Accurate bracket placement and positioning

Design advantages for treatment flexibility

Easy to open & close doors can be opened using either the U-notch or the gingival tab - no special instruments needed

Rounded slot edges designed to reduce archwire binding

Ample under tie-wing area for use of elastic chains and other ligated auxiliaries

Compatible with other active self-ligating bracket archwire sequences

Heat-treated Nickel-Cobalt alloy door for secure, reliable, and robust ligation

Full slot-width door for optimal rotational control

Intricate based features micro-etched, 80-gauge mesh pad applied at a 45° angle

for greater bond strength

Optional APC™ Adhesive Coated Appliance System

Typical properties

Adhesive Precoat: Non-coated Angulation: +10°, +5°, +7°, +9°, 0°

Arch: Lower, Upper Brands: Victory Series™

Factory ISO Certification: No ISO Certification information available

Hook Included: Yes Hook Type: Distal

Material: Stainless Steel

Prescription: Roth*

Side: Left, Left/Right, Right Slot Size (Imperial): 0.022 in Slot Size (Metric): 0.559 mm

Tooth: Anterior, Bicuspid, Central Incisor, Cuspid, First Bicuspid, Lateral Incisor,

Second Bicuspid

Torque: +12°, +8°, -11°, -17°, -1°, -22°, -2°, -7°

Units per Pack: 5

3M™ Victory Series™ Superior Fit Buccal Tubes

Details

Large compound contour base for superior fit
"Contoured funnel" design for easy wire insertion
Sleek, low profile, tapered body and flush-mount hook for patient comfort
Parallel mesial and distal surfaces for ease of gripping and handling
Micro-etched, 80-gauge mesh pad for reliable bonding

Typical properties

Adhesive Precoat: APC™ Flash-Free, APC™ II, Non-coated

Arch: Lower, Upper Brands: Victory Series™

Distal Offset: 0 Degree, 6 Degree, 7 Degree, 10 Degree

Factory ISO Certification: No ISO Certification information available

Hook Included: Yes Hook Type: Distal

Material: Stainless Steel Prescription: MBT™, Roth*

Side: Left, Right

Slot Size (Imperial): 0.018 in, 0.022 in Slot Size (Metric): 0.457 mm, 0.559 mm

Tooth: First Molar, Second Molar

Torque: -10°, -14°, -15°, -20°, -25°, -30°

Tube Style: Double Convertible, Single Convertible, Single Non-convertible

Units per Inner Pack: 4 per Pack, 5 per Pack eClass 14 Classification Group: 34110000

Dimensions and Classifications

Overall Length (Metric): 3.7 mm, 4.3 mm, 4.7 mm

3M™ RelyX™ Ceramic Primer

Details

RelyX Ceramic Primer is designed to improve the bonding of dental adhesives to porcelain and metals

Typical properties

Brands: RelyX™

Catalog Number: 2721 Delivery System: Vial

Factory ISO Certification: No ISO Certification information available

Inner Packs per Case: 10 Packs per Case

Package Contents: 1 – 5ml Vial Total Units per Case: 10 per Case Units per Inner Pack: 1 per Pack

3M™ Durelon™ Carboxylate Luting Cement Maxicap Refill, 56310

3M ID 7000054671 UPC 04035077013031

Details

Simple technique requires no conditioning and no bonding agents for reduced chair time

Thin film thickness Fluoride release

Typical properties

Brands: Durelon™

Catalog Number: 56310 Cure Type: Self-cure

Delivery System: Maxicap™ Capsule

Factory ISO Certification: No ISO Certification information available

Material: Carboxylate

Package Contents: 20 – 3M™ Maxicap™ Capsules

3M™ RelyX™ Try-In Paste Syringe Refill

Details

Customer-preferred shading system with excellent color-match to cured cement Easy to dispense, apply and seat Water-soluble for easy clean-up Compatible color-match to cured $3M^{\text{TM}}$ RelyX $^{\text{TM}}$ Veneer Cement

Complete your final veneer restoration with 3M™ RelyX™ Veneer Cement

Typical properties

Brands: RelyX™

Catalog Number: 7614TRT,7614B0.5T,7614A1T,7614A3T,7614A5T,7614WOT

Cure Type: Light Cure Delivery System: Syringe

Factory ISO Certification: No ISO Certification information available

Package Contents: 1 – 2g A3 Opaque/Yellow Opaque,1 – 2g A5 Opaque/Dark,1

- 2g B0.5/White,1 - 2g Translucent,1 - 2g A1/Light Yellow

Product Type: Try-in Paste

Shade: Multiple, TR, White Opaque

3M™ Adper™ Single Bond 2 Intro Kit, 51201

3M ID 7100156294 UPC 30605861066984

Details

Offers good bond strength and performance, providing your patients with lasting dental restorations

Convenient bottle design offers easier opening and dispensing

No shaking required thanks to stable nano-filler designed not to settle within the bottle

Reduces the risk of post-operative sensitivity, creating a better experience for your patients

Compatible with light-cured composite restoratives

Typical properties

Adhesive Type: Total Etch

Brands: Adper™

Catalog Number: 51201 Cure Type: Light Cure

Factory ISO Certification: No ISO Certification information available

Product Type: Bonding Adhesive Kit

3M™ Unitek™ Stainless Steel Primary Molar Crowns

Details

Easy to customize with parallel walls to allow for clinician fit control Pre-trimmed to optimum length and contour to save placement time Made in the U.S.A. of high quality stainless steel Shallow occlusal anatomy means less tooth reduction Thick occlusal surface to prevent bite-through Crown set is packaged in a unique, stackable tray that provides convenient, easy storage of 3M[™] ESPE[™] prefabricated crowns Wide array of shapes and sizes

Typical properties

Arch: Lower, Upper Brands: Unitek™ Catalog Number:

 $900117,900237,900116,900115,900236,900235,900114,900234,900113,900112\\,900233,900111,900232,900231,900227,900226,900225,900147,900224,90022\\3,900146,900222,900145,900221,900144,900143,900142,900141,900217,9002\\16,900215,900137,900214,900136,900213,900135,900212,900134,900211,900\\133,900132,900131,900127,900247,900126,900246,900125,900124,900245,90\\0244,900123,900122,900243,900242,900121,900241$

Dentition: Primary

Factory ISO Certification: No ISO Certification information available Mesial/Distal Width (Imperial): 0.26 in, 0.272 in, 0.276 in, 0.287 in, 0.291 in, 0.303 in, 0.307 in, 0.319 in, 0.323 in, 0.335 in, 0.339 in, 0.35 in, 0.354 in, 0.366 in, 0.374 in, 0.386 in, 0.394 in, 0.402 in, 0.413 in, 0.417 in, 0.433 in, 0.453 in

Mesial/Distal Width (Metric): 6.6 mm, 6.9 mm, 7 mm, 7.3 mm, 7.4 mm, 7.7 mm, 7.8 mm, 8.1 mm, 8.2 mm, 8.5 mm, 8.6 mm, 8.9 mm, 9 mm, 9.3 mm, 9.4 mm, 9.5 mm, 9.8 mm, 10 mm, 10.2 mm, 10.5 mm, 10.6 mm, 11 mm, 11.5 mm

Number of Crowns: 5 Patient: Pediatric Side: Left, Right Size: 4, 5

Tooth: First Molar, Second Bicuspid, Second Molar

3M™ Strip Crown Form - Pedo

Details

Automatically contours restorative material to match natural dentition

Thin interproximal walls allow for tight contacts

Ideal for use with chemical or light-cured composites

Strips off easily leaving a smooth surface

Provides sufficient strength for easy handling

Crown set is packaged in a unique, stackable tray that provides convenient, easy storage

Typical properties

Arch: Upper Brands: 3M™ Catalog Number:

914004,914014,914003,914024,914013,914002,914034,914001,914023,914012

,914033,914022,914011,914032,914021,914031

Dentition: Primary

Factory ISO Certification: No ISO Certification information available

Mesial/Distal Width (Imperial): 0.169 in, 0.201 in, 0.232 in, 0.236 in, 0.264 in,

0.291 in, 0.319 in

Mesial/Distal Width (Metric): 4.3 mm, 5.1 mm, 5.9 mm, 6 mm, 6.7 mm, 7.4 mm,

8.1 mm

Number of Crowns: 5 Side: Left, Right

Size: 1, 4

Tooth: Central Incisor, Lateral Incisor eClass 14 Classification Group: 34111100

3M™ Lava™ Ultimate CAD/CAM Restorative Intro Kit for CEREC®, 2914SK

3M ID 7000021346 UPC 30605861041134

Details

Fast and easy to mill - no firing necessary Brilliant esthetics with long-lasting polish Less wear to opposing enamel than glass ceramics Resilient, not brittle

Resin nano ceramic material delivers exceptional durability, excellent polish retention

Available in eight shades and two translucencies

Typical properties

Brands: Lava™

CAD/CAM System: CEREC®

Digital Materials: Digital Chairside Restorative Materials

Factory ISO Certification: No ISO Certification information available

Material: Resin Nano Ceramic

Restoration Type: Inlay, Onlay, Veneer

Shade: Multiple Shades

3M™ Paradigm™ Temporization Material Refill

Details

Strong and fracture resistant with high toughness
Nanofiller technology provides strength and esthetics
Available in three of the most popular shades: A1, A2, A3; designed to match
3M™ Filtek™ Supreme Ultra Flowable Restorative
Temporaries that look natural in most conditions
Three of the most popular shades: A1, A2, A3 match 3M™ Filtek™ Supreme
Ultra Flowable Restorative

Typical properties

Brands: Paradigm™

Catalog Number: 46955,46961,46958 Delivery System: Garant™ Dispenser

Factory ISO Certification: No ISO Certification information available

Material: Bis-acrylic Composite

Pre-formed: No Shade: A1, A2, A3

3M[™] Clinpro[™] Sealant & 3M[™] Adper[™] Prompt[™] L-Pop[™] Self-Etch Adhesive Trial Kit, 12626TP, 1.2 mL/Syringe, 1 Syringe/Kit

3M ID 7100156276 UPC 30605861066069

Typical properties

Brands: Clinpro™

Factory ISO Certification: No ISO Certification information available

Inner Packs per Case: 10 Kits per Case

Package Contents: Trial kit includes 1 – 1.2mL sealant syringe, 10 sealant

syringe tips, 20 3M™ Adper™ Prompt™ L Pop™ Applicators

Total Units per Case: 10 per Case Units per Inner Pack: 1 per Case

3M™ Unitek™ Miniature Twin Metal Brackets Kit

Details

A microetched base features strength comparable to mesh Ample under tie-wing area means smoother and more secure ligation Torque-in-base design with compound base contour

Typical properties

Adhesive Precoat: APC™ II, APC™ PLUS, Non-coated

Arch: Upper/Lower Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Hook Included: Yes Hook Type: Distal

Material: Stainless Steel Prescription: MBT™, Roth*

Side: Left/Right

Slot Size (Imperial): 0.018 in, 0.022 in Slot Size (Metric): 0.457 mm, 0.559 mm

Units per Pack: 12, 20

3M™ Clinpro™ Sealant Introductory Kit, 12600 Series

Details

Color change from pink to white to help visualize placement

Convenient direct-delivery syringe for easy application

Low viscosity to flow easily into pits and fissures

Contains and releases fluoride — recharges with use of fluoride-containing products

Bonds to enamel

Lower fill rate for fewer occlusal adjustments than filled sealants

Excellent stain resistance

Can be used with self-etch adhesive or phosphoric acid-etch technique

Typical properties

Brands: Clinpro™

Delivery System: Bottle, Syringe

Factory ISO Certification: No ISO Certification information available

Material: Sealant Units per Case: 10.0

Units per Inner Pack: 2 per Case

3M™ Durelon™ Carboxylate Luting Cement Pack

Details

Simple technique requires no conditioning and no bonding agents for reduced chair time

Thin film thickness Fluoride release

Typical properties

Brands: Durelon™

Catalog Number: 38267,38234,38235,38236,38237

Cure Type: Self-cure

Delivery System: Bottle, Calibrated Dispenser, Hand Mix

Factory ISO Certification: No ISO Certification information available

Material: Carboxylate

Package Contents: 3 – 20g Powder; 40ml Liquid CD (12 Dispensers),3 – 20g Powder, Regular Set; 40ml Liquid; Mixing Pads (70144),3 – 20g Powder, Regular

Set,40ml Triple Liquid,3 - 20g Powder, Fast Set,40ml Triple Liquid CD (12

Dispensers)

3M™ RelyX™ Veneer Cement Syringe Refill

Details

Light-cure formula

Available in customer-preferred delivery and shading system

Excellent color stability

Easy handling and clean-up

Easy to dispense, apply and seat

Ensure precise and convenient color-matching with 3M™ RelyX™ Try-In Paste

Typical properties

Brands: RelyX™ Catalog Number:

7614WO,8714B0.5,8714WO,7614A1,8714A3,7614TR,7614A5,8714TR,7614B0.

5,7614A3,8714A1 Cure Type: Light Cure Delivery System: Syringe

Factory ISO Certification: No ISO Certification information available

Inner Packs per Case: 10 Packs per Case

Package Contents: 1 – 2g White Opaque,1 – 3g A3 Opaque/Yellow Opaque,1 –

3g B0.5/White,1 – 2g Translucent,1 – 3g A5 Opaque/Dark

Shade: A1, A3, B0.5, Multiple, TR, White Opaque

Total Units per Case: 10 per Case Units per Inner Pack: 1 per Pack

3M™ Adper™ Single Bond Plus Adhesive Intro Kit, 51101

3M ID 7100156390 UPC 30605861066892

Details

Offers good bond strength and performance, providing your patients with lasting dental restorations

Convenient bottle design offers easier opening and dispensing

No shaking required thanks to stable nano-filler designed not to settle within the bottle

Reduces the risk of post-operative sensitivity, creating a better experience for your patients

Compatible with light-cured composite restoratives

Typical properties

Adhesive Type: Total Etch

Brands: Adper™

Catalog Number: 51101 Cure Type: Light Cure Delivery System: Vial

Factory ISO Certification: No ISO Certification information available

Package Contents: Intro kit includes 1-6g 3M[™] Adper[™] Single Bond Plus adhesive, 2-3ml phosphoric etchant syringes, mixing wells, 1-brush handle, 1-

bag of brush tips (60) and Technique Guide

Product Type: Dental Adhesive

3M™ Imprint™ Bite Registration Material Refill, 36850

3M ID 7000145314 UPC 04035077001892

Details

Shore-A final hardness of 90
Thixotropic properties
Simple, convenient, automatic mixing with the cartridge dispenser from 3M ESPE

Typical properties

Brands: Imprint™

Catalog Number: 36850

Delivery System: Garant™ Hand Dispenser

Factory ISO Certification: No ISO Certification information available

Impression Type: Bite Registration

Material: VPS

Set Type: Regular Set

3M™ Imprint™ II Monophase VPS Impression Material

Details

High-flow material

Elastic recovery and stability reduce potential for distortion, especially when

using dual-arch trays

Can be disinfected without altering detail reproduction

Void-free and homogenous mix

Hydrophilic

Typical properties

Brands: Imprint™

Catalog Number: 9378,9579,9379,9578,9371P

Delivery System: Garant™ Dispenser, Pentamix™ Mixing Unit Factory ISO Certification: No ISO Certification information available

Impression Type: Final

Inner Packs per Case: 10 Packs per Case

Material: VPS

Package Contents: 4 - 50ml cartridges; 20 - 3M™ Garant™ Mixing Tips - Yellow

Product Color: Green, Orange, Purple, Yellow

Set Time: 2.5 min, 4 min

Set Type: Quick Set, Regular Set Setting Time: 2.5 Minute, 4 Minute Total Units per Case: 10 per Case Units per Inner Pack: 1 per Pack Viscosity: Heavy, Light, Regular

Work Time: 1 min, 40 sec

3M™ Stainless Steel Permanent Molar Crowns

Details

Pre-belled, designed for quick, easy-snap fit

Superior longevity compared to multi-surface amalgam fillings

Accurately duplicates anatomy of molars

Crown set is packaged in a stackable tray for storage of 3M™ prefabricated crowns

Wide array of shapes and sizes

Bicuspid and molar crown selection to cover adult posterior applications

Made in the U.S.A

Accurately duplicates anatomy for better fit and performance

Parallel walls allow for operator fit control

High quality stainless steel

Typical properties

Arch: Lower, Upper

Brands: 3M™

Catalog Number: 6-LR-3,6-LR-2,6-LL-7,6-LL-6,6-LR-7,6-LR-6,6-LR-5,6-LR-4,6-UL-6,6-UL-5,6-UL-4,6-UL-3,6-UL-2,6-UR-7,6-UR-6,6-UR-5,6-UR-4,6-LL-5,6-UR-

3,6-LL-4,6-UR-2,6-LL-3,6-LL-2,6-UL-7

Dentition: Permanent

Factory ISO Certification: No ISO Certification information available

Mesial/Distal Width (Imperial): 0.409 in, 0.421 in, 0.425 in, 0.437 in, 0.441 in,

0.453 in, 0.457 in, 0.469 in, 0.472 in, 0.484 in, 0.488 in, 0.504 in

Mesial/Distal Width (Metric): 10.4 mm, 10.7 mm, 10.8 mm, 11.1 mm, 11.2 mm,

11.5 mm, 11.6 mm, 11.9 mm, 12 mm, 12.3 mm, 12.4 mm, 12.8 mm

Number of Crowns: 5 Patient: Pediatric Side: Left, Right Size: 2, 3, 4, 5, 6, 7 Tooth: First Molar

3M™ Unitek™ Nitinol Heat-Activated Archwire

Details

Can be cooled or chilled resulting in a softer, more pliable wire for easy engagement

Provides light continuous forces

Force activation occurs around body temperature

Available in square sizes making it excellent for early torque control

Good bend retention when bent in hot water

Typical properties

Arch: Lower, Upper

Arch Shape: LA, Ovoid, Square

Archwire Size (Imperial): 0.014 in, 0.014 x 0.025 in, 0.016 in, 0.016 x 0.022 in, 0.016 x 0.025 in, 0.017 x 0.017 x 0.017 x 0.017 in, 0.018 in, 0.018 x 0.018 in,

0.019 x 0.025 in, 0.02 in, 0.021 x 0.021 in, 0.021 x 0.025 in, 0.022 in

Archwire Size (Metric): 0.356 mm, 0.356 x 0.635 mm, 0.406 mm, 0.406 x 0.559 mm, 0.406 x 0.635 mm, 0.432 x 0.635 mm, 0.457 mm, 0.457 x 0.457 mm, 0.483

x 0.635 mm, 0.508 mm, 0.533 x 0.533 mm, 0.533 x 0.635 mm, 0.559 mm

Brands: Unitek™

Cross Section Shape: Hybrid

Factory ISO Certification: No ISO Certification information available

Material: Nitinol

Units per Inner Pack: 10 per Pack

3M™ Sof-Lex™ Pre-Polishing Spiral, 5090, 15ct

3M ID 7100063732 UPC 30605861059153

Details

Used to smooths and removes scratches Unique, flexible shape adapts to all tooth surfaces Multi-use, can be sterilized and reused

Typical properties

Application: Pre-polishing

Brands: Sof-Lex™ Catalog Number: 5090

Factory ISO Certification: No ISO Certification information available Package Contents: Refil includes 15 pre-polishing spirals beige

Product Color: Beige

Product Type: Rubberized Spiral

eClass 14 Classification Group: 34111407

Dimensions and Classifications

Overall Diameter (Imperial): 0.5 in Overall Diameter (Metric): 12.7 mm

3M™ Lava™ Ultimate CAD/CAM Restorative for CEREC®

Details

Fast and easy to mill - no firing necessary
Brilliant esthetics with long-lasting polish

Less wear to opposing enamel than glass ceramics

Resilient, not brittle

Resin nano ceramic material delivers exceptional durability, excellent polish

retention

Available in eight shades and two translucencies

Typical properties

Block Size: 12, 14L Brands: Lava™

CAD/CAM System: CEREC®

Catalog Number: 2914A3-LT,2914A2-LT,2914A3-HT,2914A2-HT

Digital Materials: Digital Chairside Restorative Materials

Factory ISO Certification: No ISO Certification information available

Material: Resin Nano Ceramic

Restoration Type: Inlay, Onlay, Veneer Shade: A1, A2, A3, A3.5, B1, BL, C2, D2

3M™ Lava™ Plus High Translucency Zirconia Shading Liquid

Details

Available in each VITA® classical shade, plus two bleach shades Brilliant esthetics with long-lasting polish

Shading liquid is completely absorbed and provides coloration throughout the restoration

Widest range of dyeing liquids

Lava™ Plus Zirconia shading system enables highly esthetic all-zirconia or traditional layered restorations

Typical properties

Brands: Lava™

CAD/CAM System: Lava

Factory ISO Certification: No ISO Certification information available

Material: Zirconia

Product Type: Liquid Dye

Shade: A1, A2, A3, A3.5, A4, B1, B2, B3, B4, C1, C2, C3, C4, D2, D3, D4, W1,

W3

3M™ Crown Scissor

Details

Designed for trimming the margins of 3M pre-formed crowns

Typical properties

Brands: 3M™

Catalog Number: 801201,801202,801203

Factory ISO Certification: No ISO Certification information available

Product Type: Crown Scissors, Curved Festooning

Tip Shape: Curved, Straight

3M[™] Applicator Handles and Disposable Applicator Brush Tips, 1919, (3 Handles/120 Tips)

3M ID 7100214043 UPC 30605861068803

Details

Handles for 3M™ adhesives application

Typical properties

Brands: 3M™

Factory ISO Certification: No ISO Certification information available

Package Contents: Includes 3 color-coded handles and 120 disposable brush

tips

Product Type: Applicator Handle

3M™ ESPE™ Intraoral Syringe Value Pack

Details

Single-use syringe provides a hygienic procedure with no disinfection or messy cleanup

Smaller, ergonomic design for easier handling

Less impression material waste

Advance prep--can be filled up to 12 hours before application

Compatible with your preferred wash impression material

Compatible with your preferred wash material

Typical properties

Brands: Garant™

Catalog Number: 71508,71506

Factory ISO Certification: No ISO Certification information available Package Contents: 50 - green syringes for VPS precision impression materials,50 - purple syringes for Polyether precision impression materials

Product Type: Syringe

3M™ Penta™ Elastomer Syringe Knurled Hub Cap, 77720

3M ID 7000054661 UPC 04035077012898

Details

Easy to use

Typical properties

Catalog Number: 77720

Factory ISO Certification: No ISO Certification information available

Package Contents: 1 - Knurled Hub Cap

Product Type: Knurled Hub Cap

3M™ Unitek™ Prestige End Cutter, 900-750, 1/Ea

Part Number 900-750 3M Product Number 900-750 3M ID 7000025677

Details

72° angle for very close cut to buccal tube 3M™ Unitek™ Prestige Series

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

3M™ Imprint™ 4 Penta™ Heavy Cartridge for Pentamix™ 3, 77784

3M ID 7000055273 UPC 04035077007924

Details

Durable, metal-reinforced cartridge
Automatic mixing and dispensing with the Pentamix 3 system
Absolutely homogeneous and void-free mixing at the touch of a button
Helps get an accurate impression on the first take--reducing the need for costly adjustments, remakes and retakes

Typical properties

Brands: Imprint™

Catalog Number: 77784

Factory ISO Certification: No ISO Certification information available

Package Contents: 1 - Penta™ Cartridge

Product Type: Cartridge

3M™ Unitek™ Sickle-Style Cement Remover and Bracket Positioner, 900-822, 1/Ea

3M ID 7000004394

Details

Lightweight, well-balanced Serrated tip to assure a positive grip on band edges and seating lugs

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

3M™ AlastiK™ B Modules

Details

Not made with natural rubber latex Injection molded 4 sizes

Typical properties

Brands: AlastiK™

Center to Center Distance (Imperial): 0.145 in, 0.175 in, 0.23 in, 0.325 in Factory ISO Certification: No ISO Certification information available

Force Rating: Heavy, Medium Ligature Size: B-1, B-2, B-3, B-4 Product Color: Gray, Transparent Units per Inner Pack: 100 per Pack

3M™ SmartClip™ SL3 Self-Ligating Brackets Kit

Details

Proprietary self-ligating clips which can save valuable chair time during wire engagement and disengagement

"Intelligent" clips that release the archwire if forces are excessive

Provide selective ligation capability, allowing clinicians to engage the archwire in a single clip

Precisely finished nickel titanium clips and smooth radiused surfaces for enhanced patient comfort

Offered in the MBT[™] System prescription, as well as Ricketts (High Torque) and Roth (Low Torque) for desired treatment philosophy options

3M™ APC™ Flash-Free and 3M™ APC™ II Adhesive coating options for bonding efficiency

Molar brackets with small, low profile tie-wings for treatment flexibility

Notched tie-wings for additional bicuspid identification

Low profile for patient comfort

Torque-in-base for level slot line up

Micro-etched, three point contact base anatomical contouring for better base to tooth fit

Vertical scribe line for accurate bracket placement

Typical properties

Adhesive Precoat: Non-coated

Arch: Upper/Lower Bracket Style: Twin Brands: SmartClip™

Factory ISO Certification: No ISO Certification information available

Hook Included: Yes Hook Type: Distal

Material: Stainless Steel

Prescription: MBT™, Ricketts*, Roth* Slot Size (Imperial): 0.018 in, 0.022 in Slot Size (Metric): 0.457 mm, 0.559 mm

Units per Pack: 20, 24

3M™ Unitek™ Beta III Titanium Archwire

Details

Processed to reduce breakage

Polished surface finish for improved sliding mechanics

Tight corner radii specifications reducing torque loss, especially in final phase of treatment

Typical properties

Arch

Lower, Upper

Arch Shape

Ovoid, Square

Archwire Size (Imperial)

0.016 x 0.022 in, 0.016 x 0.025 in, 0.017 x 0.025 in, 0.0175 x 0.0175 in, 0.019 x

0.025 in, 0.021 x 0.025 in

Archwire Size (Metric)

0.406 x 0.559 mm, 0.406 x 0.635 mm, 0.432 x 0.635 mm, 0.483 x 0.635 mm,

0.533 x 0.635 mm

Brands

Unitek™

Factory ISO Certification

No ISO Certification information available

Material

Beta III Titanium

Units per Inner Pack

10 per Pack

eClass 14 Classification Group

34110604

3M™ RelyX™ Fiber Post 3D Glass Fiber Post

Details

Higher mechanical retention to the core build-up material thanks to coronal 3D macro retentions

Strong bond in the root canal due to its microporous surface without a silanating step

Better post position control thanks to improved radiopacity that is among the highest of available glass fiber posts

Low risk of root fracture due to dentin-like elasticity

More esthetic option than metal post

Excellent light translucency for efficient light cure deep into the root canal

Typical properties

Apical Width (Metric)

0.6 mm, 0.7 mm, 0.8 mm, 0.9 mm

Brands

RelyX™

Catalog Number

56951,56950,56948,56949

Coronal Width (Metric)

1.1 mm, 1.3 mm, 1.6 mm, 1.9 mm

Factory ISO Certification

No ISO Certification information available

Material

Fiber Post

Package Contents

10 Posts - Size 1 (1.3mm dia.) Yellow, 10 Posts - Size 0 (1.1mm dia.) White, 10

Posts - Size 2 (1.6mm dia.) Red,10 Posts - Size 3 (1.9mm dia.) Blue

Product Color

Blue, Red, White, Yellow

Product Type

Fiber Post 3D

Size

0, 1, 2, 3

3M[™] Ketac[™] Nano Light-Curing Glass Ionomer Quick Mix Restorative Trial Kit, 3305TK

3M ID 7000128833 UPC 30605861012073

Details

Quick Mix Capsule for intraoral placement

The first pasts/pasts, resin-modified class innemer (PMGI) re

The first paste/paste, resin-modified glass ionomer (RMGI) restorative based on bonded nanofiller technology

Esthetic glass ionomer restorative

High fluoride release

Improved wear resistance

Light cure on demand

Typical properties

Brands

Ketac™

Catalog Number

3305TK

Cure Type

Light Cure

Delivery System

Capsule

Factory ISO Certification

No ISO Certification information available

Package Contents

Trial kit includes 10 - A2 Quick Mix Capsules; 10 - A3 Quick Mix Capsules;

Primer (3ml); Technique Guide

Total Units per Case

10 per Case

eClass 14 Classification Group

34110000

3M™ Ketac™ Cem Glass Ionomer Luting Cement Refill, 56000 Series

Details

Glass ionomer luting cement offers chemical bond to enamel and dentin Releases fluoride High biocompatibility Extremely low film thickness optimizes fit and marginal integrity Radiopaque

Typical properties

Brands
Ketac™
Catalog Number
56021,56060,56061
Cure Type
Self-cure
Delivery System
Aplicap™ Capsule, Maxicap™ Capsule
Factory ISO Certification
No ISO Certification information available
Package Contents
50 – 3M™ Aplicap™ Capsules,50 – 3M™ Maxicap™ Capsules

3M™ Victory Series™ Second Molar Bands

Details

Permanent laser ID marking Straight interproximal walls for reduced separation Available with micro-etched inner surface for increased retention Increased occlusal "edge" strength for reduced deformation

Typical properties

Arch

Lower, Upper

Band Size

128, 128+, 129, 129+, 130, 130+, 131, 131+, 132, 132+, 133, 133+, 134, 134+,

135, 135+, 136, 136+, 137, 137+, 138, 138+, 139, 139+, 140, 140+, 141, 28,

28+, 29, 29+, 30, 30+, 31, 31+, 32, 32+, 33, 33+, 34, 34+, 35, 35+, 36, 36+, 37,

37+, 38, 38+, 39, 39+, 40, 40+, 41

Brands

Victory Series™

Factory ISO Certification

No ISO Certification information available

Material

Stainless Steel

Product Type

Band

Side

Left, Right

Tooth

Second Molar

eClass 14 Classification Group

34110603

3M[™] Photac[™] Fil Quick Aplicap[™] Light-Cured Glass Ionomer Restorative, 61000 Series

Details

High-fluoride release that can recharge after appluing a topical fluoride source Self-bonding and light curable

Excellent mechanical properties, plus it's radiopaque and biocompatible Fewer steps - no glaze required

Convenience, simplicity and accuracy of the 3M[™] Aplicap[™] delivery system Conservative cavity preparation, fluoride releas eand fast, easy technique make it ideal for geriatric and pediatric patients

Dual cure provides light-cure control and the assurance of a complete cure Dual cure provides light-cure control and the assurance of a complete cure

Typical properties

Brands

Photac™

Catalog Number

61060,61070,61080,61090,61010,61030,61050

Cure Type

Light Cure

Delivery System

Capsule

Factory ISO Certification

No ISO Certification information available

Package Contents

Refill includes 50 capsules in assorted shades, Refill includes 50 capsules in shade A1, Refill includes 50 capsules in shade B2, Refill includes 50 capsules in shade A2, Refill includes 50 capsules in shade A35, Refill includes 20 capsules in shade C4, Refill includes 50 capsules in shade A3, Refill includes 20 capsules in shade DBO

Shade

A1, A2, A3, A3.5, B2, C4, DBO, Multiple Shades eClass 14 Classification Group 34110000

3M™ Filtek™ Flowable Dispensing Tips

Details

Single use for easy dispensing

Typical properties

Brands

Filtek™

Catalog Number

3700T27,3700T

Environmental claim being made

Yes

Environmental claim relates to

Product

Factory ISO Certification

No ISO Certification information available

Features

Access to Hard to Reach Cavity Preparations, Easy-to-bend Cannula, Finer

Gauge, Improved Access, More Precise Material Placement

Package Contents

Refill includes 20 blue dispensing tips (20-gauge), Refill includes 20 dispensing

tips

Product Type

Tip

eClass 14 Classification Group

34110108

3M)™ AlastiK™ Plastic Box With Dividers, 406-890, 1/Ea

3M ID 7100154417

Details

Not made with natural rubber latex Each divided section can hold up to 50 sticks Latch close

Typical properties

Brands
AlastiK™
Factory ISO Certification
No ISO Certification information available
Units per Inner Pack
1 per Pack
eClass 14 Classification Group
34110603

3M[™] Unitek[™] Ligature Wire, 390-160, Preformed Coated, Natural Colored, .012 in., 100/Pk

3M ID 7100154542

Details

Preformed ligature slips easily over all brackets Coated preformed ligature wire; ideal for esthetic brackets such as 3M™ Clarity Brand Brackets

Dead soft permachrome wire in color coded bundles of 200 Coated ligature wire reduces friction between archwire and bracket

Typical properties

Brands
Unitek™
Factory ISO Certification
No ISO Certification information available
Units per Inner Pack
100 per Pack
eClass 14 Classification Group
34110603

Dimensions and Classifications

Overall Diameter (Imperial) 0.012 in

3M™ Unitek™ Prestige Series Omega Loop Forming Pliers Replacement Tip, 900-708, 1/Ea

3M ID 7000025663

Details

Fine tip for precise bends 3M™ Unitek™ Prestige Series

Typical properties

Brands
Unitek™
Factory ISO Certification
No ISO Certification information available
eClass 14 Classification Group
21043800

3M™ Unitek™ How Utility Offset Pliers, 801-701, 1/Ea

3M ID 7100069059

Details

Make uniform Omega loops in square and rectangular wires One-piece tip construction with serrated concave beak

Typical properties

Brands
Unitek™
Factory ISO Certification
No ISO Certification information available
eClass 14 Classification Group
21043800

3M™ Unitek™ Lingual Distal End Cutter, 803-552, 1/Ea

3M ID 7100267824

Details

Light duty cutter specifically for soft wires and elastics Sharp pointed tips can reach into tight and difficult to reach areas Tips are made of Tungsten Carbide for durability

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

3M™ Unitek™ General Purpose Molar Bands

Details

Universal anatomy fits rights and lefts with minor burnishing Permanent laser ID Available in narrow or wide

Typical properties

Arch: Lower, Upper

Band Size: 128, 128+, 129, 129+, 130, 130+, 131, 131+, 132, 132+, 133, 133+, 134, 134+, 135, 135+, 136, 136+, 137, 137+, 138, 138+, 139, 139+, 140, 140+, 141, 141+, 142, 25, 26, 27, 28, 28+, 29, 29+, 30, 30+, 31, 31+, 32, 32+, 33, 33+, 34, 34+, 35, 35+, 36, 36+, 37, 37+, 38, 38+, 39, 39+, 40, 40+, 41, 41+, 42

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Material: Stainless Steel Product Type: Band

Tooth: Molar

3M[™] Forsus[™] Fatigue Resistant Device, 885-110, L-Pins Spring Module, 10/Pk

3M ID 7000004453

Details

Universal L-Pin

Typical properties

Brands: Forsus™

Factory ISO Certification: No ISO Certification information available

Module Style: L-Pin

Product Type: Fatigue Resistant Device

Units per Inner Pack: 10 per Pack

3M™ Ortholux™ Luminous Curing Light Glasses, S0360X, 1/Ea

3M ID 78000921904

Details

Orange
Can be disinfected

Typical properties

Brands: Ortholux™

Factory ISO Certification: No ISO Certification information available

3M™ CoJet™ Sand, 68411, 3 Bottles

3M ID 7000054650 UPC 04035077012829

Details

Blasting with the CoJet system embeds a ceramic layer on the surface of metal, porcelain or previously cured composite, producing dramatically enhanced bond strengths

Embedded ceramic layer enables silane coupling with composite for a true chemical bond with alloys and porcelain

Chairside, intraoral technique allows for highly esthetic repair of porcelain or composite restorations and provides an economical alternative to replacing the restoration. Plus, repairs can be made in one session

Typical properties

Brands: CoJet™

Catalog Number: 68411 Delivery System: Bottle

Factory ISO Certification: No ISO Certification information available

Product Type: Adhesive System Units per Inner Pack: 3 /BT

3M™ Astringent Retraction Paste

Details

Offers fast and precise gingival retraction

Extra-fine tip fits right into the sulcus for precise placement

Long-lasting hemostasis

Gentle on tissue with less risk of tissue trauma or hemorrhaging than cords

Provides a more comfortable experience for patients

Hygienic unit-dose capsule for up to 3 teeth

Can be used with or without cords

Get precision at every step and create more detailed impressions for a wide range of cases with 3M™ Impregum™ Impression Material or 3M™ Imprint™ 4

VPS Impression Material

Typical properties

Brands: 3M™

Catalog Number: 56944,56945

Delivery System: Capsule

Factory ISO Certification: No ISO Certification information available

Impression Type: Gingival Retraction

Package Contents: 25 - Capsules, 0.3 g each eClass 14 Classification Group: 34110100

3M™ Ortholux™ Luminous Curing Light Guide, 704-456, 1/Ea

3M ID 7100018470 UPC 00652221149374

Details

8mm diameter tip Black coated High intensity 1600 mW/cm2 blue LED Built-in light meter Rechargeable Lithium ion battery

Typical properties

Brands: Ortholux™

Catalog Number: 704-456

Factory ISO Certification: No ISO Certification information available

Product Type: Light Guide

3M™ RelyX™ Ultimate Clicker™ Adhesive Resin Cement Kit

3M ID B10234982

Details

High bond strength
Easy-to-use, two-component system
Versatile for total-etch, selective etch or self-etch techniques
Excellent color stability with natural fluorescence
Offers wear resistance
If you like 3M™ RelyX™ Ultimate Cement, you'll love 3M™ RelyX™ Universal Resin Cement

Typical properties

Brands: RelyX™ Cure Type: Dual Cure

Delivery System: Clicker™ Dispenser

Material: Adhesive Resin

3M™ Impregum™ Super Quick Polyether Impression Material

Details

Superfast polyether with 45-second working time and quick 2-minute setting time Unique moisture tolerance provides excellent detail reproduction for void-free impressions, even in wet environments

Very high level of accuracy for precise-fitting restorations

Improved patient comfort thanks to fresh, minty taste and shorter setting time Scannable for more efficient workflow

Available in a range of viscosities and delivery systems for monophase and 1step impression techniques

Control moisture at every step and gently retract gingival tissue for more accurate impressions with 3M[™] Astringent Retraction Paste

Typical properties

Brands: Impregum™

Catalog Number: 69379,69417,69406,69418,69385

Delivery System: Garant™ Dispenser, Pentamix™ Mixing Unit Factory ISO Certification: No ISO Certification information available

Impression Type: Final Material: Polyether

Package Contents: 4 cartridges – 50 ml each, 5 3M™ Intra-oral Syringes Purple,

5 3M™ Garant™ Mixing Tips purple,4 cartridges – 50 ml each, 10 3M™

Garant™ Mixing Tips purple,2 base pastes - 300 ml each, 2 catalysts - 60 ml

each

Product Color: Dark Purple, Green, Purple

Set Time: 2 min

Set Type: Super Quick Set Setting Time: 2 Minute

Viscosity: Heavy, Light, Medium, Multiple Viscosities

Work Time: 45 seconds

3M™ Transbond™ LR Light Cure Adhesive

Details

Designed specifically for ortho bonding of lingual retainers Excellent handling properties for easy application and control Efficient, accurate and low-waste delivery system available in capsules

Typical properties

Brands: Transbond™

Catalog Number: 712-033,712-038

Cure Type: Light Cure Delivery System: Capsule

Factory ISO Certification: No ISO Certification information available

Package Contents: 25 - 0.2g capsules

Product Type: LR Capsule Kit Units per Inner Pack: 25 per Pack

3M™ Garant™ Intraoral Tips, Intraoral Dispensing

Details

Allows easy intraoral application of bite registration material using the cartridge dispenser

Typical properties

Brands: Garant™

Catalog Number: 71460,71461N

Factory ISO Certification: No ISO Certification information available

Package Contents: 50 - Intraoral Tips, white,50 - Intraoral Tips, Bite, Clear

Product Type: Tip

3M™ Victory Series™ V-Slot Brackets Kit

Details

Full contoured mesh base provides precision placement and consistent bonds The low profile torque-in-base design achieves level slot line-up and easy bracket placement

Chamfered .018 x .018 V-slot allows easy insertion and removal of various auxiliaries to aid tooth alignment

Typical properties

Adhesive Precoat: APC™ II, Non-coated

Arch: Upper/Lower Brands: Victory Series™

Factory ISO Certification: No ISO Certification information available

Hook Included: Yes

Hook Type: Cuspid, Distal

Prescription: Roth* Side: Left/Right

Slot Size (Imperial): 0.018 in, 0.022 in Slot Size (Metric): 0.457 mm, 0.559 mm

Tooth: Cuspid, Cuspid/Bicuspid

Units per Pack: 5, 20

3M[™] Forsus[™] Fatigue Resistant Device, 885-200, Reactivation Spacer, 10/Pk

3M ID 7100069253

Details

Closed ring slides easy on to the push rod and rests at the bend without any special installation 10/pack

Typical properties

Brands: Forsus™

Factory ISO Certification: No ISO Certification information available

Product Type: Fatigue Resistant Device Units per Inner Pack: 10 per Pack

3M™ Filtek™ Restoratives Dispenser, 5850RD, 1 Dispenser

Part Number 5850RD 3M Product Number 5850RD 3M ID 7100196622

Details

Rotatable capsule holder is designed for easy load in, no matter the angle Features a "snap-in" function that secures capsules during dispensing and makes it easier to load/unload capsules

Easy-pull trigger and ergonomic design

Compatible with warming devices like the Bioclear® HeatSync® Composite Capsule warming device

Can be cleaned manually or with a washer-disinfector and can be high level disinfected or sterilized (tested for 200 cycles)

Typical properties

Brands: Filtek™

Environmental claim being made: Yes Environmental claim relates to: Product

Factory ISO Certification: No ISO Certification information available

Package Contents: 1 - Dispenser, Instructions for Use

Product Type: Dispenser

3M™ Penta™ Universal Cartridge for Pentamix™ 3, 71545

3M ID 7100215172

Details

Supports consistent, homogeneous and void-free mixing at the touch of a button Compatible with 3M™ Pentamix™ 3 Mixing Unit

Reusable, metal-reinforced cartridge

Helps get an accurate impression on the first take, so you can avoid needing costly adjustments, remakes and retakes

Typical properties

Brands: Penta™

Catalog Number: 71545

Factory ISO Certification: No ISO Certification information available

Package Contents: 1 - Penta™ Cartridge

Product Type: Cartridge

3M™ Victory Series™ First Molar Bands Kit

Details

Permanent laser ID marking

Straight interproximal walls for reduced separation

Reduced gingival/occlusal dimensions for less occlusal and gingival interferences

Available with micro-etched inner surface for increased retention

Curved festooned occlusal interproximal

Increased occlusal "edge" strength for reduced deformation

Typical properties

Arch: Lower, Upper Brands: Victory Series™

Factory ISO Certification: No ISO Certification information available

Material: Stainless Steel

Product Type: Band, Introduction Kit, Master Kit, Professional Kit, Trial Kit

Tooth: First Molar

3M™ Concise™ Orthodontic Bonding Paste A, 18 g, 196-002, 1/Ea

3M ID 7010386569 UPC 00652221047830

Details

Jar 18 g

Typical properties

Adhesive Type: Direct Brands: Concise™

Catalog Number: 196-002 Cure Type: Chemical Cure

Factory ISO Certification: No ISO Certification information available

Product Type: Orthodontic Adhesive Units per Inner Pack: 1 per Pack

3M™ Unitek™ Twist Wire

Details

Three (3) strand twist More continuous light force over greater working range than solid round wire of same diameter

Typical properties

Arch: Lower, Upper Arch Shape: Square

Archwire Size (Imperial): 0.0175 x 0.0175 in

Brands: Unitek[™]

Factory ISO Certification: No ISO Certification information available

Material: Stainless Steel

Units per Inner Pack: 10 per Pack

3M™ Clinpro™ Tooth Créme 0.21% Sodium Fluoride Anti-Cavity Toothpaste

Details

Help your patients improve their everyday oral health with a fluoride toothpaste that strengthens teeth and repairs early lesions

Formulated with 0.21% NaF (950 ppm fluoride ion) and an innovative calcium and phosphate ingredient

Prevents, repairs and reverses early and white spot lesions

Proven remineralization throughout the lesion, not just on the surface

Helps reverse root caries

Creates acid-resistant mineral in dentin tubules

Recharges fluoride-releasing sealants and restorations

Features a non-abrasive formula that provides gentle, effective cleaning and whitening

Typical properties

Active Ingredient: Sodium Fluoride

Brands: Člinpro™ Delivery System: Tube

Factory ISO Certification: No ISO Certification information available

Features: Anti-cavity

Flavor: Spearmint, Vanilla Mint Fluoride Content: 0.21 % Material: Sodium Fluoride

Net Weight (Imperial): 0.847 oz, 3.985 oz, 4 oz, 4.691 oz Net Weight (Metric): 24 g, 113 g, 113.398 g, 133 g

Product Form: Paste Product Type: Dentifrice

Units per Inner Pack: 1 per Pack, 2 per Pack, 24 per Pack

3M™ Protemp™ Crown Temporization Material Kit

Details

Nine preformed sizes with malleable walls for a custom fit Natural-looking esthetics Low intra-oral heat during fitting process Excellent strength Light curable Low intra-oral heat during fitting process

Typical properties

Arch: Upper/Lower Brands: Protemp™

Catalog Number: 50600TK,50600

Dentition: Permanent

Factory ISO Certification: No ISO Certification information available

Material: Composite

Number of Crowns: 16, 42

Pre-formed: Yes Side: Left/Right

Tooth: Molar & Bicuspid

3M™ ESPE™ Ketac™ Molar Easymix Glass Ionomer Filling Material Refill, 56625

3M ID 7100041170

Details

High flexural strength – reduces the risk of restoration fracture
Adhesion to enamel and dentin – ideal for Minimal Intervention Dentistry (MID)
and Atraumatic Restorative Treatment (ART)
Low acid erosion – maintains excellent marginal integrity of restorations
More durable and longer lasting

Typical properties

Brands: Ketac™

Catalog Number: 56625 Cure Type: Self-cure Delivery System: Bottle

Factory ISO Certification: No ISO Certification information available

3M™ Adper™ Scotchbond™ 1 XT Adhesive

3M ID B00007542

Details

Offers good bond strength and performance, providing your patients with lasting dental restorations

Convenient bottle design offers easier opening and dispensing

No shaking required thanks to stable nano-filler designed not to settle within the bottle

Reduces the risk of post-operative sensitivity, creating a better experience for your patients

Compatible with light-cured composite restoratives

Typical properties

Brands: Adper™

3M™ Scotchbond™ Etchant, 7423

3M ID 7000030404 UPC 30605861032699

Details

Blue-colored for precise application Flexible dispensing tip can be individually adjusted before application

Typical properties

Brands: Scotchbond™ Delivery System: Vial

Factory ISO Certification: No ISO Certification information available

Package Contents: Refill includes 9ml bottle eClass 14 Classification Group: 34110204

3M™ Penta™ Cartridge for Pentamix™ Lite, 77944

3M ID 7100215302 UPC 04054596734142

Details

For consistent, homogeneous and void-free mixing at the touch of a button Compatible with 3M™ Pentamix™ Lite Mixing Unit Reusable, metal-reinforced cartridge Helps get an accurate impression on the first take, reducing the need for costly adjustments, remakes and retakes

Typical properties

Brands: Penta™

Catalog Number: 77944

Factory ISO Certification: No ISO Certification information available

Package Contents: 1 - Penta™ Cartridge

Product Type: Cartridge

3M™ Polycarbonate Crowns, C-180, Anterior, Kit

3M ID 7000030120 UPC 30605861011151

Details

Polycarbonate material contours and crimps similar to metal crowns to yield easy adaptation to margins

Universal shade No. 62 can be adjusted using cements and liners for better esthetics

Crown set is packaged in a unique, stackable tray that provides convenient, easy storage of 3M™ ESPE™ prefabricated crowns

Typical properties

Arch: Upper/Lower Brands: 3M™

Catalog Number: C-180 Dentition: Permanent

Factory ISO Certification: No ISO Certification information available

Material: Polycarbonate Resin

Number of Crowns: 120

Side: Left/Right Size: Assorted Tooth: Assorted

3M™ Paradigm™ MZ100 Blocks for CEREC®

Details

Created from Z100 restorative material from 3M
High flexural strength and fracture toughness
Enamel-like wear characteristics are superior to that of ceramic blocks
Can be easily characterized with light-cured composite stains and glazes
An ideal material for CAD/CAM where strength and occlusal wear properties are needed
Easy to repair intraorally

Typical properties

Block Size: 14

Brands: Paradigm™

CAD/CAM System: CEREC® Catalog Number: 2714E-S

Digital Materials: Digital Chairside Restorative Materials

Factory ISO Certification: No ISO Certification information available

Material: Composite

Restoration Type: Crown, Inlay, Onlay, Veneer

Shade: A1, A2, A3, A3.5, B3, Enamel

Lava™ Frame Shading Liquid

3M ID B00008233

Details

Patented substructure system available in eight shades: Unshaded, FS1, FS2, FS3, FS4, FS5, FS6 and FS7

Typical properties

Brands: Lava™

3M™ Unitek™ Permachrome Resilient Archwire

Details

Better formability than Hi-T™ II Mid-range tensile strength

Typical properties

Arch: Lower, Upper

Arch Shape: LA, Ovoid, Square, Standard, Tapered

Archwire Size (Imperial): 0.012 in, 0.014 in, 0.016 in, 0.016 x 0.016 in, 0.016 x 0.022 in, 0.016 x 0.025 in, 0.017 x 0.022 in, 0.017 x 0.025 in, 0.017 x 0.025 in, 0.017 x 0.025 in, 0.018 x 0.025 in, 0.018 x 0.025 in, 0.018 x 0.025 in, 0.025 in

Archwire Size (Metric): 0.305 mm, 0.356 mm, 0.406 mm, 0.406 x 0.406 mm, 0.406 x 0.559 mm, 0.406 x 0.635 mm, 0.432 x 0.559 mm, 0.432 x 0.635 mm, 0.457 mm, 0.457 x 0.559 mm, 0.457 x 0.635 mm, 0.483 x 0.635 mm, 0.508 mm, 0.532 x 0.635 mm, 0.535 mm, 0.555 mm,

0.533 x 0.635 mm Brands: Unitek™

Cross Section Shape: Hybrid

Factory ISO Certification: No ISO Certification information available

Material: Stainless Steel

Units per Inner Pack: 10 per Pack, 50 per Pack eClass 14 Classification Group: 34110604

3M™ Sof-Lex™ Diamond Polishing Spiral, 5091

3M ID 7100063733 UPC 30605861059177

Details

Diamond paste-like gloss in the convenience of a rubberized system Unique, flexible shape adapts to all tooth surfaces Multi-use, can be sterilized and reused

Typical properties

Application: Polishing Brands: Sof-Lex™
Catalog Number: 5091

Factory ISO Certification: No ISO Certification information available Package Contents: Refil includes 15 diamond polishing spirals pink

Product Color: Pink

Product Type: Rubberized Spiral

eClass 14 Classification Group: 34111407

Dimensions and Classifications

Overall Diameter (Imperial): 0.5 in Overall Diameter (Metric): 12.7 mm

956-208 CNTRD STD ME U7

Typical properties

Arch: Lower, Upper Brands: Victory Series™

Factory ISO Certification: No ISO Certification information available Product Type: Band, Introduction Kit, Professional Kit, Trial Kit

Tooth: Second Molar

3M™ Unitek™ Bracket Positioning Gauge Five, 900-840, 6.0-6.5 mm, 1/Ea

Part Number 900-841 3M ID 7000004396

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

3M™ AlastiK™ Continuous C Module Chain Reel

Details

Not made with natural rubber latex 3 ft. (91 cm) per reel Die-cut chain Medium force

Typical properties

Brands: AlastiK™

Center to Center Distance (Imperial): 0.11 in, 0.14 in

Factory ISO Certification: No ISO Certification information available

Force Rating: Medium Ligature Size: C-1, CK

Product Color: Gray, Transparent

Units per Inner Pack: 4 per Pack, 40 per Pack eClass 14 Classification Group: 34110603

3M™ Unitek™ Auxiliary Ligature Tie Hook 390-145, 0.014, 100/Pk

3M ID 7100071110

Details

Slips under tie-wing, secures like ligature Tie-in mesially or distally with archwire in place

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Units per Inner Pack: 100 per Pack

eClass 14 Classification Group: 34110603

Dimensions and Classifications

Overall Diameter (Imperial): 0.014 in

3M™ Unitek™ Posterior Band Removing Pliers, 801-713, 1/Ea

3M ID 7100069081

Details

Long, narrow chisel-tip beak facilitates posterior band removal Replaceable occlusal pad

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

3M™ Unitek™ Lingual Retainer Kit, 709-060, 1 Kit

3M ID 70202089143

Details

25 each Twist wire 50 each AlastiK™ EZ3 Holders Complete Instructions

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

3M™ Unitek™ Lingual Anterior Debonding Pliers, 803-532, 1/Ea

3M ID 7100266451

Details

Gives a flush cut to the wire at the distal end of the terminal tube Tips are made of Tungsten Carbide for durability Safety hold feature secures the loose wire after the archwire is cut

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

3M™ Unitek™ General Purpose Molar Bands Kit

Details

Universal anatomy fits rights and lefts with minor burnishing Permanent laser ID Available in narrow or wide

Typical properties

Arch: Upper/Lower Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Material: Stainless Steel

Product Type: Introduction Kit, Professional Kit

Tooth: Molar

3M[™] Forsus[™] Fatigue Resistant Device, 885-100, L-Pin Spring Modules, 5/Pk

3M ID 7000004452

Details

EZ2 Module for simple, secure installation L-pin Module for flexible, more mobile installation L-pin Module for flexible, more mobile installation

Typical properties

Brands: Forsus™

Factory ISO Certification: No ISO Certification information available

Module Style: L-Pin

Product Type: Fatigue Resistant Device

Units per Inner Pack: 5 per Pack

3M™ Forsus™ Fatigue Resistant Device, Split Crimps, 885-203, 50/Pk

3M ID 7000004454

Details

Used to reactivate the spring module as treatment progresses Universal

Typical properties

Brands: Forsus™

Factory ISO Certification: No ISO Certification information available

3M[™] Unitek[™] Multi-Cure Glass Ionomer Orthodontic Band Cement Dispensing Spoon, 712-053, 1/Ea

3M ID 7100268714

Details

Allows measuring the powder component of the 3M Multi-Cure Glass Ionomer Band Cement

Double-ended spoon for dispensing different amounts of powder

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

3M[™] Unitek[™] Multi-Cure Glass Ionomer Band Cement Powder, 712-051, 35 g, 1/Btl

Part Number 712-051 3M ID 7000004385

Details

Fills the voids between the anatomy of the tooth surface and the band, minimizing wash out during the lifetime of the band Provides a sustained fluoride release during treatment

Dissoing charmy flavor

Pleasing cherry flavor

Blue color increases visibility for easy cleanup

Doctor and patient friendly, resulting in shorter banding appointments and less loose bands during treatment

Typical properties

Brands: Unitek™

Catalog Number: 712-051

Factory ISO Certification: No ISO Certification information available

Units per Inner Pack: 1 /BT

3M™ Series 5 Facebows

Details

Unless otherwise specified, Series 5 facebows are shipped with adjustable loops bent occlusally

For loops bent gingivally, indicate "special" when ordering and specify "adjustable loops bent gingivally"

Typical properties

Brands: 3M™

Factory ISO Certification: No ISO Certification information available

Hook Type: Cuspid

Inner Bow Length (Metric): 83 mm, 90 mm, 97 mm, 104 mm, 111 mm

Inner Bow Size: 1, 2, 3, 4, 5 Outer Bow Angulation: 0°, 15°

Outer Bow Size: Long, Short, Super Short

Dimensions and Classifications

Overall Diameter (Imperial): 0.045 in

3M™ Concise™ Orthodontic Bonding Paste B, 18 g, 196-003, 1/Ea

3M ID 7010315430 UPC 00652221047847

Details

Jar 18 g

Typical properties

Adhesive Type: Direct Brands: Concise™

Catalog Number: 196-003 Cure Type: Chemical Cure

Factory ISO Certification: No ISO Certification information available

Product Type: Orthodontic Adhesive Units per Inner Pack: 1 per Pack

3M™ Unitek™ Braided Archwire

Details

Three dimensional control Wide working range

Typical properties

Arch: Lower, Upper

Arch Shape: Ovoid, Square

Archwire Size (Imperial): 0.016 x 0.016 in, 0.016 x 0.022 in, 0.017 x 0.025 in,

0.018 x 0.025 in, 0.019 x 0.025 in, 0.021 x 0.025 in

Archwire Size (Metric): 0.406 x 0.406 mm, 0.406 x 0.559 mm, 0.432 x 0.635 mm,

0.457 x 0.635 mm, 0.483 x 0.635 mm, 0.533 x 0.635 mm

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Material: Stainless Steel

Units per Inner Pack: 10 per Pack

3M™ Digital Bonding Tray

Details

Designed for use with the most efficient bonding system in orthodontics, 3M[™] APC[™] Flash-Free Precoated Brackets, which eliminates the need to add adhesive and remove flash

Compatible with 3M[™] Victory Series[™] Low Profile Brackets, 3M[™] Clarity[™] Advanced Ceramic Brackets and 3M[™] SmartClip[™] SL3 Self-Ligating Brackets Reduces bonding procedure time by 30% compared to direct bonding Offers 100% satisfaction rate with bracket positioning after tray removal Two-layer design offers soft gel inlay to conform to patient's anatomy and facilitate easy tray removal. Hard outer shell adds stability during placement Easy to insert, seat and remove

Blue midline marker assists with tray alignment when seating

Typical properties

Arch: Upper/Lower Brands: 3M™

Catalog Number: 355-204,355-205,355-202,355-203,355-200,355-201 Factory ISO Certification: No ISO Certification information available

3M™ RelyX™ Universal Micro Mixing Tips

Details

Micro Mixing Tips for 3M™ RelyX™ Universal automix syringe Excellent mixing quality 80% less cement waste compared to standard automix systems Self-sealing syringe allows hygienic storage without the used tip Thin, long and easy to attach elongation tip for easy application into root canal All applications with only two mixing tip components

Typical properties

Brands: RelyX™

Catalog Number: 56975,56976

Factory ISO Certification: No ISO Certification information available Package Contents: 15x 3M™ RelyX™ Universal Micro Mixing Tips with

Elongation Tips,30x 3M™ RelyX™ Universal Micro Mixing Tips

Product Type: Tip

3M™ Unitek™ Bicuspid Bands Kit

Details

Universal anatomy fits rights and lefts Regular temper Permanent laser ID

Typical properties

Arch: Lower, Upper Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Material: Stainless Steel

Product Type: Band, Bicuspid Bands Kit, Introduction Kit, Master Kit,

Professional Kit, Trial Kit

Side: Left/Right Tooth: Bicuspid

3M™ Ketac™ Cem Maxicap™ Permanent Glass Ionomer Luting Cement Intro Kit, 56015

3M ID B5005238002

Details

Glass ionomer luting cement offers chemical bond to enamel and dentin Releases fluoride
High biocompatibility
Extremely low film thickness optimizes fit and marginal integrity
Radiopaque

Typical properties

Brands: Ketac™

Catalog Number: 56015

Delivery System: Maxicap™ Capsule

3M™ Adper™ Scotchbond™ Multi-Purpose Primer

Details

Improve patient outcomes with a total-etch system that offers excellent bond strength to a variety of dental surfaces

Creates a complete, versatile dental bonding system when used with 3M[™] Adper[™] Scotchbond[™] Multi-Purpose Adhesive

Obtain strong and reliable bonds both to wet and dry surfaces, even when complete isolation from moisture is difficult

Reduces microleakage that can contribute to post-operative sensitivity, helping to ensure a positive patient experience

Typical properties

Adhesive Type: Total Etch

Brands: Adper™

Catalog Number: 3008,7542

Container Size (Imperial): 0.271 fl oz

Container Size (Metric): 8 mL

Cure Type: Light Cure

Delivery System: Bottle, Vial

Factory ISO Certification: No ISO Certification information available

Package Contents: Refill includes 8 ml bottle

Product Type: Dental Primer

3M™ Victory Series™ 2nd Bicuspid Tubes

Details

No tie-wings to come into occlusion Ligature wire channel Chamfered tube opening to assist in wire insertion Micro-Etched base with offset option Metal Injection Molded, including integral hook

Typical properties

Adhesive Precoat: APC™ II, APC™ PLUS, Non-coated

Angulation: +2°, 0° Arch: Lower, Upper Base Type: Bondable Brands: Victory Series™

Factory ISO Certification: No ISO Certification information available

Hook Included: Yes Hook Type: Distal

Material: Stainless Steel Prescription: MBT™

Side: Left, Left/Right, Right

Slot Size (Imperial): 0.018 in, 0.022 in Slot Size (Metric): 0.457 mm, 0.559 mm

Tooth: Bicuspid, Second Bicuspid

Torque: -17°, -7°

Units per Inner Pack: 5 per Pack

3M™ AlastiK™ Color Palette, 406-825, 1/Ea

3M ID 7000004557

Details

Available replacement color dots

Typical properties

Brands: AlastiK™

Factory ISO Certification: No ISO Certification information available

Units per Inner Pack: 1 per Pack

3M™ Unitek™ Single-Ended Ligature Director, 803-003, 1/Ea

3M ID 7100154495

Details

Notched tip at one end for use on anterior brackets Notched tip at other end for easier access to posterior brackets

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Units per Inner Pack: 1 per Pack

3M™ Unitek™ Prestige Posterior Band Removing Pliers

Details

Available in long and short beak Replacement pads available 3M™ Unitek™ Prestige Series

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Pliers Type: Band Remover

Product Type: Band

3M™ Deluxe AlastiK™ Force Modules Needle Holder, 803-535, 1/Ea

3M ID 7100272870

Details

Enables .75 mm 1st and 2nd order bends intraorally without having to remove the archwire

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

3M[™] Unitek[™] Lingual Light Wire Flush Cut Distal End Cutter, 803-550, 1/Ea

3M ID 7100266452

Details

Gives a flush cut to the wire at the distal end of the terminal tube Designed for larger wires
Tips are made of Tungsten Carbide for durability
Safety hold feature secures the loose wire after the archwire is cut

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

3M™ Unitek™ Eyelet, Large, 471-003, 1/Ea

3M ID 70202082734

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

3M™ Unitek™ T Spring

Details

.188 in. (4.8 mm) close-wound retraction spring Two configurations Coil loop at distal end 90° to helix

Typical properties

Arbor Diameter (Imperial): 0.03 in, 0.036 in, 0.3 in Arbor Diameter (Metric): 0.914 mm, 7.62 mm

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Units per Inner Pack: 10 per Pack

Wire Diameter (Imperial): 0.011 in, 0.012 in

Wire Diameter (Metric): 0.305 mm

704-045 Unitek™ Disposable Spatula 100/Pk

3M ID 7100154483

Details

Can be used for mixing liquid, powder, paste components Plastic Disposable 100/pack

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

3M™ AlastiK™ QuiK-StiK™ Separators, 406-410 Series

Details

Not made with natural rubber latex

Typical properties

Brands: AlastiK™

Factory ISO Certification: No ISO Certification information available

Product Color: Gray

eClass 14 Classification Group: 34110603

Dimensions and Classifications

Overall Diameter (Imperial): 0.145 in, 0.175 in Overall Diameter (Metric): 3.7 mm, 4.4 mm

3M™ Unitek™ Prestige Lingual Lock-Forming Pliers, 900-706, 1/Ea

3M ID 7000025661

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Package Contents: 1/Ea

3M™ Unitek™ Adjustable Sliding Yoke, 310-200, 10/Pk

Part Number 310-200 3M Product Number 310-200 3M ID 70202070275

Details

Clip onto archwire Adjust length then crimp 10/pk

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

3M™ Clinpro™ Sealant Refill, 12600 Series

Details

Color change from pink to white to help visualize placement

Convenient direct-delivery syringe for easy application

Low viscosity to flow easily into pits and fissures

Contains and releases fluoride — recharges with use of fluoride-containing

products

Bonds to enamel

Lower fill rate for fewer occlusal adjustments than filled sealants

Excellent stain resistance

Can be used with self-etch adhesive or phosphoric acid-etch technique

Typical properties

Brands: Clinpro™

Delivery System: Bottle, Syringe

Factory ISO Certification: No ISO Certification information available

Inner Packs per Case: 10 Packs per Case

Package Contents: Includes 1 - 6ml bottle of sealant, Includes 1 - 1.2ml sealant

syringe and 10 sealant syringe tips Total Units per Case: 10 per Case

Units per Case: 10.0

Units per Inner Pack: 1 per Pack

3M™ Clinpro™ Clear Fluoride Treatment

Details

Water-based fluoride treatment for improved taste, smell and mouthfeel Formula has a neutral pH and is free of rosin and common allergens* *See allergen disclosure below in the Resources section Available in a fast and easy-to-use L-Pop unit dose Effective fluoride uptake with lower dosage Minimum contact time of just 15 minutes Ready-release formula contains innovative stabilized calcium phosphate

Ready-release formula contains innovative stabilized calcium phosphate technology

Suction-line friendly formula

Available in mint, watermelon or flavorless

Typical properties

Brands: Clinpro™

Catalog Number: 7200W,7210W,7210FF,7200FF,7200M,7210M

Container Size (Imperial): 0.016 fl oz Container Size (Metric): 0.5 mL Flavor: Flavorless, Mint, Watermelon

Total Units per Case: 50 per Pack, 100 per Pack

3M™ Clarity™ Retainers

Details

Digitally remove appliances to create retainers for fixed appliance and aligner cases

All-digital approach — retainers can be ordered from a new scan or final aligner stage

Available for order following 3M[™] Clarity[™] Aligners, 3M[™] Digital Bonding, or 3M[™] Filtek Matrix cases

Highly durable clear retainers

Straight trimline contributes to patient comfort and esthetics

Nearly invisible, straight and stain resistant

Direct-to-patient delivery available for your convenience

Typical properties

Brands: Clarity™

3M™ Unitek™ Bracket Positioning Gauge

Details

Precision machined handles Color coded size

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

3M™ Unitek™ Molded Retainer Cases

Details

High-impact polypropylene Hinged sections Latchable Vented, to reduce odor formation

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Product Color: Assorted, Black, Blue, Green, Hot Green, Hot Orange, Hot Pink,

Hot Red, Hot Yellow, Ivory, Purple, Red, Yellow

3M™ Unitek™ Preformed Lingual Arch, Horizontal

Details

Permachrome wire (medium-hard) Adjustable (15° bilateral loops)

Typical properties

Arch Length (Imperial): 2.25 in, 2.36 in, 2.5 in, 2.63 in, 2.75 in, 2.87 in, 3 in, 3.11

in, 3.25 in, 3.38 in, 3.5 in

Arch Length (Metric): 57.15 mm, 59.944 mm, 63.5 mm, 66.802 mm, 69.85 mm,

72.898 mm, 76.2 mm, 78.994 mm, 82.55 mm, 85.852 mm, 88.9 mm

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Material: Stainless Steel

Units per Inner Pack: 10 per Pack

Wire Diameter (Imperial): 0.03 in, 0.036 in Wire Diameter (Metric): 0.762 mm, 0.9144 mm eClass 14 Classification Group: 34110602

3M[™] Concise[™] Orthodontic Chemical Cure Adhesive Enamel Bond Resin A, 196-009, 6 ml, 1/Btl

3M ID 7000024653 UPC 00652221047885

Details

One part of two-part resin and paste system Cures chemically Reliable bond strength

Typical properties

Adhesive Type: Direct Brands: Concise™

Catalog Number: 196-009

Container Size (Imperial): 0.203 fl oz

Container Size (Metric): 6 mL Cure Type: Chemical Cure

Factory ISO Certification: No ISO Certification information available

Product Type: Orthodontic Primer

Units per Inner Pack: 1 /BT

3M™ Unitek™ Coaxial Wire

Details

Light force
Highly flexible – easy to engage
Excellent initial wire
One form for both upper and lower arches
Ends do not unravel when cut

Typical properties

Arch: Universal Arch Shape: Square

Archwire Size (Imperial): 0.0155 in, 0.0175 in, 0.0195 in, 0.0215 in

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Material: Stainless Steel

Units per Inner Pack: 10 per Pack

3M™ Typodont with Stand, 3M™ Clarity™ Advanced Ceramic Brackets, 600-236, 1/Ea

3M ID 7100015672

Typical properties

Bracket Style: Low Profile

Brands: 3M™

Factory ISO Certification: No ISO Certification information available

3M™ AlastiK™ Color ID Dots, 406-829, 1/Pk

3M ID 7100154416

Details

Adhesive backed 1 each of 25 color dots

Typical properties

Brands: AlastiK™

Factory ISO Certification: No ISO Certification information available

Product Color: Assorted

Units per Inner Pack: 1 per Pack

3M™ Unitek™ Single-Ended Ligature Director, 803-003, 1/Ea

3M ID 7100268456

Details

Made from stainless steel Small, thin, serrated tip

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Units per Inner Pack: 1 per Pack

3M™ Unitek™ Molar Band Seater, 811-001, 1/Ea

3M ID 7000004390

Details

Hardened, colbalt-based alloy tip embedded in a molded handle which serves as a bite stick

Serrated, triangular tip allows better interproximal access with a positive grip

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

3M™ Unitek™ Lingual Mathieu Pliers

Details

Fine tips Available in curved and straight tips

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available Tip Orientation: Curved, Straight eClass 14 Classification Group: 21043800

3M™ Unitek™ Rectangular Retainer Cases, 452-001, 100/Pk

3M ID 70202082452

Details

Opaque plastic Upper and lower sections are hinged and latchable Compact size Can accommodate two retainers

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

3M™ Unitek™ Clear Plastic Button, 224-800, 10/Pk

3M ID 78000921896

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

3M™ Unitek™ Steiner Rotating Spring Wedge

Details

Wings made of Alloy 718, an age-hardened nickel-based aerospace alloy for more spring back

Centered notch on rolled-edge of one spring wing stabilizes wedge under archwire

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Units per Inner Pack: 10 per Pack

3M™ RelyX™ Universal Resin Cement Kits

Details

Universal dual-cure resin cement for virtually all adhesive and self-adhesive indications

15 applications on average from 3.4g cement

Micro Mixing Tip with minimized cement waste

Fast and easy excess clean-up

Unique self-sealing syringe - more hygienic, easier to clean compared to currently used automix syringes

Outstanding self-adhesive bond strength in both self-cure and light-cure Further enhanced bond strength to all substrates by 3M[™] Scotchbond[™] Universal Plus Adhesive

High esthetics by color-stable, fluorescent shades

Thin, long elongation tip for easy root canal application

Virtually no post-op sensitivity

Simplifies workflows and reduces stock - easier to train, easier to use, saves space, time and money

Typical properties

Brands: RelyX™

Catalog Number: 56970,56968,56969

Cure Type: Dual Cure

Delivery System: Automix Syringe

Factory ISO Certification: No ISO Certification information available

Material: Universal Resin

Package Contents: 1x 3M[™] RelyX[™] Universal Resin Cement syringe (3.4g), 15x 3M[™] RelyX[™] Universal Micro Mixing Tips, 5 Elongation Tips, 1 vial 3M[™] Scotchbond[™] Universal Plus Adhesive (1.5 ml), 1x 3M[™] Scotchbond[™] Universal Etchant (3 ml), 25 blue etchant tips, 50 orange microbrushes, step-by-step guides,4x 3M[™] RelyX[™] Universal Resin Cement syringes shades TR, A1, A3O, WO (3.4g each), 60x 3M[™] RelyX[™] Universal Micro Mixing Tips, 20 Elongation Tips, 1 vial 3M[™] Scotchbond[™] Universal Plus Adhesive (5 ml), 1x 3M[™] Scotchbond[™] Universal Etchant (3 ml), 25 blue etchant tips, 50 orange microbrushes, step-by-step guides

Shade: A1, Multiple, TR

3M™ Imprint™ 4 Preliminary Penta™ Cartridge for Pentamix™ 3, 71520 Series

Details

For consistent, homogeneous and void-free mixing at the touch of a button Compatible with 3M™ Pentamix™ 3 Mixing Unit

Reusable, metal-reinforced cartridge

Helps get an accurate impression on the first take, reducing the need for costly adjustments, remakes and retakes

Typical properties

Brands: Imprint™

Catalog Number: 71527,71528

Factory ISO Certification: No ISO Certification information available

Package Contents: 1 - Penta™ Cartridge

Product Type: Cartridge

3M™ Unitek ™ Anterior Band Regular Temper Narrow Upper Left/Right Central/Lateral

3M ID 70202032390

Typical properties

Brands: Unitek™

Factory ISO Certification: No ISO Certification information available

Product Type: Band

3M™ Filtek™ Easy Match Universal Restorative

Details

Simplified shade selection for anterior and posterior single-shade restorations Just three shades (Bright, Natural and Warm) to match the classical VITA shade guide

Naturally-adaptive opacity shows Enamel-like translucency, where you need it, on the bevel and incisal edge

No blocker is needed to mask the darkness of the oral cavity in anterior restorations

Excellent polish retention and shade match

Excellent wear resistance and strength

Firm handling

Featuring the same nanotechnology behind our leading 3M™ Filtek™ Dental Restoratives

Typical properties

Brands: Filtek™

Catalog Number: 6120W,6110W,6120B,6110B,6120N,6110N,6120K,6110K

Composite Type: Universal

Delivery System: Capsule, Syringe

Inner Packs per Case: 1 Bags per Case, 3 Kits per Case, 20 Bottles per Case, 60

Kits per Case

Material Structure: Nanocomposite Net Weight (Metric): 0.2 g, 4 g

Shade: Bright, Multiple Shades, Natural, Warm

3M™ RelyX™ Universal Resin Cement Refills

Details

Universal dual-cure resin cement for virtually all adhesive and self-adhesive indications

15 applications on average from 3.4g cement

Micro Mixing Tip with minimized cement waste

Fast and easy excess clean-up

Unique self-sealing syringe - more hygienic, easier to clean compared to currently used automix syringes

Outstanding self-adhesive bond strength in both self-cure and light-cure Further enhanced bond strength to all substrates by 3M[™] Scotchbond[™] Universal Plus Adhesive

High esthetics by color-stable, fluorescent shades

Thin, long elongation tip for easy root canal application

Virtually no post-op sensitivity

Simplifies workflows and reduces stock - easier to train, easier to use, saves space, time and money

Typical properties

Brands: RelyX™

Catalog Number: 56973,56974,56971,56972,56977,56978

Cure Type: Dual Cure

Delivery System: Automix Syringe

Factory ISO Certification: No ISO Certification information available

Material: Universal Resin

Package Contents: 1x 3M™ RelyX™ Universal Resin Cement syringe (3.4g), 15x

3M™ RelyX™ Universal Micro Mixing Tips, 5 Elongation Tips

Shade: A1, A3O, TR, White Opaque

3M™ Filtek™ One Bulk Fill Restorative

Details

Fast and easy one-step, bulk fill composite placement up to 5mm — no capping layer needed

Increased opacity for improved esthetics

Innovative monomers reduce shrinkage and shrinkage stress to enable up to 5mm depth of cure

Excellent adaptation, handling and sculptability

3M's TRUE nanofiller technology provides superior wear resistance and excellent polish retention

Capsules can be safely warmed up to 70°C/158°F for up to one hour Available in capsule or syringe delivery in 5 shades: A1, A2, A3, B1, C2 Complete your posterior composite restoration using the gold standard of adhesion with 3M™ Scotchbond™ Universal Plus Adhesive

Typical properties

Brands: Filtek™ Catalog Number:

4870A2,4871A1,4870A3,4871A2,4870C2,4870B1,4870A1,4871B1,4868A3,4866 B1,4866C2,4871A3,4866A3,4867A2,4868A1,4868A2,4867A3,4867B1,4867C2,4

866A1,4868C2,4868B1,4867A1,4866A2

Composite Type: Bulk Fill Cure Type: Light Cure

Delivery System: Capsule, Syringe Environmental claim being made: Yes Environmental claim relates to: Product Factory ISO Certification: ISO 13485 Inner Packs per Case: 10 Packs per Case

Material Structure: Nanocomposite Net Weight (Metric): 0.2 g, 4 g

Package Contents: Refill includes 20 - 02g capsules in shade A1,Refill includes 20 - 02g capsules in shade B1,Refill includes 20 - 02g capsules in shade C2,Refill includes 20 - 02g capsules in shade A3,Refill includes 20 - 02g capsules in shade A2,Refill includes 1 - 4g syringe in shade A2,Refill includes 1 - 4g syringe in shade B1,Refill includes 1 - 4g syringe in shade C2,Refill includes 1 - 4g syringe in shade A1

Shade: A1, A2, A3, B1, C2

Total Units per Case: 10 per Case, 200 per Case Units per Inner Pack: 1 per Pack, 20 per Pack

Units per Pack: 1, 20

3M™ Polycarbonate Crowns

Details

Polycarbonate material contours and crimps similar to metal crowns to yield easy adaptation to margins

Universal shade No 62 can be adjusted using cements and liners for better esthetics

Made from polycarbonate resin with microglass fibers for strength and durability A low cost alternative providing esthetic temporization for patients with economic considerations

Crown set is packaged in a unique, stackable tray that provides convenient, easy storage of 3M™ ESPE™ prefabricated crowns

Typical properties

Arch: Lower, Upper, Upper/Lower

Brands: 3M™ Catalog Number:

44,50,51,52,53,54,10,11,12,13,14,16,17,18,19,200,2,5,60,61,62,63,64,20,65,21,66,22,67,23,68,24,69,25,26,27,28,29,30,31,32,33,34,35,36,37,38,39,100,101,10 2,300,103,301,302,303,40,41,42,43

Dentition: Permanent

Factory ISO Certification: No ISO Certification information available

Material: Polycarbonate Resin

Mesial/Distal Width (Imperial): 0.157 in, 0.201 in, 0.205 in, 0.209 in, 0.213 in, 0.22 in, 0.224 in, 0.228 in, 0.244 in, 0.248 in, 0.252 in, 0.256 in, 0.268 in, 0.272 in, 0.276 in, 0.283 in, 0.283 in, 0.295 in, 0.299 in, 0.303 in, 0.307 in, 0.315 in, 0.319 in, 0.323 in, 0.327 in, 0.331 in, 0.339 in, 0.35 in, 0.354 in, 0.37 in, 0.398 in Mesial/Distal Width (Metric): 4 mm, 5.1 mm, 5.2 mm, 5.3 mm, 5.4 mm, 5.6 mm, 5.7 mm, 5.8 mm, 6.2 mm, 6.3 mm, 6.4 mm, 6.5 mm, 6.8 mm, 6.9 mm, 7 mm, 7.1 mm, 7.2 mm, 7.5 mm, 7.6 mm, 7.7 mm, 7.8 mm, 8 mm, 8.1 mm, 8.2 mm, 8.3 mm, 8.4 mm, 8.6 mm, 8.9 mm, 9 mm, 9.4 mm, 10.1 mm

Number of Crowns: 5

Side: Left, Left/Right, Right

Size: 10, 100, 101, 102, 103, 11, 12, 13, 14, 15, 16, 17, 18, 19, 2, 20, 200, 21, 22, 23, 24, 25, 26, 27, 28, 29, 30, 300, 301, 302, 303, 31, 32, 33, 34, 35, 36, 37, 38, 39, 40, 41, 42, 43, 44, 50, 51, 52, 53, 54, 60, 61, 62, 63, 64, 65, 66, 67, 68, 69

Tooth: Central Incisor, Cuspid, First Bicuspid, Incisor, Lateral Incisor, Second Bicuspid

3M™ Vitremer™ Core Buildup/Restorative, 3303

Details

Versatile system can be light-cured or self-cured to meet your application needs Unique dark-cure chemistry allows for bulk placement which eliminates the need for layering and saves time

Displays a true glass ionomer cure for long-term fluoride release and adhesion to tooth structure which provides better patient protection and increased doctor confidence

As a restorative material, it provides good handling properties and extended working time in the mouth for easier use

Available in a variety of shades to better match existing tooth structure for esthetic restorations. Pedo shade is available for pedodontic applications and Blue shade for core buildup applications

Typical properties

Brands: Vitremer™ Catalog Number:

3303A3.5,3303PEDO,3303C4,3303B3,3303A3,3303A4,3303BLUE,3303C2,3303

B2,3303L

Cure Type: Light Cure

Delivery System: Bottle, Hand Mix

Factory ISO Certification: No ISO Certification information available

Inner Packs per Case: 10 Packs per Case

Package Contents: 9g powder in shade B3, Spoon,9g powder in shade Blue, Spoon,9g powder in shade A4, Spoon,9g powder in shade C2, Spoon,9g powder in shade B2, Spoon,9g powder in shade A35, Spoon,9g powder in shade C4, Spoon,Refill includes 8ml liquid,9g powder in shade A3, Spoon,9g powder in shade Pedo, Spoon

Shade: A3, A3.5, A4, B2, B3, Blue, C2, C4, PEDO

Total Units per Case: 10 per Case Units per Inner Pack: 1 per Pack

3M™ RelyX™ Unicem Self-Adhesive Universal Resin Cement Maxicap™ Capsule Refill

Details

Saves time by eliminating etching, priming and bonding steps
Virtually no post-operative sensitivity
Less discoloration
Strong and moisture-tolerant
If you like 3M™ RelyX™ Unicem Self-Adhesive Resin Cement, you'll love 3M™
RelyX™ Universal Resin Cement

Typical properties

Brands: RelyX™ Catalog Number:

56820,56831,56832,56833,56816,56817,56818,56819,56834,56815

Cure Type: Dual Cure

Delivery System: Maxicap™ Capsule

Factory ISO Certification: No ISO Certification information available

Material: Adhesive Resin, Self-adhesive Resin

Package Contents: 20 – 3M[™] Maxicap[™] Capsules A2 Universal,50 – 3M[™] Aplicap[™] Capsules Translucent; 10 Elongation Tips,50 – 3M[™] Aplicap[™] Capsules A2 Universal; 10 Elongation Tips,20 – 3M[™] Maxicap[™] Capsules Translucent,20 – 3M[™] Aplicap[™] Capsules White Opaque,20 – 3M[™] Aplicap[™] Capsules A1,20 – 3M[™] Maxicap[™] Capsules A3 Opaque,50 – 3M[™] Aplicap[™] Capsules A3 Opaque; 10 Elongation Tips,20 – 3M[™] Maxicap[™] Capsules in Assorted Shades: 10 – A2 Universal, 6 – Translucent, 4 – A3 Opaque,50 – 3M[™] Aplicap[™] Capsules in Assorted Shades: 24 – A2 Universal, 16 – Translucent, 10 – A3 Opaque; 10 Elongation Tips

Shade: A1, A2, A3, Multiple, TR, White Opaque

3M™ ESPE™ Stainless Steel Primary Molar Crowns

Details

Pre-belled, designed for quick, easy-snap fit
Superior longevity compared to multi-surface amalgam fillings
Accurately duplicates anatomy of primary molars
Crown set is packaged in a stackable tray for storage of 3M™ prefabricated crowns
Wide array of shapes and sizes
High quality stainless steel

Typical properties

Arch: Lower, Upper Brands: 3M™

Catalog Number: D-LL-3,D-LL-2,D-LR-5,D-LR-4,D-LR-3,D-LR-2,D-LL-7,D-LL-6,D-LL-5,D-LR3,D-LL-4,D-LR-7,D-LR-6,D-UR-2,D-UL-7,D-UL-6,D-UL-5,D-UR-6,D-UR-5,D-UR-4,D-UR-3,D-UL-3,D-UL-2,D-UR-7,E-UL-2,E-LR-7,E-UL-3,E-UL-4,E-UL-5,E-LR-3,E-UL-6,E-LR-4,E-UL-7,E-LR-5,E-UR-2,E-LR-6,E-UR-3,E-LL-5,E-UR-4,E-LL-6,E-UR-5,E-LL-7,E-UR-6,E-LR-2,E-UR-7,E-LL-2,E-LL-3,E-LL-4

Dentition: Primary

Factory ISO Certification: No ISO Certification information available Mesial/Distal Width (Imperial): 0.283 in, 0.287 in, 0.299 in, 0.303 in, 0.315 in, 0.319 in, 0.331 in, 0.335 in, 0.346 in, 0.35 in, 0.362 in, 0.366 in, 0.37 in, 0.378 in, 0.386 in, 0.394 in, 0.402 in, 0.409 in, 0.417 in, 0.425 in, 0.433 in, 0.441 in, 0.449 in

Mesial/Distal Width (Metric): 7.2 mm, 7.3 mm, 7.6 mm, 7.7 mm, 8 mm, 8.1 mm, 8.4 mm, 8.5 mm, 8.8 mm, 8.9 mm, 9.2 mm, 9.3 mm, 9.4 mm, 9.6 mm, 9.8 mm, 10 mm, 10.2 mm, 10.4 mm, 10.6 mm, 10.8 mm, 11 mm, 11.2 mm, 11.4 mm

Number of Crowns: 5 Patient: Pediatric Side: Left, Right Size: 2, 3, 4, 5, 6, 7

Tooth: First Molar, Second Molar

Units per Case: 5.0

3M™ Ketac™ Conditioner Refill, 37470

3M ID 7000054695 UPC 04035077013192

Details

Removes the dentin smear layer generated during the tooth preparation

Typical properties

Brands: Ketac™

Catalog Number: 37470 Delivery System: Vial

Factory ISO Certification: No ISO Certification information available

Package Contents: Refill includes 10ml liquid eClass 14 Classification Group: 34110000

3M™ Imprint™ 4 VPS Bite Registration Material, 71529

3M ID 7100196326 UPC 04054596648333

Details

60-second setting time

Convenient working time reduces chair time for patients and stress for dentists Great flow properties and no-slump consistency for accurate impressions Suitable hardness, low flexibility and easy to cut and trim after setting Can be scanned for occlusal recording and registration with CAD/CAM systems For successful indirect procedures, capture your precision impressions with 3M[™] Imprint[™] 4 VPS Impression Material or 3M[™] Impregum[™] Impression Material

Typical properties

Brands: Imprint™

Catalog Number: 71529

Delivery System: Garant™ Dispenser

Factory ISO Certification: No ISO Certification information available

Impression Type: Bite Registration

Material: VPS

Package Contents: 2 - 50ml cartridges; 5 - 3M™ Garant™ Mixing Tips - Green; 5

- 3M™ Garant™ Intra-oral Tips - Bite

Product Color: Sky Blue

Set Time: 1 min

Set Type: Regular Set Setting Time: 1 Minute Work Time: 20 Seconds

3M™ Pentamix™ 3 Wall Mount Kit, 77601

3M ID 7100215176 UPC 04054596736696

Typical properties

Brands: Pentamix™ Catalog Number: 77601

Factory ISO Certification: No ISO Certification information available

Package Contents: 1 - wall mount kit eClass 14 Classification Group: 34110100

3M™ Unitek™ Stainless Steel Permanent Bicuspid Crowns, 902600, Kit

3M ID 7000054529 UPC 30605861012851

Details

Easy to customize with parallel walls to allow for clinician fit control Pre-trimmed to optimum length and contour to save placement time Made in the U.S.A. of high quality stainless steel

Typical properties

Arch: Upper/Lower Brands: Unitek™

Catalog Number: 902600 Dentition: Permanent

Factory ISO Certification: No ISO Certification information available

Number of Crowns: 84 Patient: Pediatric Size: Assorted

Tooth: Assorted, Bicuspid

3M[™] Filtek[™] Z250 XT Nano Hybrid Universal Restorative Introductory Kit, 1470IK, Syringe

3M ID 7100156383 UPC 30605861066960

Details

Visible light-activated nanohybrid composite
Used to permanently bond the composite to the tooth structure
Available in 12 shades, two of which are Opaque
All shades are radiopaque and fluorescent
Packaged in syringes and capsules

Typical properties

Brands: Filtek™

Catalog Number: 1470IK Composite Type: Universal Delivery System: Syringe

Factory ISO Certification: No ISO Certification information available

3M™ Penta™ Mixing Tips Refill, Red (50), 77949

3M ID 7100215286

Details

Unique design reduces dispensing forces by up to 50 percent over previous mixing tip to facilitate easier extrusion of impression materials Homogeneous mixing for void-free impressions
Less waste--dispenses only amount of material needed
Can be used with polyether, VPS and alginate replacement materials

Typical properties

Brands: Penta™

Catalog Number: 77949

Factory ISO Certification: No ISO Certification information available

Package Contents: 50 - Mixing Tips - Red

Product Type: Tip Units per Case: 50.0

3M™ Elipar™ DeepCure-S Li Ion Battery, 76985

Part Number 76985 3M Product Number 76985 3M ID 7100196365

Details

A new fully-charged battery provides 60 minutes of cure time Indicator LED shows level of battery's charge Battery can be fully recharged in 90 minutes

Easy tool-free battery exchange

Simple two-button operation with one cure mode for all indications

Tack cure function allows for predictable curing with a reproducible short light output

Silent operation -- innovative design eliminates noisy fan, and includes option to turn off beep signals

Two light guides, both offering consistently high and focused light output

Typical properties

Brands: Elipar™

Catalog Number: 76985

Factory ISO Certification: No ISO Certification information available

Package Contents: 1 - Elipar DeepCureS ion Battery

Product Type: Battery

3M™ Vanish™ XT Extended Contact Varnish, 12148, 10 g Clicker™ Dispenser

Part Number 12148 3M Product Number 12148 3M ID 7000128831

Details

Delivers fluoride, calcium and phosphate

Creates a barrier against demineralization and acid erosion* *Based on in vitro

Helps repair demineralized enamel* *Based on in vitro data

Provides protection beyond the coating

Coating remains for 6 months or longer

Moisture-tolerant

Typical properties

Brands: Vanish™

Factory ISO Certification: No ISO Certification information available

Inner Packs per Case: 10 Packs per Case Material: Resin-modified Glass Ionomer

Total Units per Case: 10 per Case Units per Inner Pack: 1 per Pack

3M™ Express™ VPS Bite Registration Material Putty, 6160J

3M ID 7000053382 UPC 30605861010871

Details

Quick setting time High rigidity Mixes easily

Typical properties

Brands: Express™
Catalog Number: 6160J
Delivery System: Hand Mix

Factory ISO Certification: No ISO Certification information available

Impression Type: Bite Registration

Material: VPS

Package Contents: 1 - jar of base, 305ml; 1 - jar of catalyst, 305 ml

Product Color: Green
Set Time: 1.5 Minutes
Set Type: Quick Set
Setting Time: 1.5 Minute

Work Time: 30 sec

3M™ ESPE™ Express™ XT Putty Soft VPS Impression Material Refill, 36972

3M ID 7000042952 UPC 04035077003926

Details

Hydrophilicity and flow behavior helps result in void-free impressions
Elastic recovery helps prevent distortion for a more accurate impression
Resists tearing upon removal thanks to material toughness
Recovery from deformation helps result in fewer retakes
Available in a full range of viscosities, deliveries and setting times
Control moisture at every step and gently retract gingival tissue for more accurate impressions with 3M™ Astringent Retraction Paste

Typical properties

Brands: Express™

Catalog Number: 36972

Factory ISO Certification: No ISO Certification information available

Material: VPS

Set Type: Regular Set

Viscosity: Multiple Viscosities

3M™ Scotchbond™ Dispensing Tips, 41274

3M ID 7100156110 UPC 04054596539679

Details

The tip can be individually adjusted before application

Typical properties

Brands: Scotchbond™

Factory ISO Certification: No ISO Certification information available

Product Type: Tip

3M™ Polyether Adhesive Refill, 30600

3M ID 7100196367 UPC 04054596648616

Details

Ensures adhesion between the polyether impression material and the impression tray

Typical properties

Brands: 3M™

Catalog Number: 30600 Cure Type: Air Dry Delivery System: Vial

Factory ISO Certification: No ISO Certification information available

Material: Polyether

Package Contents: 1 - 17ml vial Product Type: Tray Adhesive

По вопросам продаж и поддержки обращайтесь:

Алматы (7273)495-231 Ангарск (3955)60-70-56 Архангельск (8182)63-90-72 Астрахань (8512)99-46-04 Барнаул (3852)73-04-60 Белгород (4722)40-23-64 Благовещенск (4162)22-76-07 Брянск (4832)59-03-52 Владивосток (423)249-28-31 Владикавказ (8672)28-90-48 Владимир (4922)49-43-18 Волгоград (844)278-03-48 Вологда (8172)26-41-59 Воронеж (473)204-51-73 Екатеринбург (343)384-55-89 Иваново (4932)77-34-06 Ижевск (3412)26-03-58 Иркутск (395)279-98-46 Казань (843)206-01-48 Калининград (4012)72-03-81 Калуга (4842)92-23-67 Кемерово (3842)65-04-62 Киров (8332)68-02-04 Коломна (4966)23-41-49 Кострома (4942)77-07-48 Краснодар (861)203-40-90 Красноярск (391)204-63-61 Курск (4712)77-13-04 Курган (3522)50-90-47 Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13 Москва (495)268-04-70 Мурманск (8152)59-64-93 Набережные Челны (8552)20-53-41 Нижний Новгород (831)429-08-12 Новокузнецк (3843)20-46-81 Ноябрьск (3496)41-32-12 Новосибирск (383)227-86-73 Омск (3812)21-46-40 Орел (4862)44-53-42 Оренбург (3532)37-68-04 Пенза (8412)22-31-16 Петрозаводск (8142)55-98-37 Псков (8112)59-10-37 Пермь (342)205-81-47 Ростов-на-Дону (863)308-18-15 Рязань (4912)46-61-64 Самара (846)206-03-16 Санкт-Петербург (812)309-46-40 Саратов (845)249-38-78 Севастополь (8692)22-31-93 Саранск (8342)22-96-24 Симферополь (3652)67-13-56 Смоленск (4812)29-41-54 Сочи (862)225-72-31 Ставрополь (8652)20-65-13 Сургут (3462)77-98-35 Сыктывкар (8212)25-95-17 Тамбов (4752)50-40-97 Тверь (4822)63-31-35

Тольятти (8482)63-91-07 Томск (3822)98-41-53 Тула (4872)33-79-87 Тюмень (3452)66-21-18 Ульяновск (8422)24-23-59 Улан-Удэ (3012)59-97-51 Уфа (347)229-48-12 Хабаровск (4212)92-98-04 Чебоксары (8352)28-53-07 Челябинск (351)202-03-61 Череповец (8202)49-02-64 Чита (3022)38-34-83 Якутск (4112)23-90-97 Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +7(7172)727-132

Киргизия +996(312)96-26-47